For use in conjunction with the FastPack® IP and FastPack® Total PSA Immunoassay and the FastPack® IP System Analyzer

**CAUTION: United States Federal law restricts this device to sale and distribution by or on the order of a physician, or to a clinical laboratory; and use is restricted to, by or on the order of a physician.**

### INTENDED USE

The FastPack® Total PSA Calibrator is intended to calibrate the FastPack® IP System when used for the quantitative determination of total PSA in human serum and plasma.

### SUMMARY AND PRINCIPLE

Quantitative assay calibration is the process by which a set of samples with known analyte concentrations (i.e. assay calibrators) is tested to measure the response. The mathematical relationship between the measured responses and the known analyte concentrations establishes the calibration curve. The FastPack® IP Analyzer must be calibrated by the user to ensure that it is properly adjusted for the particular lot of FastPacks that is being used. The FastPack® Total PSA Calibrator is used for this purpose.

### PRODUCT INFORMATION

FastPack® Total PSA Calibrator - Cat. No. 25000002

- Provided ready to use.
- Mix contents by gently inverting before use. Avoid bubble formation.
- Calibrator: 5 mL/vial. Liquid. Contain components of human origin prepared in a Tris Buffer with protein stabilizers to yield predetermined concentrations.  
Calibrator B
- Calibrator Card: 1
- Preservatives: 0.1% Sodium Azide

### WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Do not mix calibrators from different lots.
- After opening, calibrators are stable until the expiration date on the label when stored and handled as directed. Do not use calibrators beyond the expiration date.
- Avoid microbial contamination of reagent when removing aliquots from the bottles.
- Refer to the FastPack® IP Procedure Manual for calibration procedures.
- Discard unused or expired calibrator material, in stoppered vial, into a Biohazard container.
- Reagents in this kit contain sodium azide as a preservative, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build up.
- The components containing Sodium Azide are classified per applicable European Economic Community (EEC) Directives as: Very toxic and dangerous to the environment (T+ N). The following are the appropriate Risk (R) and Safety (S) phrases:

R28	Very toxic if swallowed.
R32	Contact with acids liberates very toxic gas.
R50/53	Very toxic to aquatic organisms may cause long-term adverse effects in the aquatic environment.
S28	After contact with skin, wash immediately with plenty of soap-suds.
S45	In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
S60	This material and its container must be disposed of as hazardous waste.
S61	Avoid release to the environment. Refer to special instructions/safety data sheets.
- **Human source material. Treat as potentially infectious.**

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© 2000 Qualigen, Inc. All rights reserved. Qualigen and FastPack are trademarks or registered trademarks of Qualigen, Inc.  
All other trademarks are the property of their respective owners.

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