## AMYLASE LIQUID SELECT

**CATALOGUE NUMBER:** 80-5383-00 **SIZE:** 10 x 5 mL  
80-5451-00 1 x 1 L

### INTENDED USE

For the quantitative measurement of  $\alpha$ -amylase ( $\alpha$ 1, 4-glucoan-4-glucohydrolase, E.C.3.2.1.1) activity in serum and plasma.

### TEST SUMMARY


For many years, the levels of serum  $\alpha$ -amylase in patients have provided needed evidence for the diagnosis of acute pancreatitis<sup>1-3</sup>. Early assay techniques were based on either a change in the absorption maxima of the complex between starch and iodine as the  $\alpha$ -amylase degraded the starch; or a measurement of the increase in reducing groups as the starch was hydrolyzed by the  $\alpha$ -amylase<sup>4</sup>. These methods are not as reliable and easy to quantitate as spectrophotometric methods using a defined substrate<sup>5</sup>.

Some methods are based on the production of NADH proportionate to the activity of the  $\alpha$ -amylase. A defined substrate, such as maltotetraose, is degraded by  $\alpha$ -amylase to produce glucose which can be measured in a coupled enzyme assay. However, this method necessitates the removal of endogenous glucose which would give a high background to the assay<sup>5</sup>.

More recent methods are based on the production of p-nitrophenol from defined oligosaccharide substrates with blocking groups attached on the terminal sugar. The action of the  $\alpha$ -amylase on the oligosaccharide yields a variety of chain lengths after hydrolysis. These methods then use a variety of coupling enzymes to hydrolyze the resulting short chain oligosaccharides to produce p-nitrophenol<sup>6</sup>. The coupling enzymes contain residual  $\alpha$ -amylase activity that may significantly reduce the stability of the reagent.

### TEST PRINCIPLE

The Amylase Liquid Select assay involves the use of a chromogenic substrate, 2-chloro-4-nitrophenol linked with maltotriose<sup>7</sup>.

As shown above,  $\alpha$ -amylase hydrolyzes the 2-chloro-4-nitrophenyl- $\alpha$ -D-maltotriose (CNPG<sub>3</sub>) to release 2-chloro-4-nitrophenol (CNP) and form 2-chloro-4-nitrophenyl- $\alpha$ -D-maltoside (CNPG<sub>2</sub>), maltotriose (G<sub>3</sub>), and glucose (G). The rate of formation of the 2-chloro-4-nitrophenol can be detected spectrophotometrically at 405 nm to give a direct measurement of  $\alpha$ -amylase activity in the sample. The reaction is not readily inhibited by endogenous factors.

### REAGENTS

Component	Ingredients	Concentration
Amylase Liquid Select Reagent	2-Chloro-4-Nitrophenyl- $\alpha$ -D-Maltotriose (CNPG <sub>3</sub> )	2.25 mM
	Sodium Chloride	350 mM
	Calcium Acetate	6 mM
	Potassium Thiocyanate	900 mM
	Sodium Azide	0.1 %
	Buffer	pH 6.0

### WARNINGS AND PRECAUTIONS FOR USE

#### WARNING

H319 – Causes serious eye irritation

#### Prevention

P264 – Wash thoroughly after handling.

P280 – Wear eye protection/face protection.

#### Response

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

P337 + P313 If eye irritation persists, get medical advice/attention.

For In Vitro Diagnostic Use.

See Material Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

The Amylase Liquid Select reagent is supplied ready to use.

Unopened reagent is stable until the expiration date printed on the label when stored at 2-8°C.

After opening, the reagent is stable for 60 days when properly capped immediately after each use and stored at 2-8°C.

Reagent is stable for 14 days on board the Roche/Hitachi 911 Analyzer.

### REAGENT DETERIORATION

The following results indicate deterioration:

- Inability to recover control values.
- Presence of turbidity.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Serum, sodium heparinized plasma or lithium heparinized plasma are the recommended sample types. Other anti-coagulants such as EDTA or citrate should not be used.

Refer to CLSI GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture<sup>8</sup> and CLSI GP44-A4 Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests<sup>9</sup> for guidance.

If not analyzed promptly, internal studies have shown that samples may be stored ambient (20-25°C) for 14 days, refrigerated (2-8°C) for 21 days or frozen (-20°C) for 10 weeks. Samples may be frozen (-20°C), in a non-frost (frost free) freezer, and thawed for up to two freeze/thaw cycles.

### ANALYTICAL SPECIFICITY

- Refer to the work of Young et al<sup>10</sup> for a review of drug effects on amylase levels.
- Amylase enzyme activity is temperature dependent.<sup>4</sup>
- Macroamylase has been shown to cause hyperamylasemia which may lead to overdiagnosis of acute pancreatitis when using oligosaccharide substrates.<sup>11</sup>
- No significant interference was detected in the Amylase Liquid Select assay up to and including the concentrations stated using the Roche/Hitachi 911 Analyzer.

<b>Substance Tested</b>	<b>Concentration with no significant (<math>\pm 10\%</math>) Interference</b>
Liposyn	3000 mg/dL (1%)
Triglyceride	3000 mg/dL
Ascorbic Acid	50 mg/dL
Bilirubin (unconjugated)	50 mg/dL
Bilirubin (conjugated)	50 mg/dL
Hemoglobin	500 mg/dL
Glucose	2000 mg/dL

Samples with hemoglobin interference higher than the upper limit may be diluted 1 part sample with 1 part physiological saline. Multiply the result by two to correct for the dilution.

#### **ANALYTICAL PROCEDURE**

#### **MATERIALS PROVIDED**

The Amylase Liquid Select Reagent is packaged in various formats. Any of the following items may be included in the package you receive.

<b>Description</b>	<b>Configuration</b>	<b>Catalog Number</b>
Amylase Liquid Select Reagent	1 x 1 L	80-5451-00
Amylase Liquid Select Reagent	10 x 5 mL	80-5383-00

#### **MATERIALS REQUIRED (BUT NOT PROVIDED)**

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
2. Calibration material. (If applicable.)
3. Quality Control materials.

#### **TEST CONDITION**

For the data presented in this insert, studies were performed on an automated analyzer using a rate test mode, with a sample to reagent ratio of 1:40 and a wavelength of 415 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at 800-565-0265. Outside Canada and the U.S. please contact our local distributor.

#### **CALIBRATION**

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

#### **QUALITY CONTROL**

Reliability of test results should be routinely monitored with quality control materials or serum pools that reasonably represent performance on patient specimens. Controls or serum pools should be run with each assay to ensure that the reagents are functioning properly and that correct procedures have been followed. Quality control materials are intended for use only as monitors of accuracy and precision. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range confirm procedures were performed correctly and follow normal troubleshooting measures. If the problem persists call Sekisui Diagnostics Technical Services at 800-565-0265. Outside Canada and the U.S., please contact your local distributor.

Quality control requirements should be established in accordance with local, state and/or Federal regulations or accreditation requirements.

#### **RESULTS**

The analyzer automatically calculates the  $\alpha$ -amylase activity of each sample.

#### **REFERENCE INTERVALS**

The expected range for Amylase is 20 – 104 U/L.<sup>12</sup> It is recommended that each laboratory establish the normal range for its population.

#### **PERFORMANCE CHARACTERISTICS**

##### **REPORTABLE RANGE**

The Amylase Liquid Select Reagent is linear to 2000 U/L on the Roche/Hitachi 911 Analyzer. If a sample exceeds 2000 U/L, it should be diluted with an equal volume of normal saline and re-assayed. Multiply the value from the resulting calculation by 2 to correct for the dilution.

##### **PRECISION STUDIES**

Within run precision of the Amylase Liquid Select reagent was determined using four levels of frozen serum pools analyzed 20 times in a single run on the Roche/Hitachi 911 Analyzer.

<b>Serum Pool</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
n	20	20	20	20
Mean (U/L)	61	272	902	1509
% CV	1.34	0.61	0.51	0.62

Between run precision of the Amylase Liquid Select reagent was determined using four levels of frozen serum pools. Each pool was run in duplicate for 20 times over at least five days on the Roche/Hitachi 911 Analyzer.

<b>Serum Pool</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
n	40	40	40	40
Mean (U/L)	60	273	917	1507
% CV	1.7	0.8	0.9	0.6

#### **ACCURACY**

In comparative performance studies, 50 samples ranging in activity from 28-304 U/L were evaluated by Amylase Liquid Select Reagent and another commercially available  $\alpha$ -amylase method on the Roche/Hitachi 911 analyzer. The correlation coefficient for the study was  $r = 0.9985$ . The regression equation was  $y = 0.90x - 2.50$ .

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

#### **TRADEMARKS**

The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

All trademarks, brands, product names and trade names are the property of their respective companies.

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©2011 Sekisui Diagnostics, LLC - All rights reserved.
©2011 Sekisui Diagnostics, LLC - Tous droits réservés.