(Diagnostic Reagent Grade)


PHOSPHOLIPASE D

[PLDP (T-39), PLDPV (T-138)]

(Glycerophospholipid specific)

from *Streptomyces* sp.

(Phosphatidylcholine phosphatidohydrolase, EC 3.1.4.4)

Preparation and Specification

PLDP (T-39)

Appearance : Light brownish amorphous powder, lyophilized
 Specific activity : More than 100 U/mg solid

PLDPV (T-138)

Appearance : White to brownish lyophilized powder
 Specific activity : More than 30 U/mg solid
 ※ Animal derived material free.

Properties

| | | |
|-----------------------------|-------------------------------------|----------|
| Substrate specificity | : See Table 1 | |
| Molecular weight | : 46 kDa (gel filtration) | |
| Isoelectric point | : pH 4.2 | |
| Optimum pH | : 5.5-6.0 | Figure 1 |
| pH stability | : 4.2-8.5 (37°C, 60 min, 0.05% BSA) | Figure 2 |
| Storage stability | : At least one year at -20°C | Figure 3 |
| Effect of various chemicals | : See Table 2 and Table 3 | |

Transphosphatidylation Catalyzed by Phospholipase D

Figure 4
 Figure 5
 Figure 6
 Figure 7

Table 1. Substrate specificity


| Substrate | Relative activity (%) |
|---------------|-----------------------|
| Lecithin | 100 |
| Lysolecithin | 3.4 |
| Sphingomyelin | 0.03 |

Table 2. Effect of detergents on PLDP activity


| Additive | Concentration (%) | Relative activity (%) |
|--------------|-------------------|-----------------------|
| None | 0 | 0 |
| Triton X-100 | 0.1 | 28 |
| | 0.5 | 100 |
| | 1.0 | 78 |
| | 1.0 | 15 |
| SDS | 0.1 | 8 |
| Deoxycholate | 0.1 | 8 |

Table 3. Effect of metal ions on PLDP activity


| Metal ion | Concentration (mM) | Relative activity (%) |
|-------------------|--------------------|-----------------------|
| None | 0 | 100 |
| CaCl ₂ | 10 | 99 |
| MgCl ₂ | 10 | 101 |
| MnCl ₂ | 10 | 132 |
| ZnCl ₂ | 1 | 99 |
| CoCl ₂ | 1 | 106 |
| BaCl ₂ | 1 | 101 |
| CuCl ₂ | 1 | 65 |
| EDTA | 1 | 99 |
| | 10 | 100 |

- : 3,3-Dimethylglutarate-NaOH buffer
- : Phosphate buffer
- △ : Tris-HCl buffer
- : Glycine-NaOH buffer

- 37°C, 60 min.
- : 3,3-Dimethylglutarate-NaOH buffer
 - : Phosphate buffer
 - △ : Tris-HCl buffer
 - : Glycine-NaOH buffer

-20°C

Transphosphatidyltion Catalyzed by Phospholipase D




Fig.4 The Influence of Organic Solvents on the Yield of 5'-Phosphatidylcytidine by PLDP-Catalyzed Transphosphatidyltion


Organic phase : Solvent (0.5ml) Contained 1mg Dipalmitoylphosphatidylcholine
 Aqueous phase : 0.25M Cytidine-HCl buffer pH4.5 (0.1ml) and 0.05ml of water contained 0.9 Unit of PLDP
 Reaction Temp : 35°C, Reaction : 90 min.


Fig.5 Optimum pH of PLDP-Catalyzed

Donor : Uridine 250 mM pH6-7 : MES-NaOH buffer
pH4-6 : Acetate buffer

pH3.5-4 : Glycine-HCl buffer Cytidine, Deoxycytidine : 250mM
pH4-5.5 : Acetate buffer

pH3-4 : Glycine-HCl buffer Adenosine : 50mM
pH4-6 : Acetate buffer Deoxyadenosine : 50mM

Phospholipid acceptor : Dipalmitoylphosphatidylcholine in Chloroform (1ml)
Buffer : 0.5ml contained 1.8 Units of PLDP
Reaction Temp : 35°C
Reaction period : 90 min.

Fig.6 The Influence of Reaction Temperature on the Yield of 5'-Phosphatidylcholine by PLDP-Catalyzed Transphosphatidylation (Donor:dipalmitoylphosphatidylcholine, Acceptor: Cytidine, pH4.5)



Fig.7 The Effect of Acceptor-Donor Ratio on the Yield of 5'-Phosphatidylcytidine by PLDP-Catalyzed Transphosphatidylation




Organic solvent : 10 mM Dipalmitoylphosphatidylcholine in Chloroform (1ml)
Aqueous solvent : 0.1 M Glycine-HCl (pH4.5) 0.5ml containing 1.8 Units of PLDP
Reaction : at 35°C, with continuous stirring, for 180min.

Assay

Principle

The assay is based on the increase in absorbance at 500 nm as the formation of quinoneimine dye proceeds in the following reactions:

Unit definition

One unit is defined as the amount of enzyme which hydrolyzes 1 μmole of phosphatidylcholine to phosphatidic acid and choline per minute at 37°C under the conditions specified in the assay procedure.

■ Reagents

1. Reaction mixture for the first reaction

| | |
|--|---------|
| 0.2 M DMG-NaOH buffer pH 5.5 | 0.10 ml |
| 10 mM CaCl ₂ | 0.05 ml |
| 3% (W/V) Triton X-100 solution | 0.05 ml |
| 10 mM Substrate solution ¹⁾ | 0.10 ml |
| Distilled water | 0.15 ml |

DMG: 3, 3'-Dimethylglutarate
1): 10 mM Substrate solution
Dissolve 78.6 mg of L- α -phosphatidylcholine, dioleoyl with 10 ml of 5% Triton X-100 (W/V).
2. Reaction mixture for the second reaction

| | |
|-----------------------------|---------|
| 15 mM 4-AA | 0.05 ml |
| 0.2% (W/V) Phenol | 0.05 ml |
| 1 M Tris-HCl buffer, pH 8.0 | 0.05 ml |
| 50 U/ml POD ²⁾ | 0.05 ml |
| 50 U/ml COD ³⁾ | 0.05 ml |
| Distilled water | 0.25 ml |

2): 50 U/ml POD
Dissolve 500 U (PPU) of POD with 10 ml of distilled water.
- 3): 50 U/ml COD
Dissolve 500 U of COD with 10 ml of 10 mM Tris-HCl buffer, pH 8.0.
3. Reaction stopper
1 M Tris-HCl buffer containing 10 mM EDTA and 1% (W/V) Cetyltrimethylammonium chloride
EDTA: Ethylenediamine tetraacetic acid
4. Reaction dilution solution
1% (W/V) Triton X-100
5. Enzyme dilution buffer
10 mM DMG-NaOH buffer, pH 5.5 containing 0.1% Triton X-100
6. Reagents:
DMG : Tokyo Kasei Kogyo Co., Ltd. #D1322
Triton X-100: The Dow Chemical Company
1,2-Dioleoyl-sn-glycero-3-phosphocholine:
Sigma Chemical Co. #P-6354
EDTA (2Na · 2H₂O) : KISHIDA CHEMICAL Co., Ltd.
#060-29133
COD: Asahi Kasei Pharma Corporation #T-05
4-AA: NACALAI TESQUE, INC. Special grade #01907-52
Cetyltrimethylammonium chloride:
Wako Pure Chemical Industries, Ltd. #087-06032
POD: Sigma Chemical Co. Type II #P-8250

■ Enzyme solution

Weigh about 20 mg of test sample exactly and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

■ Procedure

1. Pipette 0.45 ml of reaction mixture for the first reaction into a small test tube and preincubate them at 37°C.
2. After 5 min, add exactly 50 μ l of enzyme solution and mix to start the first reaction.
※ In the case of a test blank, add 50 μ l of enzyme dilution buffer in place of enzyme solution at this point.
3. After 10 min, add 0.50 ml of reaction stopper and mix. On stopping the first reaction, add 0.50 ml of the reaction mixture for the second reaction immediately to start the second reaction at 37°C.
4. After 20 min, add 1.5 ml of reaction dilution solution and mix.
5. After 10 min, measure the absorbance at 500 nm.
Absorbance sample: A_s
blank : A_b
 $\Delta A = (A_s - A_b) \leq 0.40$ Abs

■ Calculation

$$\text{Activity (U/mg)} = \frac{\Delta A/10}{12.2 \times 1/2} \times \frac{1}{2} \times \frac{3.00}{0.05} \times \frac{1}{X}$$

- 12.2 : millimolar extinction coefficient of quinoneimine dye at 500 nm (cm² / μ mole)
2 : the multiplier derived from the fact that 1 mole of phosphatidyl choline produces 2 mole of H₂O₂
1/2 : the multiplier derived from the fact that 2 mole of H₂O₂ produces 1 mole of quinoneimine dye
10 : reaction time (min)
3.00 : final volume (ml)
0.05 : volume of enzyme solution (ml)
X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition (Figure 3).

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