HDL Ultra Cholesterol Calibrator

CATALOGUE NUMBER: 6272-3 SIZE: 3 x 1 mL

INTENDED USE

For the calibration of the HDL Ultra Cholesterol assay and Ultra N-geneous® HDL Cholesterol assay.

TEST SUMMARY

Refer to the HDL Ultra Cholesterol Reagent and Ultra N-geneous® HDL Cholesterol package insert.

REAGENTS

The HDL Ultra Cholesterol Calibrator is a preparation of lyophilized human serum containing lipoproteins from the various lipoprotein classes including high-density lipoproteins.

Note: The value of the calibrator was assigned by procedures traceable to the CDC HDL cholesterol reference method.^{1,2}

WARNINGS AND PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

Rx ONLY

- Do not use beyond the expiration date.

Warning

Contains: Sodium Azide

H302 – Harmful if swallowed.

H412 – Harmful to aquatic life with long lasting effects.

Prevention – P264 – Wash thoroughly after handling.

Avoid release to the environment.

P270 – Do not eat, drink or smoke when using this product.

Response – P314 – Get medical advice/attention if you feel unwell.

If swallowed: Call a poison center/doctor if you feel unwell. Rinse mouth.

Disposal – P501 – Dispose of contents/container in accordance with local/regional/national/international regulations.

- For *In Vitro* Diagnostic Use.
- Do not pipette by mouth.
- Do not use the calibrator after the expiration date printed on the box.
- Warning:** Human source material. Treat as potentially infectious. Each plasma donor unit used in the preparation of this product has been tested by an FDA-approved method and found nonreactive for HBsAg, HCV, HIV 1 & 2 and HIV-1 Antigen. Because no known test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good

laboratory practices using appropriate precautions.

- Caution:** Contains sodium azide, which may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush drain with a large volume of water to prevent buildup. Dispose of in accordance with local, state and federal regulations.
- HDL Ultra Cholesterol Reagent and Ultra N-geneous® HDL Cholesterol Reagent must be used with HDL Ultra Cholesterol Calibrator.

See Safety Data Sheet for additional information.

PRODUCT PREPARATION, STORAGE AND STABILITY

Allow the vial of calibrator to equilibrate to room temperature before reconstitution. Lyophilized serum calibrator should be reconstituted by adding 1.0 mL of distilled or deionized water. Close the vial and let stand for 20 minutes. Dissolve the contents of the vial by swirling gently to avoid the formation of foam.

DO NOT SHAKE

Unopened calibrator is stable at 2-8°C until the expiration date printed on the vial.

After reconstitution, calibrator is stable for 14 days at 2-8°C. Reconstitution stability of the calibrator may be extended by aliquotting and freezing the reconstituted calibrator preparation at less than -70°C for up to 4 weeks.

PRODUCT DETERIORATION

Turbidity would indicate deterioration.

DISPOSAL

Product must be disposed of in accordance with all Federal, Provincial, State and local regulations.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Description	Configuration	Catalog Number
HDL Ultra Cholesterol Calibrator	3 x 1 mL	6272-3

MATERIALS REQUIRED (BUT NOT PROVIDED)

Description	Configuration	Catalog Number
HDL Ultra Cholesterol Reagent	R1 1 x 60 mL R2 1 x 20 mL	6121
HDL Ultra Cholesterol Reagent	R1 1 x 250 mL R2 1 x 80 mL	6122
Ultra N-geneous® HDL Cholesterol Reagent R1	R1 1 x 750 mL	80-6283-00
Ultra N-geneous® HDL Cholesterol Reagent R2	R2 1 x 250 mL	80-6277-00

1. Class A volumetric pipettes.
2. Distilled, deionized, Type II water or equivalent.

CALIBRATION

The calibrator should be used to analyze patient samples in accordance with the instructions outlined in Sekisui Diagnostics HDL Ultra Cholesterol Reagent and Ultra N-geneous® HDL Cholesterol Reagent package insert. The value of the HDL Ultra Cholesterol Calibrator, which can be found on the calibrator vial label, is assigned by procedures traceable to the CDC HDL cholesterol reference method.^{1,2} Refer to the instrument manufacturer's recommendation for calibration frequency. If assistance is required please call Sekisui Diagnostics Technical Services at 800-565-0265. Quality Control values should be within the expected range.

All trademarks, brands, product names and trade names are the property of their respective companies.

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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In vitro diagnostic medical device
Appareil médical de diagnostic *in vitro*.

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Caution, consult accompanying document
Prudence, consulter les documents joints

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For use by or on the order of a physician only (applicable to USA classification only)

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1. National Reference System for Cholesterol. CRMLN HDL Cholesterol Protocol, November 2002.
2. Kimberly MM, Leary ET, Cole TG, Waymack PW. Selection, validation, standardization, and performance of a designated comparison method for HDL cholesterol for use in the Cholesterol Reference Method Laboratory Network. *Clinical Chemistry* 1999; 45:1803-12.


Refer to the HDL Ultra Cholesterol Reagent and the Ultra N-geneous® HDL Cholesterol package insert for additional references.

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N-geneous® is a registered trademark of Sekisui Diagnostics, LLC. ©2011 Sekisui Diagnostics, LLC - All rights reserved.

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The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany
Phone: (+49)-511-6262 8630
Fax: (+49)-511-6262 8633

IN6272-3
March 27, 2019