# EN

## RPR CALIBRATOR SET

### INTENDED USE

For the calibration of the RPR assay. The RPR Calibrator Set can only be used with the Sekisui Diagnostics RPR Reagent.

### SUMMARY AND PRINCIPLE

See the RPR Reagent package insert.

### REAGENTS/COMPOSITION

The RPR calibrator Set is a set of 5 calibrators based on human serum and bovine serum albumin (BSA) with assigned concentrations of syphilitic anti-lipid antibodies.

#### Precuations and Warnings

1. For In Vitro Diagnostic Use.
2. Do not use the calibrators beyond the expiration date printed on the labels.
3. **Warning:** Human source material. Treat as potentially infectious. Each plasma donor unit used in the preparation of this product has been tested by an FDA- approved method and found non-reactive for HBsAg, anti-HCV, anti-HIV 1 & 2 and HIV-1 antigens. Because no known test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good laboratory practices using appropriate precautions.<sup>1</sup>
4. **Caution:** RPR Calibrators contain 0.1% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.
5. Disposal of all waste material should be in accordance with local guidelines.
6. **Warning:** This product contains syphilis positive serum and should be handled as a potentially biohazardous material in accordance with universal/ standard precautions.

#### Preparation of Calibrator

Calibrators: Liquid, ready to use.

Invert to mix before use. Avoid the formation of foam.

#### Storage and Stability

Unopened calibrator is stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened and capped, the calibrators are stable up to 4 weeks at 2 - 8°C.

DO NOT FREEZE.

#### Indications of Deterioration

Presence of turbidity or microbial growth may indicate deterioration.

### PROCEDURE

#### Materials Provided

Description	Configuration	Catalog Number
RPR Calibrator Set	5 levels x 1mL	486623

#### Materials Required but not Provided

Description	Configuration	Catalog Number
RPR Reagent 1	1 x 60mL	486616
RPR Reagent 2	1 x 20mL	
RPR Control Set	Positive 2 x 1mL Negative 2 x 1mL	486630

- Volumetric pipette or calibrated single channel pipette.
- Analyzer capable of running two-reagent chemistries.

#### Calibration

When using this calibrator set, treat it in exactly the same manner as a patient sample.

The RPR assay requires a multipoint calibration curve. The values of the RPR Calibrators, which can be found on the calibrator vial labels, are traceable to an in-house standard.

Refer to the instrument operator's manual for analyzer specific calibration procedures and for guidance in determining calibration frequency.

Quality Control values should be within the expected ranges.

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1. Richardson JH and Barkley WE, eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC: 1984.
2. Data on file at SEKISUI.

#### Definitions for Symbols/ Definition Der Symbole/ Définition des symboles/ Definizioni dei simboli/ Definiciones de los símbolos/ Definice symbolů/ Определение символов

##### REF

Catalog number  
Bestellnummer  
Référence  
Numero di catalogo  
Número de catálogo  
Katalogové číslo  
Номер по каталогу

Temperature limitation  
Temperaturbegrenzung  
Limites de température  
Limite di temperatura  
Límite temperature  
Teplotní omezení  
Температурное ограничение

##### IVD

For in vitro diagnostic use  
Für die In-Vitro-Diagnostik bestimmt  
Pour diagnostic in vitro  
Per uso diagnostico in vitro  
Para uso en diagnóstico in vitro  
K diagnostickému použití in vitro  
Только для in vitro диагностики

Manufactured by  
Hergestellt von  
Fabriqué par  
Prodotto da  
Fabricante  
Výrobce  
Изготовлено

Use by  
Zu verwenden bis  
Utiliser avant  
Usare entro il  
Fecha de caducidad  
Datum spotřeby  
Использовать до

Consult instructions for use  
Gebrauchsanweisung beachten  
Consulter le mode d'emploi  
Consultare le istruzioni per l'uso  
Consulte las instrucciones de uso  
Prostudujte si návod k použití  
Обратитесь к инструкции по применению

##### LOT

Batch code  
Chargenkennung  
Code de lot  
Codice del lotto  
Código del lote  
Kód šarže  
Код партии

Caution, consult accompanying documents  
Achtung, Begleitunterlagen beachten  
Attention, consulter la documentation jointe  
Attenzione, consultare i documenti di accompagnamento  
Precaución, consulte los documentos anexos  
Upozornění, prostudujte si příloženou dokumentaci  
Внимание! Обратитесь к сопровождающей документации

##### EC REP

Authorized Representative in the European Community  
Autorisierte Vertretung in der Europäischen Gemeinschaft  
Représentant agréé dans la Communauté européenne  
Rappresentante autorizzato nell'Unione Europea  
Representante autorizado en la Comunidad Europea  
Autorizovaný zástupce pro Evropské společenství  
Уполномоченный представитель в Европейском сообществе

Biological risks  
Biologisches Risiko  
Risques biologiques  
Rischi biologici  
Riesgos biológicos  
Biologická rizika  
Биологический риск

SEKISUI MEDICAL CO., LTD.  
1-3 , Nihonbashi 2 -chome, Chuo-ku  
Tokyo, 103-0027, Japan

**EC REP**

Medical Device Safety Service GmbH (MDSS)  
Schiffgraben 41  
30175 Hannover, Germany

**SEKURE™**  
CHEMISTRY

**CE**

April 2019 KI486623.03

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

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Das Wort SEKURE und das SEKURE-Logo sind Markenzeichen von Sekisui Diagnostics, LLC.
Le terme SEKURE et le logo Sekure sont des marques déposées de Sekisui Diagnostics, LLC.
La parola SEKURE e il logo Sekure sono marchi commerciali registrati di Sekisui Diagnostics, LLC.
La palabra SEKURE y el logotipo de SEKURE son marcas comerciales registradas de Sekisui Diagnostics, LLC.
Slovo SEKURE a logo Sekure jsou registrované ochranné známky společnosti Sekisui Diagnostics, LLC.