

Amylase Direct Reagent

FOR THE QUANTITATIVE MEASUREMENT OF α -AMYLASE

METHOD: DIRECT (CNPG₃); KINETIC

α -Amylase measurements are used to diagnose pancreatitis and other pancreatic diseases.

The SEKISUI Direct Amylase assay involves the use of a chromogenic substrate, 2-chloro-4-nitrophenol (CNPG₃). The methodology enables direct measurement of α -amylase without the use of auxiliary enzymes. It is intended for the quantitative measurement of α -amylase activity in serum and plasma.

Features

- One part stable liquid; ready to use reagent
- Short reaction time with no lag phase
- Precise
- No significant interferences from triglyceride, glucose, ascorbic acid, bilirubin or hemoglobin

Benefits

- Easy to use, no additional reagent preparation required
- Faster throughput on automated analyzers
- Confidence in results
- Reduces repeat testing and the need for sample dilutions

Performance Characteristics

Precision

- Within-Run: $\leq 1.4\%$
- Total Precision: $\leq 1.7\%$

Accuracy^(a)

- Slope: 0.9
- Intercept: -2.50 U/L
- Correlation Coefficient: 0.9985

Linearity

- Up to 2000 U/L

(a) The performance of this method (y) was compared with the performance of another α -amylase method (x) on a Roche/Hitachi® 911 analyzer.

No Significant Interferences Up to Levels Indicated

- Liposyn: 3000 mg/dL (1%)
- Triglyceride: 3000 mg/dL (33.9 mmol/L)
- Ascorbic Acid: 50 mg/dL (2840 μ mol/L)
- Unconjugated Bilirubin: 50 mg/dL (855 μ mol/L)
- Conjugated Bilirubin: 50 mg/dL (593 μ mol/L)
- Hemoglobin: 500 mg/dL (77.5 μ mol/L)
- Glucose: 2000 mg/dL (111 mmol/L)

Reference Range

- 20-104 U/L

Amylase Direct Reagent


Ordering Information

	Configuration	Catalog Number
AMYLASE DIRECT REAGENT	10 x 5 mL	80-5383-00
AMYLASE DIRECT REAGENT	1L	80-5451-00

To learn more about our clinical chemistry reagent products scan the QR code or go to [Clinical Chemistry Reagents – SEKISUI Diagnostics](#).



Sam



Skyler



Suzy

Securing accurate and reliable Clinical Chemistry results!



THE AMERICAS

SEKISUI Diagnostics, LLC
One Wall Street
Burlington, MA 01803
Phone: 800 332 1042

Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent, ME16 0LQ, UK
Phone: +44 1622 607800

Email: info@sekisui-dx.com

SEKISUI
DIAGNOSTICS

Because every result matters™

sekisuidiagnostics.com