

**ENZYMES**

# Reverse Transcriptase

**ORIGIN** *E.coli* (Recombinant)

**CAT# TRT Range**  
**EC# 2.7.7.7**

**SPECIFICATION**

**Concentration/Activity\*** 100 U/ $\mu$ L  
**Ribonuclease Activity** None detected

\*One unit is defined as the amount of enzyme required to incorporate 1 nmole of dTTP into an acid-insoluble material in 10 min at 42°C.

**PRODUCT FORMAT**

Formulated in 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 10 mM DTT, 0.01% Nonidet P-40, 50% Glycerol and store at -20°C.

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CONTENT/FORMAT</b>
TRT-101	ReverTra Ace (Reverse Transcriptase) Kit	0.1mL (10,000U) Reverse Transcriptase 1mL 5x RTase Buffer
TRT-179	ReverTra Ace (Reverse transcriptase) 500KU	5mL (500KU) Reverse Transcriptase
TRT-1B	5x RTase Buffer 1mL	1mL 5x Reverse Transcriptase Buffer

**DESCRIPTION AND APPLICATION**


ReverTra Ace is a high efficient M-MLV (Moloney Murine Leukemia Virus) reverse transcriptase that has been genetically modified to remove RNase H activity and increase reaction efficiency. It is the preferred enzyme for applications requiring full-length cDNAs and high product yields from total RNA, mRNA, rRNA, etc.

**CHARACTERISTICS**

**Features:**

- RNase minus M-MLV RTase with improved performance.
- Enables the synthesis of longer cDNAs ( $\geq 14$  kb) than the WT-enzyme.
- Exhibits excellent reaction efficiency at high temperatures.

**Application Data:** cDNAs were synthesized with oligo (dT)30 primers and 100 U enzyme/ poly (A)<sup>+</sup> RNA mixture (1.35-9.49 kb, 0.4 mg) as templates for 30 min at various temperatures. cDNAs were labeled with (<sup>32</sup>P-dCTP) during the reaction. The synthesized cDNAs were separated by 1% denatured agarose gel electrophoresis, and detected. The results suggested that our Reverse Transcriptase **can elongate efficiently at 42-55°C** compared to other RNase H minus RTases from other companies.



**Application Data:** G3PDH genes (500 bp) were amplified by PCR using cDNA templates that were synthesized with various RNase H minus RTases from G3PDH mRNA (10<sup>2</sup>-10<sup>5</sup> copies/ reaction). The RTase reaction was performed with specific reverse primers and 100 U enzyme at 42°C for 20 min. The results indicated that our Reverse Transcriptase is suitable for RT-PCR amplifications that require sensitivity.



**Application Data:** cDNA was synthesized by ReverTra Ace using a specific primer for the 3'-end of dystrophin mRNA at 42°C for 30 min. The 5' region at a distance of 14 kb from the 3' end of the dystrophin gene was amplified by PCR. The result indicated that ReverTra Ace can elongate cDNA of ≥14 kb.



**THE AMERICAS**

SEKISUI Diagnostics, LLC  
 One Wall Street,  
 Burlington, MA 01803  
 Phone: 800 332 1042  
 Fax: 800 762 6311

**INTERNATIONAL**

SEKISUI Diagnostics (UK) Limited  
 Liphook Way, Allington  
 Maidstone, Kent, ME16 0LQ, UK  
 Phone: +44 1622 607800  
 Fax: +44 1622 607801

MANUFACTURED BY:



DISTRIBUTED BY:



*Because every result matters™*