EN

PHOSPHORUS-SL ASSAY

CATALOGUE NUMBER: 117-30 **SIZE:** R1: 3 x 100 mL, R2: 1 x 75 mL
CAL: 1 x 15 mL

Note: Changes are highlighted.

INTENDED USE

For the *in vitro* quantitative measurement of inorganic phosphorus in serum.

TEST SUMMARY

Increased serum inorganic phosphorus concentrations (hyperphosphatemia) are usually the result of Vitamin D overdose, hypoparathyroidism, or renal failure. Decreased serum concentrations (hypophosphotemia) usually result from rickets, hyperparathyroidism, or Fanconi Syndrome (a defect in absorption of phosphorus and other metabolites from the glomerular filtrate).⁽¹⁾

The majority of methods for determining inorganic phosphorus are based on the reduction of a phosphorus-molybdate complex with a reducing agent resulting in the subsequent formation of molybdenum blue. In 1972, Daly and Ertingshausen introduced a method for inorganic phosphorus which measured the unreduced phosphorus-molybdate complex. This procedure is a modification of the Daly and Ertingshausen technique.⁽²⁾

TEST PRINCIPLE

Inorganic Phosphorus + Ammonium Molybdate $\xrightarrow{\text{H}_2\text{SO}_4}$ unreduced Phosphomolybdate complex

Inorganic phosphorus reacts with ammonium molybdate in the presence of sulfuric acid to produce an unreduced phosphomolybdate complex. The increase in absorbance at 340 nm is directly proportional to the concentration of inorganic phosphorus in the sample.

REAGENTS

Phosphorus Blank Reagent (R1): An acidic solution (pH < 1 at 25°C) containing a surfactant.

Phosphorus Molybdate Reagent (R2): An acidic solution (pH < 1 at 25°C) containing 2.8 mmol/L ammonium molybdate.

Phosphorus Calibrator: 1 x 15 mL of a solution containing 6.2 mg/dL (2 mmol/L) phosphorus. (For use with manual assays only.)

WARNINGS & PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

R_x ONLY

PHOSPHORUS-SL Blank Reagent (R1)

Warning

Contains: Sulphuric acid (CAS No) 7664-93-9

Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

Precautionary statements

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

PHOSPHORUS-SL Molybdate Reagent (R2)

Warning

Contains: Sulphuric acid (CAS No) 7664-93-9

Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

Precautionary statements

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

PHOSPHORUS-SL Calibrator

Warning

Contains: ethanediol (CAS No) 107-21-1

Hazard statements

H373 - May cause damage to organs through prolonged or repeated exposure

Precautionary statements

P260 - Do not breathe dust/fume/gas/mist/vapors/spray.

P314 - Get medical advice/attention if you feel unwell.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are provided in a ready to use format.

Supplied reagent is stable at 18-26°C until expiry date stated on the labels. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. Serum should be separated promptly from the red blood cells because erythrocytes contain phosphate concentrations several times greater than those found in sera.⁽³⁾

SAMPLE STORAGE

Phosphate is stable in serum that has been separated from the clot for 1 year at -20 °C, 4 days at 4-8 °C and 1 day at 20-25 °C.⁽⁴⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁵⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, and hemolysis were evaluated for this phosphorus method on a Roche/Hitachi® 704 analyzer using a significance criterion of > 10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
mg/dL	mmol/L			
3.2	1.04	Hemoglobin	400 mg/dL	62.0 µmol/L
3.4	1.10	Bilirubin	40 mg/dL	684 µmol/L
3.2	1.03	Intralipid	600 mg/dL	1800 mg/dL (20 mmol/L) Simulated Triglycerides

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenström's Macroglobulinemia may produce unreliable results.

Samples containing eltrombopag should not be used.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁶⁾

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' Phosphorus-SL reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:50 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800) 565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the inorganic phosphorus concentration of each sample.

TEST LIMITATIONS

A sample with an inorganic phosphorus concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹⁾

Adults 2.7-4.5 mg/dL (0.87-1.45 mmol/L)
Children 4.0-7.0 mg/dL (1.29-2.26 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

RESULTS

Inorganic phosphorus concentration is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁵⁾

The linearity of the procedure described is 20 mg/dL (6.46 mmol/L). The lower limit of detection of the procedure described is 0.1 mg/dL (0.03 mmol/L). This data results in a reportable range of 0.1 to 20 mg/dL (0.03 to 6.46 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁵⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
3.3	1.07	0.06	0.02	1.7	3.3	1.07	0.05	0.02	1.5
7.5	2.42	0.10	0.03	1.4	7.7	2.49	0.04	0.01	0.6

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁵⁾

The performance of this method (y) was compared with the performance of a similar phosphorus method (x) on a Roche/Hitachi® 704 analyzer. Forty-four patient serum samples ranging from 2.3-20.4 mg/dL (0.74-6.58 mmol/L) were tested and gave a correlation coefficient of 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 0.9926 (\text{reference method}) + 0.00 \text{ mg/dL (0.0 mmol/L)}.$$

All other trademarks are the property of their respective owners.

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REFERENCES / BIBLIOGRAFÍA

- Burtis, C.A., Ashwood, E.R., *Tietz Textbook of Clinical Chemistry, 2nd ed.* W.B. Saunders Co., Toronto, pp 1950-1952, p. 2202 (1994).
- Daly, J.A. and Ertingshausen, G., *Direct Method for Determining Inorganic Phosphorus in Serum with the Centrifichem*, Clin. Chem. 18:263, (1972).
- Kaplan, L., Pesce, A.J., *Clinical Chemistry, Theory, Analysis, Correlation, Third Edition*, Mosby, p. 552 (1996).
- World Health Organization. *Use of Anticoagulants in Diagnostic Laboratory Investigations*. Geneva: World Health Organization; 2002:39
- CLSI Method Evaluation Protocols, *Clinical and Laboratory Standards Institute, Wayne, PA*.
- Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington, 1990.

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

IN11710-20
November 29, 2022

The Americas
SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisui-dx.com
techservices@sekisui-dx.com

International
SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisui-dx.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com