EN

Calcium Assay

CATALOGUE NUMBER: 140-20 **SIZE:** 4 x 125 mL

Note: Changes are highlighted.

INTENDED USE

For the *in vitro* quantitative measurement of total calcium in serum.


TEST SUMMARY

The majority of calcium in the body is present in bones. The remainder of the calcium is present in serum and has various functions. For example, calcium ions decrease neuromuscular excitability, participate in blood coagulation, and activate some enzymes.

Hypercalcemia can result from hyperparathyroidism, hypervitaminosis D, multiple myeloma, and some neoplastic diseases of bone.⁽¹⁾ Hypocalcemia can result from hypoparathyroidism, steatorrhea, nephrosis, nephritis, and pancreatitis.⁽¹⁾

Calcium has traditionally been difficult to measure accurately and precisely, and a large variety of methods have been developed. Among these are flame photometry, oxalate precipitation with titration, atomic absorption spectrophotometry, EDTA chelation, and more recently calcium dye complexes which are measured spectrophotometrically. Examples of calcium dyes are o-cresolphthalein complexone and Arsenazo III, the latter being the dye used for calcium measurement in this method.

TEST PRINCIPLE

Arsenazo III [2,2'-(1,8-dihydroxy-3,6-disulfonaphthylene-2,7-bisazo) bisbenzenearsonic acid] reacts with calcium in an acid solution to form a blue-purple complex. The color developed has a maximum absorbance at 650 nm and is proportional to the calcium concentration in the sample.

REAGENTS

Calcium Reagent: A solution containing ~ 0.15 mmol/L arsenazo III, a buffer, and a surfactant.

Calcium Calibrator for Manual Assays: 1 x 15 mL of a solution containing 10.0 mg/dL (2.50 mmol/L) calcium and a preservative.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use.

Rx ONLY

Calcium Reagent

Contains: 4-nonylphenol, branched, ethoxylated (CAS No) 127087-87-0,

Formaldehyde (CAS No) 50-00-0

Hazard statements

H412 - Harmful to aquatic life with long lasting effects

Precautionary statements

P273 - Avoid release to the environment.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use. Supplied reagent is stable at 18-26°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, State and local regulations.

SPECIMEN

Fresh clear serum.

SAMPLE STORAGE

Serum calcium is stable for one day at 4°C and for one year when frozen.

GLASSWARE PREPARATION

All glassware used in calcium assays should be decontaminated by soaking overnight in 1 N HCl or for 6 hours in 6 N HCl. The glassware should be well rinsed in deionized water before use. The syringes on dispensers should also be rinsed in 6 N HCl followed by deionized water before being used for a calcium assay.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽²⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown. Environmental contamination of calcium assays causes inaccurate results unless the previously mentioned decontamination procedures are used.

Interferences from icterus, lipemia and hemoglobin were evaluated for this Calcium method on a Roche/Hitachi® 911 analyzer using a significance criterion of > 10% variance from control.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
9.6 mg/dL	2.4 mmol/L	Hemoglobin	1000 mg/dL	155 µmol/L
10.1 mg/dL	2.5 mmol/L	Bilirubin	40 mg/dL	684 µmol/L
9.3 mg/dL	2.3 mmol/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽³⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' Calcium reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:100, and a wavelength reading of 650 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the calcium concentration of each sample.

TEST LIMITATIONS

A sample with a calcium concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS

Serum:⁽⁴⁾

Premature Infants	6.0-10.0 mg/dL	1.5-2.5 mmol/L
Full-term Infants	7.3-12.0 mg/dL	1.8-3.0 mmol/L
1 to 2 years	10.0-12.0 mg/dL	2.5-3.0 mmol/L
Adults	8.0-10.5 mg/dL	2.0-2.6 mmol/L

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 911 automated analyzer unless otherwise stated.

RESULTS

Calcium concentration is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽²⁾

The linearity of the procedure described is 15.0 mg/dL (3.75 mmol/L) in serum. The lower limit of detection of the procedure described is 0.04 mg/dL (0.01 mmol/L). This data results in a reportable range of 0.04 to 15.00 mg/dL (0.01 to 3.75 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽²⁾

Precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV (%)	Within Run SD		Within Run CV (%)
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	
9.5	2.4	0.09	0.02	1.0	0.08	0.02	0.8
11.9	3.0	0.11	0.03	0.9	0.09	0.02	0.7

ACCURACY (CLSI EP9)⁽²⁾

The performance of this method (y) was compared with the performance of another calcium method (x). One hundred and eighteen samples ranging from 6.5 to 13.5 mg/dL (1.6 to 3.4 mmol/L) were tested and gave a correlation coefficient of 0.9889. Linear regression analysis gave the following equation:

$$\text{This method} = 0.875 (\text{reference method}) + 0.91 \text{ mg/dL} (0.23 \text{ mmol/L}).$$

All other trademarks are the property of their respective owners.

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Manufacturer
Fabricant
Fabricante

Consult instructions for use
Consulter les directives d'utilisation
Consulte las instrucciones de uso

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1. Tietz, N.W. (Ed.), *Fundamentals of Clinical Chemistry*, W.B. Saunders Co., Toronto, 636-638, 937 (1970).
2. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
3. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition, 1990.
4. Kaplan, Lawrence A. and Pesce, Amadeo J., (Eds.), *Clinical Chemistry: Theory, Analysis, Correlation*, Mosby-Year Book Ins., St. Louis (1996).

IN14020-26
November 2, 2023

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com