EN

MAGNESIUM REAGENT

CATALOGUE NUMBER: 175-12 SIZE: 4 x 125 mL

Note: Changes are highlighted

INTENDED USE

The SEKISUI Magnesium Assay is intended for the quantitative measurement of magnesium in human serum and plasma (Lithium Heparin) on automated clinical chemistry analyzers. Magnesium measurements are used in the diagnosis and treatment of hypomagnesemia and hypermagnesemia. This device is intended for professional use and *in vitro* diagnostic use only.

TEST SUMMARY

Magnesium measurements are clinically important. Magnesium, along with potassium, is the most abundant intracellular cation. Magnesium is a required coenzyme in the metabolism of carbohydrates, lipids, and proteins.

Decreased magnesium concentrations, hypomagnesemia, may cause tremors, muscle irritability, increased blood pressure and heart rate.

Elevated magnesium concentrations, hypermagnesemia, may lead to impaired respiration, coma, and cardiac arrest if left uncorrected.⁽¹⁾

The method of choice for magnesium measurement is atomic absorption spectrophotometry.⁽²⁾ This method is not practical for clinical laboratories so dye binding methods using methylthymol blue, titan yellow, calmagite, and xylydyl blue-1 are used instead. This method employs xylydyl blue-1 which allows the use of a single reagent which is rapid and precise.

TEST PRINCIPLE

Xylydyl blue-1 + Mg⁺⁺ → Mg-xylydyl blue complex (red)

Under the alkaline conditions of the Magnesium reagent, magnesium present in the sample, when mixed with the reagent, forms a red complex with the diazonium salt of xylydyl blue. The concentration of magnesium in the sample is directly proportional to the amount of xylydyl blue – magnesium complex formed and can be measured spectrophotometrically by the decrease in absorbance measured at 660 nm.

REAGENTS

Magnesium Reagent: A solution containing buffer (pH 11.35 at 25 °C), ~0.14 mmol/L xylydyl blue-1, 0.1 mmol/L EGTA, and a surfactant.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

R₀ ONLY

MAGNESIUM REAGENT

Danger

Contains: Ethanolamine CAS 141-43-5, Boric Acid CAS 10043-35-3

Hazard statements

H314 - Causes severe skin burns and eye damage
 H318 - Causes serious eye damage
 H335 - May cause respiratory irritation
 H360 - May damage fertility or the unborn child

Precautionary statements

P201 - Obtain special instructions before use.
 P202 - Do not handle until all safety precautions have been read and understood.
 P260 - Do not breathe dust/fume/gas/mist/vapors/spray.
 P261 - Avoid breathing dust/fume/gas/mist/vapors/spray.
 P264 - Wash hands, forearms and face thoroughly after handling.
 P271 - Use only outdoors or in a well-ventilated area.
 P280 - Wear protective gloves/protective clothing/eye protection/face protection.
 P301+P330+P331 - If swallowed: rinse mouth. Do NOT induce vomiting.
 P303+P361+P353 - If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
 P304+P340 - If inhaled: Remove person to fresh air and keep comfortable for breathing.
 P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P308+P313 - If exposed or concerned: Get medical advice/attention.

P310 - Immediately call a poison center or doctor.
 P312 - Call a poison center or doctor if you feel unwell.
 P321 - Specific treatment (see supplemental first aid instruction on this label).
 P363 - Wash contaminated clothing before reuse.
 P403+P233 - Store in a well-ventilated place. Keep container tightly closed.
 P405 - Store locked up.
 P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

Supplied reagent is stable at 2-8°C until the expiry date stated on the label. Stability claims are based on accelerated stability studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum and plasma. Lithium Heparin specimens were tested and found to be acceptable. Use only suitable tubes for collection and storage. Various manufacturers of sample collection systems may contain differing materials that may affect results. Follow all instructions provided by the collection system manufacturer.

SAMPLE STORAGE⁽³⁾

Magnesium is stable in serum for 7 days at 4-8°C and for 12 months at -20°C.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁴⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

EGTA prevents calcium interference by preferential chelation of the calcium.

Interferences from hemolysis, icterus, ascorbic acid and lipemia were evaluated for this magnesium method on a Roche/Hitachi® analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum. Plasma is expected to be similar.

Concentration of Analyte		Substance Tested	Concentration of Interferent where Interference is Insignificant	
mg/dL	mmol/L			
2.1	0.86	Hemoglobin	1000 mg/dL	155 µmol/L
4.2	1.73		1000 mg/dL	155 µmol/L
7.9	3.25		1000 mg/dL	155 µmol/L
2.0	0.82	Conjugated Bilirubin	40 mg/dL	474 µmol/L
4.3	1.77		40 mg/dL	474 µmol/L
6.8	2.79		40 mg/dL	474 µmol/L
2.0	0.82	Unconjugated Bilirubin	40 mg/dL	684 µmol/L
4.2	1.73		40 mg/dL	684 µmol/L
6.8	2.79		40 mg/dL	684 µmol/L
2.3	0.95	Ascorbic Acid	3000 µg/dL	170 µmol/L
4.4	1.81		3000 µg/dL	170 µmol/L
6.8	2.79		3000 µg/dL	170 µmol/L
1.9	0.78	Intralipid	800 mg/dL	2400 mg/dL (27.1 mmol/L) Simulated Triglycerides
4.1	1.69		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides
6.9	2.84		1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁵⁾

ANALYTICAL PROCEDURE MATERIALS PROVIDED

SEKISUI Diagnostics' Magnesium reagent.
MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Calibration material.
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:117, and wavelength readings of 660/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800) 565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the magnesium concentration of each sample.

TEST LIMITATIONS

A sample with a magnesium concentration exceeding the linearity limit should be diluted 1:1 with 0.9% saline and reassayed, incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹⁾

1.6-2.6 mg/dL (0.66-1.07 mmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

RESULTS

Magnesium concentration is reported as mg/dL (mmol/L).

To convert magnesium results to mEq/L, use the following conversion factor:⁽⁶⁾
 $\text{mmol/L} \times 2 = \text{mEq/L}$ ($\text{mEq/L} \times 0.5 = \text{mmol/L}$)

REPORTABLE RANGE (CLSI EP6)⁽⁴⁾

The linearity of the procedure described is 8.0 mg/dL (3.29 mmol/L). The limit of blank and limit of detection of the procedure described are 0.2 mg/dL (0.08 mmol/L) and 0.3 mg/dL (0.12 mmol/L), respectively. The limit of quantitation of the procedure described is 0.3 mg/dL (0.12 mmol/L). This data results in a reportable range of 0.3 mg/dL (0.12 mmol/L) to 8.0 mg/dL (3.29 mmol/L). The reportable range will depend on the sample to reagent ratio used.

PRECISION STUDIES (CLSI EP5)⁽⁴⁾

Total precision data was collected on three concentrations of serum in 40 runs conducted over 20 days.

Concentration		Total SD			Total CV %	Concentration		Within Run SD		With-in Run CV %
mg/dL	mmol/L	mg/dL	mmol/L	mg/dL		mmol/L	mg/dL	mmol/L		
0.8	0.33	0.05	0.021	6.1	0.9	0.37	0.06	0.024	7.0	
2.0	0.82	0.07	0.029	3.7	2.0	0.82	0.04	0.016	2.2	
4.7	1.93	0.15	0.062	3.2	4.9	2.01	0.08	0.033	1.6	

Within run precision data was collected on three concentrations of serum each run 20 times in a single assay.

ACCURACY (CLSI EP9)⁽⁴⁾

The performance of this method (y) was compared with the performance of a similar xylidyl blue method (x) on a Roche/Hitachi® 717 analyzer. One hundred patient serum samples ranging from 0.5 mg/dL to 7.3 mg/dL (0.21-3.00 mmol/L) gave a correlation coefficient of 0.985. Linear regression analysis gave the following equation:

$$\text{This method} = 0.981 (\text{reference method}) - 0.0 \text{ mg/dL} (0.02 \text{ mmol/L}).$$

The performance of this method (plasma) was compared with the performance of this method (serum) on a Roche/Hitachi® 717 analyzer. Fifty-four (54) patient serum and plasma samples ranging from 0.5 mg/dL to 5.4 mg/dL (0.21-2.22 mmol/L) gave a correlation coefficient of 0.9903. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 1.009 (\text{serum}) - 0.1 \text{ mg/dL} (0.04 \text{ mmol/L}).$$

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

All other trademarks are the property of their respective owners.

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Manufacturer
Fabricant
Fabricante

Consult instructions for use
Consulter les directives d'utilisation
Consulte las instrucciones de uso

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1. Tietz, N.W., (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1983) p.338, 339.
2. Ratge, D., Kohse, K.P., and Wisser, H.; Measurement of Magnesium in Serum and Urine with a Random Access Analyzer by use of a Modified Xylidyl Blue-1 Procedure, Clinical Chimica Acta, 159 (1986) 197-203.
3. World Health Organization; Use of Anticoagulants in Diagnostic Laboratory Investigations and Stability of Blood, Plasma and Serum Samples. WHO/Lab/99.1 Rev. 2, 2002.
4. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition, 1990.
6. Tietz, N.W. (Editor) Clinical Guide to Laboratory Tests, W.B. Saunders Company, Philadelphia (1995) p. 418.

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

IN17512-8
July 16, 2024

The Americas
SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International
SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com