EN

Total Protein

CATALOGUE NUMBER: 200-55

SIZE: 4 x 125 mL

Note: Changes are highlighted.

INTENDED USE

For the *in vitro* quantitative measurement of total protein in serum.

TEST SUMMARY

The use of the Biuret reaction as a method for the estimation of protein in plasma was first introduced by Reigler.⁽¹⁾ Gornall, et. al.⁽²⁾ modified the procedure by adding sodium potassium tartrate which acts as a complexing agent to form a stable copper protein complex. This method uses the biuret reaction in which protein reacts with the reagent at an alkaline pH to form a blue-violet colored complex.

Total protein measurements are used in the diagnosis and treatment of a variety of diseases involving the liver, kidney, or bone marrow as well as other metabolic or nutritional disorders.

TEST PRINCIPLE

Protein + Cu⁺⁺ → blue-violet complex

At an alkaline pH, the protein reacts with the copper in the biuret reagent causing an increase in absorbance.

The increase in absorbance at 540 nm due to the formation of the colored complex is directly proportional to the concentration of protein in the reaction.

REAGENTS

Total Protein Reagent: A solution containing 31.9 mmol/L sodium potassium tartrate, 12.0 mmol/L copper sulfate pentahydrate, 30.1 mmol/L potassium iodide and 0.20 mol/L sodium hydroxide.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

R_oONLY

TOTAL PROTEIN REAGENT

Warning

Contains: sodium hydroxide; caustic soda (CAS No) 1310-73-2

Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

H402 - Harmful to aquatic life

H412 - Harmful to aquatic life with long lasting effects

Precautionary statements

P264 - Wash hands, forearms and face thoroughly after handling.

P273 - Avoid release to the environment.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid ingestion.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use. Supplied reagent is stable at 18 to 26°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE

Samples should be stored at 2 to 8°C.

ANALYTICAL SPECIFICITY (CLSI EP7)⁽³⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemoglobin were evaluated for this total protein method on an Olympus AU640[®] analyzer using a significance criterion of >10% from control. Hemolysed specimens should be avoided as hemoglobin falsely elevates total protein values.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
6.1 g/dL	61 g/L	Hemoglobin	200 mg/dL	31 µmol/L
6.8 g/dL	68 g/L	Bilirubin	40 mg/dL	684 µmol/L
6.6 g/dL	66 g/L	Intralipid	100 mg/dL	300 mg/dL (3.39 mmol/L) Simulated Triglycerides

Samples containing the following should not be used: Hydroxocobalamin.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁴⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' Total Protein reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

TEST CONDITION

For data presented in this insert, studies using this reagent were performed using an endpoint test mode, with a sample to reagent ratio of 1:50 and a wavelength reading of 540 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the total protein concentration of each sample.

TEST LIMITATIONS

A sample with a total protein concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁵⁾

6.0-8.2 g/dL (60-82 g/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected using automated procedures unless otherwise stated.

RESULTS

Total protein concentration is reported as g/dL (g/L).

REPORTABLE RANGE

The procedure is linear to 10.0 g/dL (100 g/L).

PRECISION STUDIES

Total precision data was collected on two control sera using a single lot of reagent in two runs conducted over 10 days.

Total Concentration		Total SD		Total CV (%)	Within Run Concentration		Within Run SD		Within Run CV (%)
g/dL	g/L	g/dL	g/L		g/dL	g/L	g/dL	g/L	
4.8	48	0.08	0.8	1.7	4.8	48	0.03	0.3	0.6
6.3	63	0.10	1.0	1.6	6.3	63	0.08	0.8	1.2

Within run precision data was collected by assaying two control sera 24 times.

ACCURACY

The performance of this method (y) was compared with the performance of another total protein method (x). Fifty-eight patient serum samples ranging from 4.6 to 7.7 g/dL (46 to 77 g/L) gave a correlation coefficient of 0.947. Linear regression analysis gave the following equation:

This method = 1.037 (reference method) - 0.36 g/dL (3.6 g/L)

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

TRADEMARKS

All other trademarks are the property of their respective owners.

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Manufacturer
Fabricant
Fabricante

Consult instructions for use
Consulter les directives d'utilisation
Consulte las instrucciones de uso

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1. Riegler, E., A Colorimetric Method for the Determination of Plasma Albumin, *Z. Anal. Chem.* 53, 242 (1914).
2. Gornall, A.G., Bardawill, C.J., David, M.M., Determination of Serum Proteins by Means of the Biuret Reaction, *J. Biol. Chem.* 177, 751 (1949).
3. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
4. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).
5. Kaplan, A., Szabo, L.L., *Clinical Chemistry: Interpretation and Techniques*, 2nd Ed., (1983) Lea and Febiger, Philadelphia, p. 405.

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

IN20070-13
June 20, 2024

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com