EN

TRIGLYCERIDE-SL ASSAY

CATALOGUE NUMBER: 236-60 SIZE: 2 x 100 mL
236-99 1 x 1000 mL

Note: Changes are highlighted.

INTENDED USE





For the *in vitro* quantitative measurement of triglyceride in serum.

TEST SUMMARY

Triglycerides are synthesized in intestinal epithelial cells by the esterification of glycerol and fatty acids. Once formed, they combine with other lipids and lipoproteins in complex assemblages known as chylomicrons, which are released from the intestinal epithelium into the circulatory system. Elevated serum triglyceride concentrations are associated with increased risk of coronary artery disease, pancreatitis, and numerous genetic lipoprotein disorders⁽¹⁾.

Both chemical and enzymatic methods have been described for serum triglyceride measurement. The chemical methods are complex, requiring extraction and purification of the sample prior to quantitation⁽²⁾. A fully enzymatic, single reagent procedure has been described by Fossati and Prencipe⁽³⁾. The method employs glycerol phosphate oxidase in a modified Trinder reaction, and is both rapid and sensitive. This triglyceride assay is a modification of that method.

TEST PRINCIPLE

Serum triglycerides are hydrolyzed to glycerol and free fatty acids by lipase. In the presence of ATP and glycerol kinase (GK), the glycerol is phosphorylated to glycerol-1-phosphate. Glycerol-1-phosphate is then oxidized by glycerol phosphate oxidase (GPO) to yield hydrogen peroxide (H₂O₂). The hydrogen peroxide causes oxidative coupling of p-chlorophenol and 4-aminoantipyrine, producing a red colored quinoneimine dye complex. The increase in absorbance at 520 nm due to the formation of the quinoneimine dye is directly proportional to the concentration of triglycerides in the sample.

REAGENTS

Triglyceride Reagent: A buffered solution containing 0.4 mmol/L 4-aminoantipyrine, 2.3 mmol/L adenosine triphosphate, 3.0 mmol/L p-chlorophenol, >2400 U/L glycerol phosphate oxidase (microbial), >1000 U/L lipoprotein lipase (microbial), >540 U/L peroxidase (botanical), >400 U/L glycerol kinase (microbial), stabilizers, and preservatives.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

R₁ ONLY

Avoid contact with skin and eyes. See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use. The reagent included is stable until the expiry date stated on the label at 2-8°C. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. Specimens should be collected following a 12 hour fast and the serum should be separated from the clot as soon as possible. The patient should be seated for at least 5 minutes before collection.

SAMPLE STORAGE

Specimens may be stored for up to seven days at 4-8°C or for several weeks at -20°C to 0°C⁽⁴⁾. Repeated freezing and thawing of specimens should be avoided. Specimens should not be stored at room temperature for prolonged periods.

GLASSWARE PREPARATION

All glassware, cuvettes, and tubing should be thoroughly rinsed and glycerol free. Avoid the use of tubes with glycerol coated stoppers. Glycerol contamination will falsely elevate results.⁽⁵⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Grossly lipemic samples have been found to adversely affect the reaction for GPO triglyceride methods. Sera with extremely elevated triglycerides may falsely indicate a normal result⁽⁷⁾. It is recommended to dilute grossly lipemic sera 1 part serum to 4 parts saline prior to testing for triglyceride; this manual dilution will require the final result to be multiplied by a factor of 5.

Interferences from icterus and hemoglobin were evaluated for this method on a Roche/Hitachi 704 analyzer using a significance criterion of > 10% variance from control data was collected in serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
123 mg/dL	1.39 mmol/L	Hemoglobin	500 mg/dL	78 µmol/L
127 mg/dL	1.43 mmol/L	Bilirubin	20 mg/dL	342 µmol/L

Sample contamination with detergents and elevated free glycerol concentrations will cause positive interference in this assay.⁽⁵⁾

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

Samples containing the following should not be used: N-acetylcysteine (NAC), Methyl-4-aminoantipyrine (MAA).

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' Triglyceride-SL reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:75 and a wavelength reading of (primary/secondary) 505/660 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the triglyceride concentration of each sample.

TEST LIMITATIONS

A sample with a triglyceride concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

This assay is not a blanked assay. For the majority of fresh, refrigerated samples, the contribution to the triglyceride measurement caused by free glycerol is relatively low and constant and is of minimal clinical importance and does not generally exceed 10 mg/dL (0.11 mmol/L).⁽²⁾

REFERENCE INTERVALS

Expected values for serum triglycerides vary with age, sex, diet, ethnic origin, and geographic location. The upper limit of reference intervals for men and women have been reported as follows:⁽²⁾

Age (Years)	Males		Females	
	mg/dL	mmol/L	mg/dL	mmol/L
20-29	44-201	0.50-2.27	37-144	0.41-1.63
30-39	51-321	0.56-3.62	39-176	0.44-1.99
40-49	55-327	0.62-3.70	45-214	0.51-2.42
50-59	58-286	0.62-3.23	52-262	0.59-2.62

ATP III Classification of Serum Triglycerides (mg/dL)⁽⁹⁾

<150	Normal
150-199	Borderline High
200-499	High
≥500	Very High

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on an automated analyzer unless otherwise stated.

RESULTS

Triglyceride concentration/activity is reported as mg/dL (mmol/L).

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

The linearity of the procedure described is 1000 mg/dL (11.3 mmol/L). The lower limit of detection of the procedure described is 3.0 mg/dL (0.03 mmol/L). This data results in a reportable range of 3.0 – 1000 mg/dL (0.03-11.3 mmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Within Run SD		Within Run CV %
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	
56.5	0.64	0.78	0.009	1.4	0.32	0.004	0.6
111.9	1.26	1.46	0.020	1.3	0.69	0.008	0.6

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi 704. Fifty patient serum samples ranging from 48 mg/dL (0.54 mmol/L) to 369 mg/dL (4.17 mmol/L) gave a correlation coefficient of 0.9932. Linear regression analysis gave the following equation:

$$\text{This method} = 0.9946 (\text{reference method}) - 0.8 \text{ mg/dL (0.01 mmol/L)}$$

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REFERENCES / BIBLIOGRAFÍA

1. Fredrickson, D.S., et.al., New Eng J. Med. 276:32, 1976.
2. Burtis, C., and Ashwood, E.R., Editors, Tietz Textbook of Clinical Chemistry, 3rd edition, W.B. Saunders Co., Toronto, 1999.
3. Fossati, P. and Prencipe, L., Clin. Chem. 28:2077, 1982.
4. Heil, W, Koberstein, R. and Zawta, B., Reference Ranges for Adults and Children, Roche Diagnostics GmbH, 2004.
5. Sampson, M., Ruddel, M., and Elin, R.J., Effects of Specimen Turbidity and Glycerol Concentration on Nine Enzymatic Methods for Triglyceride Determination, Clin. Chem. 40:12, 221-228, 1994.
6. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
7. Shephard MDS, Whiting MJ Clin. Chem. 36(2):325-329, 1990.
8. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Third Edition, Washington, 1990.
9. ATP III Guidelines, National Cholesterol Education Program, National Institutes of Health, 2001.

IN23660-24
July 17, 2024

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com