EN

## ETHANOL L3K<sup>®</sup> ASSAY

CATALOGUE NUMBER: 273-30      SIZE: 6 x 30 mL

**Note:** Changes are highlighted

### INTENDED USE


For the *in vitro* quantitative measurement of ethanol concentration in serum, plasma, and urine.

### TEST SUMMARY

Ethanol (ethyl alcohol) is the most common toxic substance involved in medical-legal cases. Ethanol consumption is often a factor in all types of accidents and ethanol poisoning can be fatal. Ethanol testing of comatose patients aids in the differential diagnosis of the cause of the coma.<sup>(1)</sup>

Early methods for ethanol analysis depended on separation of the ethanol from the specimen. The method described here is an enzymatic procedure which does not require separation and is rapid and easily automated. This method is based on the first enzymatic method described by Bonnichsen and Theorelle<sup>(2)</sup> and later modified by other authors.<sup>(3)</sup> However, the current method is unique in that an aldehyde trap is no longer necessary in order to drive the reaction. This is accomplished by the use of a co-factor with a high oxidizing potential.

### TEST PRINCIPLE

The enzyme alcohol dehydrogenase (ADH) catalyzes the conversion of ethanol to acetaldehyde and NADH analog. The high oxidizing potential of the NAD<sup>+</sup> analog drives the reaction forward. The rate of increase in absorbance at 340 nm or 380 nm due to the reduction of NAD analog to NADH analog is proportional to the amount of ethanol present in the sample.

### REAGENTS

Ethanol L3K<sup>®</sup> Reagent: A solution containing a buffer (pH 7.3 at 25°C), 1.35 mmol/L NAD analog, > 203 KU/L alcohol dehydrogenase (yeast), stabilizers, and a preservative.

Ethanol Calibrator: 1 x 5 mL of a solution containing 92 mg/dL (20 mmol/L) ethanol and a preservative.

### WARNINGS & PRECAUTIONS FOR USE

For *in vitro* diagnostic use

**R<sub>o</sub>ONLY**

### ETHANOL L3K<sup>®</sup> CALIBRATOR

Warning

Contains: reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (CAS No) 55965-84-9

#### Hazard statements

H315 – Causes skin irritation  
H317 – May cause an allergic skin reaction  
H319 – Causes serious eye irritation  
H411 – Toxic to aquatic life with long lasting effects

#### Precautionary statements

P261 – Avoid breathing dust/fume/gas/mist/vapors/spray.  
P264 – Wash hands, forearms and face thoroughly after handling.  
P272 – Contaminated work clothing must not be allowed out of the workplace.  
P273 – Avoid release to the environment.  
P280 – Wear protective gloves/protective clothing/eye protection/face protection.  
P302+P352 – If on skin: Wash with plenty of water.  
P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P321 – Specific treatment (see supplemental first aid instruction on this label).  
P332+P313 – If skin irritation occurs: Get medical advice/attention.  
P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.  
P337+P313 – If eye irritation persists: Get medical advice/attention.  
P362+P364 – Take off contaminated clothing and wash it before reuse.  
P363 – Wash contaminated clothing before reuse.  
P391 – Collect spillage.

P501 – Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready to use.

Supplied reagent is stable at 2-8°C until expiry date.

Stability claims are based on real time studies.

The ethanol calibrator should be tightly sealed when not in use as exposure to the atmosphere causes evaporation of the ethanol.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with the Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum, plasma, or urine. Specimens should be kept sealed. The anticoagulants EDTA, heparin, and citrate have been tested and may be used with this assay.

Ethanol should not be used to swab the venipuncture site or to clean or sterilize glassware or other equipment used to collect the specimen or perform the assay. Avoid the use of test tubes, cuvettes, or pipettes which have been cleaned or sterilized with ethanol.

### SAMPLE STORAGE

Samples should be stored at 2-8°C for up to 3 days. If necessary, they may be frozen but repeated freeze-thaw cycles should be avoided.

Samples should be at room temperature (18-26°C) for testing.

Urine samples should be stored with as little dead air space as possible. The sample container should be gently inverted several times before opening.

Always tightly cap sample tubes to prevent evaporation of ethanol.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(4)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences for LDH (sample containing a specific concentration of lactate), icterus, lipemia (Intralipid), and hemoglobin were evaluated for this ethanol method on Beckman CX<sup>®</sup> analyzer using a significance criterion of >10% from the control. Interference data was collected in serum. Plasma data is similar.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
100 mg/dL	22 mmol/L	Bilirubin	40 mg/dL	684 µmol/L
100 mg/dL	22 mmol/L	Hemoglobin	600 mg/dL	93 µmol/L
100 mg/dL	22 mmol/L	LDH (in the presence of 200 mg/dL lactate)	10,000 U/L	10,000 U/L
100 mg/dL	22 mmol/L	Intralipid	1000 mg/dL	3000 mg/dL (34 mmol/L) Simulated Triglyceride

Solutions of potentially interfering compounds were prepared and analyzed. The concentrations of the solvents and the assay results, expressed as units of ethanol, are listed below. The results indicate the assay is designed to detect ethanol and not other alcohols.

Solvent	Solvent Concentration (mg/dL)	Units of Ethanol		% Reactivity
		mg/dL	mmol/L	
Acetone	100	1.3	0.3	1.3
	2000	1.2	0.3	0.1
Butanol	100	5.3	1.2	5.3
	200	12.4	2.7	6.2
Ethylene Glycol	100	2.8	0.6	2.8
	2000	4.7	1.0	0.2
Methanol	100	1.0	0.2	1.0
	2000	0.4	0.1	0.0
Isopropanol	100	0.3	0.1	0.3
	2000	6.4	1.4	0.3

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(5)</sup>

## ANALYTICAL PROCEDURE

### MATERIALS PROVIDED

SEKISUI Diagnostics' Ethanol L3K<sup>®</sup> reagent and calibrator.

### MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- 2) Quality Control materials.

### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a rate test mode, with a sample to reagent ratio of 1:100 and a primary wavelength reading of 340 nm or 380 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

### CALIBRATION

An ethanol calibrator is included with the reagents and is used as directed to calibrate the procedure. The ethanol calibrator is aqueous based and values were assigned using NIST traceable ethanol standards. The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

### CALCULATIONS

The analyzer automatically calculates the ethanol concentration of each sample.

### TEST LIMITATIONS

A sample with an ethanol concentration exceeding the linearity limit should be diluted with 0.9% saline and reanalyzed incorporating the dilution factor in the calculation of the value.

### REFERENCE INTERVALS

The legal plasma or serum limit for drivers is less than 80 mg/dL (17.4 mmol/L) (Canada)<sup>(6)</sup>.

### PERFORMANCE CHARACTERISTICS

Data presented was collected on a Beckman CX<sup>®</sup> analyzer unless otherwise stated.

### RESULTS

Ethanol concentration is expressed as mg/dL (mmol/L).

### REPORTABLE RANGE (CLSI EP6)<sup>(4)</sup>

The linearity of the procedure is 600 mg/dL (130.3 mmol/L). The lower limit of detection of the procedure described is 10 mg/dL (2.2 mmol/L). This data results in a reportable range of 10-600 mg/dL (2.2 to 130.3 mmol/L).

### PRECISION STUDIES (CLSI EP5)<sup>(4)</sup>

Total precision was collected on three concentrations of control materials in forty runs conducted over twenty days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
mg/dL	mmol/L	mg/dL	mmol/L		mg/dL	mmol/L	mg/dL	mmol/L	
40.3	8.8	1.8	0.4	4.6	41.5	9.0	1.4	0.3	3.4
98.1	21.3	2.6	0.6	2.6	97.6	21.2	1.4	0.3	1.5
245.1	53.2	6.1	1.3	2.5	240.0	52.1	2.4	0.5	1.0

Within run precision data was collected on three concentrations of control materials; each control material was run twenty times in a single assay.

### ACCURACY (CLSI EP9)<sup>(4)</sup>

The performance of this method (y) was compared with the performance of another rate ethanol method (x) on a Beckman CX<sup>®</sup>. Forty patient serum samples ranging from 22.5 to 651.9 mg/dL (4.9 to 141.5 mmol/L) gave a correlation coefficient of 0.9994. Linear regression analysis gave the following equation:

$$\text{This method} = 1.06 (\text{reference method}) - 2.45 \text{ mg/dL} (0.5 \text{ mmol/L}).$$

The performance of this method (y) was compared with the performance of another rate ethanol method (x) on a Beckman CX<sup>®</sup>. Forty patient plasma samples ranging from 17.3 to 580.8 mg/dL (3.8 to 125.9 mmol/L) gave a correlation coefficient of 0.9994. Linear regression analysis gave the following equation:

$$\text{This method} = 1.03 (\text{reference method}) - 5.86 \text{ mg/dL} (1.3 \text{ mmol/L}).$$

The performance of this method (y) was compared with the performance of another rate ethanol method (x) on a Beckman CX<sup>®</sup>. Forty patient urine samples ranging from 17.6 to 583.4 mg/dL (3.8 to 126.6 mmol/L) gave a correlation coefficient of 0.9998. Linear regression analysis gave the following equation:

$$\text{This method} = 1.01 (\text{reference method}) + 0.59 \text{ mg/dL} (0.13 \text{ mmol/L}).$$

The information presented above is based on results from SEKISUI Diagnostics P.E.I. Inc. studies and is current at the date of publication.

### TRADEMARKS

L3K is a registered trademark of SEKISUI.

All other trademarks are the property of their respective owners.

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1. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, 1171-1175 (1994).
2. Bonnichsen, R.K., Theorell, H., *An Enzymatic Method for the Micro Determination of Ethanol*, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
3. Lundquist, F., *The Determination of Ethyl Alcohol in Blood and Tissues*. Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
4. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, 1990.
6. *Criminal Code*, R.S.C. 1985, c. C-46, s. 253(1)(b).

## PATENT NUMBER

US Patent No. 6,380,380 B1

MDSS GmbH  
 Schiffgraben 41,  
 30175 Hannover, Germany  
**Phone:** (+49)-511-6262 8630  
**Fax:** (+49)-511-6262 8633

IN27330-23  
 July 12, 2024

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.  
 ©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

### The Americas

SEKISUI Diagnostics P.E.I. Inc.  
 70 Watts Avenue  
 Charlottetown, PE C1E 2B9  
 Canada  
 Phone: 800-565-0265  
 Fax: 902-628-6504  
 Email: questions@sekisuidiagnostics.com  
 techservices@sekisuidiagnostics.com

### International

SEKISUI Diagnostics (UK) Limited  
 Liphook Way  
 Allington, Maidstone  
 KENT, ME16 0LQ, UK  
 Email: info@sekisuidiagnostics.com

**SEKISUI**  
 DIAGNOSTICS

sekisuidiagnostics.com