EN

AMMONIA L3K® ASSAY

CATALOGUE NUMBER: 293-10 SIZE: 2 x 20 mL

NOTE: Changes are highlighted.

INTENDED USE

For the *in vitro* quantitative measurement of ammonia concentration in plasma.

TEST SUMMARY


Circulating ammonia in normal individuals is relatively low, despite the fact that ammonia is continuously produced from dietary and amino acid metabolism. Blood ammonia measurements have been used in the diagnosis of coma associated with hepatic dysfunction caused by cirrhosis and neoplasms. The measurement of ammonia is very useful in the diagnosis and prognosis of Reye's Syndrome.

The assay of blood ammonia has always been a tedious and time consuming process. Ammonia assays have generally been based on two approaches: the diffusion of ammonia from an alkaline medium with trapping in acid⁽¹⁾ or separation of ammonia from a sample using ion exchange resin.⁽²⁾

This assay is an enzymatic method⁽³⁾ which requires no sample preparation and employs glutamate dehydrogenase and a stabilized NADPH analog⁽⁴⁾ which is easy to use and applicable to routine instrumentation.

TEST PRINCIPLE

GLDH

Ammonia reacts with α -ketoglutarate and reduced cofactor to form L-glutamate and the cofactor. The reaction is catalyzed by glutamate dehydrogenase. The decrease in absorbance due to the oxidation of the reduced cofactor can be monitored at 340 nm or 380 nm and is proportional to the ammonia concentration.

REAGENTS

Ammonia L3K® Reagent: A solution containing a buffer (pH 8.0 at 25°C), 10 mmol/L α -ketoglutarate, ≥ 24 KU/L GLDH (microbial), ~ 0.2 mmol/L NADPH analog, stabilizers, a preservative and a detergent.

Ammonia Calibrator: 1 x 15 mL of a solution of ammonium sulfate (ammonia) 5.0 $\mu\text{g/mL}$ (294 $\mu\text{mol/L}$).

WARNINGS & PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

R₁ ONLY

AMMONIA L3K® Reagent

Contains: tertiary-octylphenoxypoly(ethoxyethanol) (CAS No) 9036-19-5

Hydrochloric acid (CAS No) 7647-01-0

Hazard statements

H401 - Toxic to aquatic life

H412 - Harmful to aquatic life with long lasting effects

Precautionary statements

P273 - Avoid release to the environment.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE & STABILITY

Reagents are ready for use. Supplied reagent is stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

SPECIMEN

Fresh, clear, nonhemolysed / non-lipemic, EDTA plasma or lithium heparinized plasma are acceptable samples. All specimens must be centrifuged using adequate force and time to ensure platelets removal and plasma clarity. Serum and ammonium heparinized plasma are not an acceptable sample types. The accuracy of ammonia measurement is extremely dependent on specimen collection.

Plasma collected into an EDTA or heparin (not ammonium heparin) evacuated tube is recommended. Release the residual vacuum immediately, place the sample on ice, and deliver to the lab as quickly as possible. Separate the plasma from the sample without delay. Do not use hemolyzed or lipemic samples.

SAMPLE STORAGE

Samples should be analyzed within 15 minutes. If this is not possible, samples may be tightly stored at 2-8°C for 2 hours.^(5,6)

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁷⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Avoid ammonia contamination from the laboratory environment. Heavy metals will interfere in the reaction by inhibiting GLDH. Hemolysed samples should not be used as erythrocytes contain larger amounts of ammonia than are found in plasma.⁽⁸⁾

Interferences from icterus, ascorbic acid, plasma pyruvate and plasma lactate were evaluated for this Ammonia method on an ADVIA® 1650 analyzer using a significance criterion of >10% variance from control.

Hemolyzed or lipemic samples should not be used.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
0.87 $\mu\text{g/mL}$	51.1 $\mu\text{mol/L}$	Conjugated Bilirubin	40 mg/dL	474 $\mu\text{mol/L}$
3.65 $\mu\text{g/mL}$	214.2 $\mu\text{mol/L}$	Unconjugated Bilirubin	40 mg/dL	684 $\mu\text{mol/L}$
2.44 $\mu\text{g/mL}$	143.3 $\mu\text{mol/L}$	Ascorbic Acid	3000 $\mu\text{g/dL}$	170 $\mu\text{mol/L}$
2.16 $\mu\text{g/mL}$	126.8 $\mu\text{mol/L}$	Plasma Pyruvate	6.6 mg/dL	0.75 mmol/L
2.72 $\mu\text{g/mL}$	159.7 $\mu\text{mol/L}$	Plasma Lactate	200 mg/dL	22.2 mmol/L

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

Samples containing the following should not be used: Sulfapyridine, Sulfasalazine, Temozolomide, (5-(3-methyltriazin-1-yl)imidazole-4-carboxamide.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' Ammonia L3K® reagent and calibrator.

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Alternative calibration material (if kit calibrator is not appropriate for use).
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:12.7 and a primary wavelength of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The kit calibrator or another calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the ammonia concentration of each sample.

TEST LIMITATIONS

A sample with an ammonia concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

Do not use samples which are not clear; including hemolysed samples, lipemic samples, and samples with visible turbidity.

REFERENCE INTERVALS

0.20 - 0.80 µg/mL (12 - 47 µmol/L)⁽³⁾

These values are suggested ranges based on existing literature. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on an ADVIA® 1650 analyzer unless otherwise stated.

RESULTS

Ammonia concentration is reported as µg/mL (µmol/L). To convert µg/mL to µmol/L, multiply by 58.71.

REPORTABLE RANGE (CLSI EP6)⁽⁷⁾

The linearity of the procedure described is 20.00 µg/mL (1174.2 µmol/L). The lower limit of detection of the procedure described is 0.07 µg/mL (4.1 µmol/L). The limit of quantitation of the procedure described is 0.15 µg/mL (8.8 µmol/L). This data results in a reportable range of 0.15 to 20.00 µg/mL (8.8 - 1174.2 µmol/L).

PRECISION STUDIES (CLSI EP5)⁽⁷⁾

Total precision was collected on three aqueous standards in forty runs conducted over twenty days.

Concentration		Total SD		Total CV%	Concentration		Within Run SD		Within Run CV%
µg/mL	µmol/L	µg/mL	µmol/L		µg/mL	µmol/L	µg/mL	µmol/L	
0.50	29.4	0.02	1.1	3.7	0.48	28.2	0.02	1.1	3.7
1.58	92.8	0.02	1.2	1.3	2.38	139.7	0.03	1.8	1.3
5.08	298.2	0.04	2.1	0.7	5.08	298.3	0.03	1.6	0.5

Within run precision data was collected on three aqueous standards, each standard was run twenty times in a single assay.

ACCURACY (CLSI EP9)⁽⁷⁾

The performance of this method (y) tested on an ADVIA® 1650 was compared with the performance of another commercially available ammonia method (x) tested on a Roche/Hitachi 911. Forty patient plasma samples ranging from 0.16-18.12 µg/mL (9.4-1063.8 µmol/L) gave a correlation coefficient of 0.9997. Linear regression analysis gave the following equation:

This method = 0.997 (similar Ammonia method) + 0.15 µg/mL (8.8 µmol/L).

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

All other trademarks are the property of their respective owners.

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REFERENCES/ BIBLIOGRAFÍA

1. Conway, E.J., Biochem J. 29:27 (1935).
2. Kingsley, G.R., and Tager, H.S., Standard Methods of Clinical Chemistry 6:115, Washington, D.C. 1970, American Assoc. of Clinical Chemistry.
3. Ratcliff, C.R., and Hall, F.F., Selected Methods of Clinical Chemistry 9:85, Edited by Willard R. Faulkner and Samuel Meites, American Association for Clinical Chemistry, Washington, D.C. (1982).
4. U.S. Patent No. 5,801,006.
5. Tietz, N.W., Tietz Textbook of Clinical Chemistry, 3rd Edition, W.B. Saunders, Philadelphia, PA (1999).
6. Heil, W., et. al. "Reference Ranges for Adults and Children", Roche Diagnostics, 2004.
7. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
8. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Washington, Third Edition, (1990).

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

IN29310-19
July 15, 2024

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com