EN

## CARBON DIOXIDE L3K<sup>®</sup> ASSAY

**CATALOGUE NUMBER:** 299-30      **SIZE:** 6 x 30 mL  
 299-17      5 x 100 mL  
 299-50      4 x 500 mL  
 299-55      1 x 500 mL

**Note:** Changes are highlighted.

### INTENDED USE

For the *in vitro* quantitative measurement of carbon dioxide in serum and plasma.

### TEST SUMMARY


Diagnostic determination of serum or plasma carbon dioxide (CO<sub>2</sub>) is used in conjunction with other clinical and laboratory testing for the evaluation of acid-base balance. Blood pH is buffered by the presence of the weak carbonic acid (H<sub>2</sub>CO<sub>3</sub>) and conjugate base bicarbonate (HCO<sub>3</sub><sup>-</sup>). As CO<sub>2</sub> concentration in blood increases, so does hydrogen ion (H<sup>+</sup>) concentration. Respiration rate, which is controlled by CO<sub>2</sub> concentration in the blood, compensates by increasing or decreasing in an attempt to maintain equilibrium by promoting CO<sub>2</sub> elimination or retention respectively.

High CO<sub>2</sub> content is often indicative of respiratory acidosis or metabolic alkalosis. Low CO<sub>2</sub> content is often indicative of respiratory alkalosis or metabolic acidosis. Further laboratory testing is necessary to differentiate between metabolic and respiratory conditions.<sup>(1)</sup>

Classic techniques for the measurement of carbon dioxide (CO<sub>2</sub>) involve the addition of acid to liberate the carbon dioxide and the measurement of carbon dioxide thus released by manometric, volumetric, or titrimetric techniques. These procedures are both time consuming and cumbersome. SEKISUI Diagnostics' Carbon Dioxide L3K<sup>®</sup> assay is an enzymatic procedure, employing phosphoenolpyruvate carboxylase (PEPC)<sup>(2)</sup> and a stabilized NADH analog,<sup>(3)</sup> which is easy to use and applicable to routine laboratory instrumentation.

### TEST PRINCIPLE

PEPC catalyses the first reaction which produces oxaloacetate. In the presence of MDH, the reduced cofactor is oxidized by oxaloacetate. The decrease in concentration of the reduced cofactor is monitored between 405 and 415 nm and is proportional to the total carbon dioxide concentration in the sample.

PEPC is specific for the bicarbonate ion (HCO<sub>3</sub><sup>-</sup>) and its action disturbs the following equilibrium which results in conversion of the CO<sub>2</sub> to HCO<sub>3</sub><sup>-</sup>.

### REAGENTS

Carbon Dioxide Reagent: A solution containing buffer (pH 7.5 at 25°C), 12.5 mmol/L PEP, >400 U/L PEPC (microbial), >4100 U/L malate dehydrogenase (mammalian), 0.6 mmol/L NADH analog, activators, stabilizers, a surfactant, and a preservative.

### WARNINGS & PRECAUTIONS FOR USE

**IVD**

For *in vitro* diagnostic use

**R<sub>1</sub> ONLY**

Carbon Dioxide L3K<sup>®</sup> Reagent

### Warning

Contains: ethanediol; ethylene glycol (CAS No) 107-21-1

### Hazard statements

H373 - May cause damage to organs through prolonged or repeated exposure

### Precautionary statements

P260 - Do not breathe dust/fume/gas/mist/vapors/spray.

P314 - Get medical advice/attention if you feel unwell.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

When using, do not eat or drink.

Avoid contact with skin and eyes.

May cause damage to organs (kidney) through prolonged or repeated exposure. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection.

See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.

Supplied reagent stable at 2-8°C until expiry date.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparinized plasma. The specimen should be promptly separated from the clot or cells and stored tightly sealed to prevent the loss of carbon dioxide.<sup>(6)</sup>

### SAMPLE STORAGE

Samples should be stored at 2-8°C.<sup>(7)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(4)</sup>

Cross Contamination studies have not been performed on automated instruments. Certain reagent / instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemolysis were evaluated for this carbon dioxide method on a Roche/Hitachi<sup>®</sup> 717 analyzer using a significance criterion of >10% variance from control. Data for plasma is similar to that presented for serum.

| Concentration of Analyte |            | Substance Tested | Concentration of Interferent Where Interference is Insignificant |  |
|--------------------------|------------|------------------|--|--|
| Conventional Units       | SI Units   |                  |  |  |
| 26.3 mmol/L              | 26.3 mEq/L | Hemoglobin       | 600 mg/dL  | 93 µmol/L                                      |
| 27.1 mmol/L              | 27.1 mEq/L | Bilirubin        | 40 mg/dL   | 684 µmol/L                                     |
| 24.6 mmol/L              | 24.6 mEq/L | Intralipid       | 1000 mg/dL   | 3000 mg/dL (34 mmol/L) Simulated Triglycerides |

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

The information presented above is based on results from SEKISUI Diagnostic's studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(5)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

SEKISUI Diagnostics' Carbon Dioxide L3K<sup>®</sup> Reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
- Calibration material.
- Quality Control materials.

### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a fixed rate test mode, with a sample to reagent ratio of 1:100, and a wavelength reading of 405 nm or 415 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800) 565-0265. Outside Canada and the U.S., please contact your local distributor.

## CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used.

## QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the carbon dioxide concentration of each sample.

## TEST LIMITATIONS

A sample with a carbon dioxide concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS

22-29 mmol/L (22-29 mEq/L)<sup>(6)</sup>

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on Roche/Hitachi<sup>®</sup> 717 analyzer unless otherwise stated.

## RESULTS

Carbon dioxide concentration is reported as mmol/L (mEq/L).

## REPORTABLE RANGE (CLSI EP6)<sup>(4)</sup>

The linearity of the procedure described is 50.0 mmol/L (mEq/L). The lower limit of detection of the procedure described is 2.9 mmol/L (mEq/L). This data results in a reportable range of 2.9 to 50.0 mmol/L (mEq/L).

## PRECISION STUDIES (CLSI EP5)<sup>(4)</sup>

Total precision data was collected on two concentrations of control sera in triplicate over five runs, for a total of 15 sample values.

| Concentration  | Total SD       | Total CV% | Concentration  | Within Run SD  | Within Run CV% |
|----------------|----------------|-----------|----------------|----------------|----------------|
| mmol/L (mEq/L) | mmol/L (mEq/L) |           | mmol/L (mEq/L) | mmol/L (mEq/L) |                |
| 14             | 0.24           | 1.7       | 14             | 0.18           | 1.3            |
| 25             | 0.27           | 1.1       | 24             | 0.40           | 1.7            |

Within run precision data was collected on two concentrations of control sera, each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(4)</sup>

The performance of this method (y) was compared with the performance of a similar carbon dioxide method (x) on a Cobas Mira<sup>®</sup> analyzer.

Forty-five serum samples ranging from 15.4 to 43.7 mmol/L (mEq/L) were tested and gave a correlation coefficient of 0.9861. Linear regression analysis gave the following equation:

$$\text{This method} = 1.00(\text{reference method}) - 0.42 \text{ mmol/L (mEq/L)}$$

The performance of this method (y) with plasma samples (y) was compared to this method with serum samples (x) on an Advia<sup>®</sup> 1650 analyzer.

Twenty-five patient serum and plasma samples ranging from 17.1 to 25.3 mmol/L (mEq/L) gave correlation coefficient of 0.9636. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 0.950[\text{this method (serum)}] + 1.35 \text{ mmol/L (mEq/L)}$$

## TRADEMARKS

L3K is a registered trademark of SEKISUI. All other trademarks are the property of their respective owners.

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1. Tietz, N.W. (Ed.), *Fundamentals of Clinical Chemistry*, W.B. Saunders Co., Toronto, 625-636, 937 (1970).
2. Norris, K.A., Atkinson, A.R., Smith, W.G., Clin. Chem. 21 (1975).
3. US Patent #5,801,006
4. *CLSI Guidelines and Standards*, Clinical and Laboratory Standards Institute, Wayne, PA.
5. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington (1990).
6. Tietz, N.W. (Editor), *Textbook of Clinical Chemistry*, W.B. Saunders Co., Philadelphia (1986).
7. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations. Geneva: World Health Organization; 2002:39.

**PATENT NUMBER**  
US Patent No. 5,801,006

MDSS GmbH  
Schiffgraben 41,  
30175 Hannover, Germany  
**Phone:** (+49)-511-6262 8630  
**Fax:** (+49)-511-6262 8633

The word SEKISUI and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

IN29930-19  
September 23, 2022

**The Americas**  
**SEKISUI** Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada  
Phone: 800-565-0265  
Fax: 902-628-6504  
Email: [questions@sekisuidiagnostics.com](mailto:questions@sekisuidiagnostics.com)  
[techservices@sekisuidiagnostics.com](mailto:techservices@sekisuidiagnostics.com)

**International**  
**SEKISUI** Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK  
Email: [info@sekisuidiagnostics.com](mailto:info@sekisuidiagnostics.com)

**SEKISUI**  
DIAGNOSTICS

[sekisuidiagnostics.com](http://sekisuidiagnostics.com)