EN

ALANINE AMINOTRANSFERASE-SL ASSAY

CATALOGUE NUMBER: 318-30 SIZE: R1: 3 x 100 mL, R2: 1 x 75 mL

Note: Changes are highlighted

INTENDED USE



For the *in vitro* quantitative measurement of alanine aminotransferase in serum.

TEST SUMMARY

Henley and Pollard⁽¹⁾ and Wroblewski and Ladue⁽²⁾ developed kinetic procedures for assaying alanine aminotransferase, formerly called glutamate pyruvate transaminase. The procedures were based on the oxidation of NADH by lactate dehydrogenase (LDH). This procedure incorporates significant modifications including optimization of substrate concentrations and replacement of phosphate with Tris as the buffer. These modifications meet the recommendations outlined by the International Federation of Clinical Chemistry (IFCC) for ALT determinations.⁽³⁾

Elevated concentrations of alanine aminotransferase result from hepatic disease. In cases of hepatic necrosis, elevations of ALT concentrations occur prior to the onset of clinical symptoms such as jaundice. Marginally elevated concentrations of ALT may be observed with myocardial infarction. To a lesser extent elevated concentrations are indicative of diseases of some internal organs.⁽⁴⁾

TEST PRINCIPLE

The ALT catalyzes the conversion of L-alanine and alpha-ketoglutarate to pyruvate and L-glutamate. In reaction II, LDH catalyzes the oxidation of NADH to NAD⁺. The rate of decrease in absorbance of the reaction mixture at 340 nm, due to the oxidation of NADH is directly proportional to the ALT activity.

REAGENTS

ALT-SL Enzyme Reagent (R1): A solution containing a buffer (pH 7.55 at 25°C), 700 mmol/L L-alanine, ≥2300 U/L lactate dehydrogenase (microbial), and a preservative.

ALT-SL Substrate Reagent (R2): A solution containing 85 mmol/L 2-oxoglutarate, 1 mmol/L NADH, and a preservative.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

Rx ONLY

ALT-SL Substrate Reagent (R2)

Warning

Contains: KETOGLUTARIC ACID (CAS No) 328-50-7

Hazard statements

H319 - Causes serious eye irritation

Precautionary statements

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 - If eye irritation persists: Get medical advice/attention.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE, AND STABILITY

Reagents are ready for use. A single working reagent may be prepared by mixing four parts R1 with one part R2. The working reagent is stable for four weeks at 2-8°C. Supplied reagent is stable at 2-8°C until expiry date. Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum.

SAMPLE STORAGE

ALT is stable for 24 hours at 20-25°C and 7 days at 2-8°C.⁽⁵⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁶⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

This procedure measures total ALT. Red blood cells contain high concentrations of ALT; therefore, hemolysis can elevate results.

Interferences from icterus and lipemia were evaluated for this method on an Roche/Hitachi[®] analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
62 U/L	Bilirubin	40 mg/dL	684 µmol/L
100 U/L	Intralipid	400 mg/dL	1200 mg/dL (13.6 mmol/L) Simulated Triglycerides

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

Samples containing the following should not be used: Sulfasalazine, Temozolomide, and Hydroxocobalamin.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁷⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' ALT-SL reagent.

MATERIALS REQUIRED (BUT NOT PROVIDED)

- 1) Automated analyzer capable of accurately measuring absorbance at an appropriate wavelength as per instrument application.
- 2) Calibration material (if applicable).
- 3) Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:14:3.5 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada

and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory.

CALCULATIONS

The analyzer automatically calculates the alanine aminotransferase (ALT) concentration of each sample.

TEST LIMITATIONS

A sample with an ALT concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽⁴⁾

Newborn: 7-40 U/L (37°C)
Adult: 10-35 U/L (37°C)

ALT activity in men are generally higher than levels in women.

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

RESULTS

ALT activity is reported as U/L.

REPORTABLE RANGE (CLSI EP6)⁽⁶⁾

The linearity of the procedure described is 600 U/L. The lower limit of detection of the procedure described is 10 U/L. This data results in a reportable range of 10 to 600 U/L.

PRECISION STUDIES (CLSI EP5)⁽⁶⁾

Data was collected on two levels of a control sera using a single lot of reagent in 40 runs conducted over 20 days.

Concentration U/L	Total SD U/L	Total CV %	Within Run SD U/L	Within Run CV %
46	1.7	3.6	1.1	2.4
144	3.5	2.5	1.1	0.7

ACCURACY (CLSI EP9)⁽⁶⁾

The performance of this method (y) was compared with the performance of a similar ALT method (x) on a Roche/Hitachi® analyzer. Forty-nine patient serum samples ranging from 12 to 304 U/L were tested and gave a correlation coefficient of 0.9945. Linear regression analysis gave the following equation:

$$\text{This method} = 1.06 (\text{reference method}) - 2.8 \text{ U/L.}$$

All other trademarks are the property of their respective owners.

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REFERENCES / BIBLIOGRAFÍA

1. Henley, K.S., Pollard, H.M. *A New Method for the Determination of Glutamic Oxalacetic and Glutamic Pyruvic Transaminase in Plasma*, J. Lab. and Clin. Med., 46, 785 (1955).
2. Wroblewski, F., Ladue, J.S. *Serum Glutamic Pyruvic Transaminase in Cardiac and Hepatic Disease*, Proc. Soc. Exp. Biol. Med. 91, 569 (1956).
3. Bergmeyer, H.U., and Horder, M. *IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, Part 3 IFCC Method For Alanine Aminotransferase*, J. Clin. Chem. Clin. Biochem. 18, 521 (1980).
4. Burtis, C.A., Ashwood, E.R. (Ed.), *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, p. 790-791, 797 (1994).
5. Kaplan, Lawrence A. and Pesce, Amadeo J., (Ed), *Clinical Chemistry: Theory, Analysis, Correlation, 3rd ed.* Mosby-Year Book Inc., St. Louis, Missouri (1996), p. 518.
6. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
7. Young, D.S. *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Washington, Third Edition (1990).

IN31810-20
June 25, 2024

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com