## CREATINE KINASE-SL ASSAY

**CATALOGUE NUMBER:** 326-10      **SIZE:** R1: 1 x 100 mL  
R2: 1 x 25 mL

**Note:** Changes are highlighted.

### INTENDED USE

For the *in vitro* quantitative measurement of creatine kinase activity in serum.

### TEST SUMMARY

Elevated concentrations of serum creatine kinase are found in all types of muscular dystrophy and disease involving the destruction of muscle tissues. In the case of progressive muscular dystrophy, elevated CK concentrations may be present before clinical symptoms appear. Elevated CK concentrations are observed after myocardial infarction and in cases of acute cardiovascular disease and Reye's Syndrome. Neurogenic muscle diseases such as multiple sclerosis and polio do not produce elevated CK concentrations. CK activity is in inverse relation to thyroid activity, thus, in cases of hypothyroidism, elevated concentrations of CK would be observed.<sup>(1)</sup>

The method of assaying CK using creatine phosphate and adenosine diphosphate (ADP) as substrates instead of creatine and adenosine triphosphate (ATP) was first described by Oliver.<sup>(2)</sup> This method employs the optimized conditions as jointly developed by the Scandinavian Committee on Enzymes and the German Society for Clinical Chemistry.<sup>(3, 4, 5)</sup> In the procedure N-acetylcysteine (NAC) is the thiol activator, and adenosine monophosphate (AMP) and p<sup>1</sup>, p<sup>3</sup>-di(adenosine-5') pentaphosphate (Ap5A) are added to inhibit interference caused by adenylate kinase activity.

### TEST PRINCIPLE

The reactions proceed as follows:

The rate of increase in absorbance at 340 nm due to the formation of NADPH is directly proportional to the creatine kinase activity.

### REAGENTS

**Creatine Kinase-SL Buffer Reagent (R1):** A solution containing a buffer (pH 6.0), 60 mmol/L Imidazole, 27 mmol/L glucose, 27 mmol/L NAC, 14 mmol/L magnesium acetate, 2 mmol/L EDTA-Na<sub>2</sub>, 2.7 mmol/L NADP, ≥5 KU/L hexokinase (yeast).

**Creatine Kinase-SL Substrate Reagent (R2):** A solution containing 160 mmol/L Imidazole, 11 mmol/L ADP, 28 mmol/L AMP, 55 μmol/L Ap5A, ≥14KU/L G-6-PDH (microbial), 2 mmol/L EDTA-Na<sub>2</sub>, 160 mmol/L creatine phosphate.

### WARNINGS & PRECAUTIONS FOR USE

For *in vitro* diagnostic use

### CREATINE KINASE-SL Buffer Reagent (R1)

Danger

Contains: imidazole (CAS No) 288-32-4

Hazard statements

H360 - May damage fertility or the unborn child

Precautionary statements

P201 - Obtain special instructions before use.

P202 - Do not handle until all safety precautions have been read and understood.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P308+P313 - If exposed or concerned: Get medical advice/attention.

P405 - Store locked up.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

### CREATINE KINASE-SL Substrate Reagent (R2)

Danger

Contains: imidazole (CAS No) 288-32-4

Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

H360 - May damage fertility or the unborn child

Precautionary statements

P201 - Obtain special instructions before use.

P202 - Do not handle until all safety precautions have been read and understood.

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 - If exposed or concerned: Get medical advice/attention.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

P405 - Store locked up.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

The reagents are provided in a ready to use format. A working reagent may be prepared by mixing 4 parts CK-SL Buffer Reagent (R1) with 1 part CK-SL Substrate Reagent (R2).

The reagents included are stable until the expiry date stated on the labels at 2-8°C. Protect the reagent from sunlight. The working reagent is stable for 21 days at 2-8°C.

Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum.

### SAMPLE STORAGE

The CK activity in serum decreases rapidly with increasing temperature and the sample should be separated and cooled to 2-8°C as quickly as possible. Macro, or mitochondrial CK, is the most stable isoenzyme followed by CK-3 (MM), CK-2 (MB) and CK-1 (BB). The average stabilities, defined as less than 5% loss of activity, are as follows:<sup>(5, 1)</sup>

CK-MM (CK-3): 24 hours at 18-26°C, 1 week at 4°C, 1 month at -20°C


CK-MB (CK-2): 12 hours at 18-26°C, 3 days at 4°C, 1 month at -20°C

CK-BB (CK-1): less than CK-MB (CK-2)

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(6)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Adenylate kinase interferes with the determination of CK activity by generating additional ATP (see reaction below) and thus increasing the apparent creatine kinase activity.

The addition of AMP and Ap5A to the reagent in this procedure inhibits the adenylate kinase activity.<sup>(4)</sup>

Interference from icterus, lipemia and hemolysis were evaluated for this creatine kinase method on a Roche/Hitachi® analyzer using a significance criterion of > 10% variance from control. Interference data was collected in serum.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
198 U/L	Hemoglobin	200 mg/dL	31 µmol/L
205 U/L	Bilirubin	40 mg/dL	684 µmol/L
188 U/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(7)</sup>

## ANALYTICAL PROCEDURE

### MATERIALS PROVIDED

SEKISUI Diagnostics' Creatine Kinase-SL reagents.

### MATERIALS REQUIRED (BUT NOT PROVIDED)

- Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- Calibration material. (If applicable.)
- Quality Control materials.

### TEST CONDITION

For data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:35 and a wavelength reading of 340 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

### CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

### QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

### CALCULATIONS

The analyzer automatically calculates the creatine kinase activity of each sample.

### TEST LIMITATIONS

A sample with a creatine kinase activity exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value. Note that dilution of highly active sera causes a progressive increase in activity per unit volume.<sup>(8)</sup>

### REFERENCE INTERVALS<sup>(1)</sup>

Males: 38-174 U/L (37°C)

Females: 26-140 U/L (37°C)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

### PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

### RESULTS

Creatine kinase activity is reported as U/L.

### REPORTABLE RANGE (CLSI EP6)<sup>(6)</sup>

Reportable range using automated procedures will depend on the sample to reagent ratio used. The automated procedure described gives a reportable range from 2-1500 U/L.

### PRECISION STUDIES (CLSI EP5)<sup>(6)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Activity	Total SD	Total CV %	Within Run SD	Within Run CV %
	U/L		U/L	
288	7.1	2.5	1.3	0.4
748	15.7	2.1	2.3	0.3

### ACCURACY (CLSI EP9)<sup>(6)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Fifty-five patient serum samples ranging from 31-1480 U/L gave a correlation coefficient of 0.9993. Linear regression analysis gave the following equation:

$$\text{This method} = 1.02 (\text{reference method}) + 1.3 \text{ U/L.}$$

The information presented above is based on results from SEKISUI Diagnostics'

studies and is current at the date of publication.

All other trademarks are the property of their respective owners.

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
Para uso en diagnóstico *in vitro*

CREATINE KINASE-SL Buffer Reagent (R1)

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#### CREATINE KINASE-SL Substrate Reagent (R2)

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### REFERENCES/ BIBLIOGRAFÍA

- Burtis, C.A., Ashwood, E.R. (Ed.), Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Toronto, (1994).
- Oliver, I.T., A Spectrophotometric Method for the Determination of Creatine Phosphokinase and Myokinase, Biochem. J.61, 116 (1955).
- Gerhart, W., Waldenstrom, J., Gruber, W., Scand, J. Clin Lab. Invest. 39, 737-742 (1979).
- Szasz, G., Gerhardt, W., Gruber, W., Clin Chem. 23, 1888-1892 (1972).
- Recommendation of the German Chemical Society Standardization of Methods for the Estimation of Enzyme Activities in Biological Fluids, J. Clin. Chem., Clin Biochem. 15, 225-260 (1977).
- CLSI Guidelines and Standards, Clinical and Laboratory Standards Institute, Wayne, PA.
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Third Edition, Washington, 1990.
- Graig, F.A., Smith, I.C., Folds, F.F., Clin. Chem. Acta. 15, 107 (1967).

MDSS GmbH  
Schiffgraben 41,  
30175 Hannover, Germany  
**Phone:** (+49)-511-6262 8630  
**Fax:** (+49)-511-6262 8633

The word SEKURE and the Sekure logo are trademarks of SEKISUI Diagnostics, LLC.

©2024 SEKISUI Diagnostics P.E.I. Inc. - All rights reserved.

IN32610-21  
June 25, 2024

The Americas  
SEKISUI Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada  
Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
techservices@sekisuidiagnostics.com

**International**  
SEKISUI Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK  
Email: info@sekisuidiagnostics.com

**SEKISUI**  
DIAGNOSTICS

sekisuidiagnostics.com