EN

ALKALINE PHOSPHATASE-SL ASSAY

CATALOGUE NUMBER: 328-30 SIZE: R1: 3 x 100 mL
R2: 1 x 75 mL

Note: Changes are highlighted.

INTENDED USE


For the IN VITRO quantitative measurement of alkaline phosphatase in serum.

TEST SUMMARY

Elevated alkaline phosphatase activity in serum is of interest in the diagnosis of several general disease conditions including hepatobiliary disease and bone disease associated with increased osteoblastic activity. Alkaline phosphatase activity in serum may be elevated due to obstructive jaundice, occlusion of the common bile or hepatic duct, and cirrhosis.⁽¹⁾

Alkaline phosphatase (ALP) activity was first measured by Kay.⁽²⁾ Since that time many substrates such as glycerol phosphate and phenyl phosphate have been used. Bessey, Lowry, and Brock⁽³⁾ introduced a more sensitive substrate p-nitrophenyl phosphate (p-NPP). Several recommendations have been made on the optimum conditions for the measurement of ALP in serum. These suggestions have been put forth by the German Society for Clinical Chemistry⁽⁴⁾ as well as by the committee on enzymes of the Scandinavian Society for Clinical Chemistry⁽⁵⁾. This method is a modification of the recommendations of the International Federation of Clinical Chemistry (IFCC) for ALP measurement in serum.⁽⁶⁾

TEST PRINCIPLE

Alkaline phosphatase hydrolyzes p-NPP to form the yellow chromogen p-nitrophenol according to the equation.

The rate of increase in absorbance of the reaction mixture at 415 nm, due to the formation of p-nitrophenol, is proportional to the alkaline phosphatase activity.

REAGENTS

R1: ALP-SL Buffer Reagent, R2: ALP-SL Substrate Reagent

Concentrations in the test: 0.67 mol/L 2-amino-2-methyl-1-propanol (pH 10.4 at 25°C), 1.6 mmol/L magnesium acetate, 0.8 mmol/L zinc sulphate, 1.6 mmol/L HEDTA, 10 mmol/L p-nitrophenylphosphate, 10 mmol/L phenol and a preservative.

WARNINGS & PRECAUTIONS FOR USE

IVD

For *In Vitro* Diagnostic Use.

Rx ONLY

Avoid contact with skin and eyes.
See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE, AND STABILITY

Reagents are provided in a ready to use format. An ALP working reagent can be prepared by mixing four parts ALP-SL Buffer Reagent (R1) with one part ALP-SL Substrate Reagent (R2).

Supplied reagent is stable at 2-8°C until expiry date. Working reagent is stable at 2-8°C for 28 days or 18-26°C for 7 days.

Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with the Federal, Provincial, State and local regulations.

SPECIMEN

Fresh, clear, unhemolysed serum. EDTA, oxalate, and citrate inhibit the action of ALP, therefore, these anticoagulants should be avoided.

SAMPLE STORAGE

Specimens stored at room temperature should be analyzed within two to three days. Samples held for longer than two to three days may be stored at -25°C for up to one month. It has been found that alkaline phosphatase activity increases significantly as samples are warmed following refrigeration or freezing.⁽⁷⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽⁸⁾

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia, ascorbic acid, and hemolysis were evaluated for this method on a Roche/Hitachi® 717 analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
135 U/L	Hemoglobin	200 mg/dL	31 µmol/L
134 U/L	Bilirubin	40 mg/dL	684 µmol/L
65 U/L	Ascorbic Acid	1200 µg/dL	68 µmol/L
121 U/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁹⁾

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Sekisui Diagnostics' ALP-SL reagent

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material. (If applicable)
3. Quality Control materials.

TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:56:14 and a wavelength reading of (primary/secondary) 415/660 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory. The best results were obtained if the alkaline phosphatase assay is performed at the same time interval following reconstitution of the control serum each day.^(10,11)

CALCULATIONS

The analyzer automatically calculates the alkaline phosphatase concentration of each sample.

TEST LIMITATIONS

A sample with a alkaline phosphatase activity exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹²⁾

Less than 103 U/L (30°C)
Less than 138 U/L (37°C)

These values are suggested guidelines. Alkaline phosphatase activities have been found to vary with age and sex. It is recommended that each laboratory establish its own expected range.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

RESULTS

Alkaline phosphatase activity is reported as U/L.

REPORTABLE RANGE (CLSI EP6)⁽⁸⁾

The linearity of the procedure described is 2000 U/L. The lower limit of detection of the procedure described is 4 U/L. This data results in a reportable range of 4 to 2000 U/L.

PRECISION STUDIES (CLSI EP5)⁽⁸⁾

Data was collected on two concentrations of control sera using a single lot of reagent in 40 runs conducted over 20 days.

Concentration	Total SD	Total CV %	Within Run SD	Within Run CV %
47 U/L	2.1	4.6	1.1	2.4
214 U/L	7.4	3.5	2.1	1.0

ACCURACY (CLSI EP9)⁽⁸⁾

The performance of this method (y) was compared with the performance of a similar alkaline phosphatase method (x) on a Roche/Hitachi® 717 analyzer. Forty patient serum samples ranging from 18-1835 U/L were tested and gave a correlation coefficient of 1.0000. Linear regression analysis gave the following equation:

$$\text{This method} = 1.0276 (\text{reference method}) - 1.0 \text{ U/L.}$$

All other trademarks are the property of their respective owners.

Manufactured by:

SEKISUI
DIAGNOSTICS

The Americas

Sekisui Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

Sekisui Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK

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This product fulfills the requirements of the European Directive for In Vitro Diagnostic Medical Devices.

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
REFERENCES/ REFERENCIAS

1. Burtis, C.A. and Ashwood, E.R. (Eds), *Tietz Textbook of Clinical Chemistry, Second Edition*, W.B. Saunders Co., Philadelphia, (1994).
2. Kay, H.D., *Plasma Phosphatase, 1:Method Determination*, J.Biol. Chem. 89, 235 (1930).
3. Bessey, O.A., Lowry, O.H., Brock, M.J., *A Method for the Rapid Determination of Alkaline Phosphatase with Five Cubic Milliliters of Serum*, J. Biol. Chem. 164, 321 (1946).
4. Recommendations of the German Society of Clinical Chemistry, *Standardization of Methods for the Estimation of Enzyme Activity in Biological Fluids*, J. Clin. Chem. Clin. Biochem. 10, 182-192 (1972).
5. The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology: *Recommended Methods for the Determination of Four Enzymes in Blood*, Scand. J. Clin. Lab. Invest. 33, 291-306 (1974).
6. Tietz, N.W., Rinker, A.D.U., Shaw, L.M. (1983). *IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 5. IFCC Method for Alkaline Phosphatase*, J. Clin. Chem. Biochem. 21, 371-478 (1983).
7. Kaplan, L.A. and Pesce, A.J., *Clinical Chemistry: Theory, Analysis, Correlation, Third Edition*, Mosby, NY, p. 1072 (1996).
8. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
9. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, Third Edition, Washington, 1990.
10. Szasz, G., *Increase of Alkaline Phosphatase Activity in Commercial Reference Sera After Reconstitution*, Scand. J. Clin. Lab. Invest. Suppl. 126, 29, 12.7 (1972) Abstract.
11. Massion, C.G., Frankenfeld, J.K., *Alkaline Phosphatase Lability in Fresh and Frozen Human Serum and in Lyophilized Control Material*, Clin. Chem. 18, 366 (1972).
12. Pesce, A.J., Kaplan, L. A., *Methods In Clinical Chemistry*, Toronto, C.V. Mosby Co., 1078 (1987).

MDSS GmbH
Schiffgraben 41,
30175 Hannover, Germany
Phone: (+49)-511-6262 8630
Fax: (+49)-511-6262 8633

IN32830-21
July 12, 2023

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.