EN

## GGT Assay

CATALOGUE NUMBER: 334-10      SIZE: R1: 1 x 100 mL, R2: 1 x 20 mL

**Note:** Changes are highlighted

### INTENDED USE

For the *in vitro* quantitative measurement of  $\gamma$ -Glutamyltransferase (GGT) activity in serum and plasma.

### TEST SUMMARY


Although renal tissue has the highest concentration of GGT, blood concentrations primarily originate from the hepatobiliary system. GGT is a valuable tool for uncovering otherwise silent liver disease. Normally present in low concentrations in the blood, GGT concentrations rise when acute damage to the liver or bile duct occurs. Measurements of its activity may also indicate cholestasis, cirrhosis, primary/secondary liver tumors, hepatitis or hepatotoxic drug damage.<sup>(1)</sup>

$\gamma$ -Glutamyltransferase (GGT) was initially purified and characterized by Szewczuk and Baranowski.<sup>(2)</sup> GGT catalyzes the transfer of the  $\gamma$ -glutamyl group from a  $\gamma$ -glutamyl peptide to an amino acid of another peptide.

The measurement of GGT using L- $\gamma$ -glutamyl-P-nitroanilide as a substrate was introduced by Orlowski and Meizler<sup>(3)</sup> and adapted for measurement by Szasz.<sup>(4)</sup> The method was modified by Persijn and Vanderslik<sup>(5)</sup> to use a more soluble L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide substrate. Reagents in this assay use the carboxylated substrate for detection of GGT following a modification of the International Federation of Clinical Chemistry rapid kinetic procedure for serum GGT.<sup>(6)</sup>

### TEST PRINCIPLE

#### GGT

GGT catalyzes the transfer of the glutamyl group from L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide to glycylglycine, forming L- $\gamma$ -glutamylglycylglycine and 5-amino-2-nitrobenzoate. The 5-amino-2-nitrobenzoate absorbs strongly at 405 nm. Formation of this product is proportional to GGT activity and is measured kinetically at 405 nm.

### REAGENTS

GGT Buffer Reagent (R1): A solution containing 100 mmol/L buffer (pH 8.2 at 25°C), 150 mmol/L glycylglycine and a preservative.

GGT Substrate Reagent (R2): A solution containing 25 mmol/L buffer (pH 6.0 at 25°C), 25 mmol/L L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide (Glupa-C), a stabilizer and a preservative.

### WARNINGS & PRECAUTIONS

#### IVD

For *in vitro* diagnostic use

## R<sub>X</sub> ONLY

### GGT Buffer Reagent (R1)

#### Warning

Contains: potassium hydroxide (CAS No) 1310-58-3

#### Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

#### Precautionary statements

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

Avoid ingestion.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

### REAGENT, PREPARATION, STORAGE & STABILITY

Reagents are ready for use.

The unopened reagents included are stable until the expiry date stated on the labels at 2-8°C. Stability claims are based on accelerated stability studies.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

Fresh, clear, unhemolysed serum or lithium heparin plasma.

### SAMPLE STORAGE

Samples may be stored at 18-26°C for 2 days, 0-4°C for 1 week or at -25°C for 1 month.<sup>(7)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(8)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

It has been reported that some antiepileptic drugs (phenytoin, barbituates) produce false elevation of GGT activity.<sup>(9),(10)</sup>

Interferences from icterus, lipemia, hemolysis and ascorbic acid were evaluated for this GGT method on a Roche/Hitachi® 911 using a significance criterion of >10% variance from control. Interference data was collected in serum.

Activity of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant	
39.0 U/L	Bilirubin	40 mg/dL	684 $\mu$ mol/L
51.0 U/L	Hemoglobin	600 mg/dL	93 $\mu$ mol/L
52.5 U/L	Ascorbic Acid	3000 $\mu$ g/dL	170 $\mu$ mol/L
74.5 U/L	Intralipid	1000 mg/dL (33.9 mmol/L) Simulated Triglyceride	3000 mg/dL (33.9 mmol/L) Simulated Triglyceride

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(9)</sup>

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

SEKISUI Diagnostics' GGT reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

- Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
- Calibration material. (If applicable.)
- Quality Control materials.

### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using a kinetic test mode, with a sample to reagent ratio of 1:43 and a wavelength reading of 415 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800) 565-0265. Outside Canada and the U.S., please contact your local distributor.

### CALIBRATION

The frequency of calibration, if necessary, using an automated system is dependent on the system and the parameters used.

## QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state, and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the GGT activity of each sample.

## TEST LIMITATIONS

A sample with a GGT activity exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS<sup>(4)</sup>

Females: 8.8-22.0 U/L at 37°C  
Males: 10.4-33.8 U/L at 37°C

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 911 analyzer unless otherwise stated.

## RESULTS

GGT activity is reported as U/L.

## REPORTABLE RANGE (CLSI EP6)<sup>(8)</sup>

The linearity of the procedure described is 7.0 to 1200.0 U/L. The lower limit of detection of the procedure described is 3.2 U/L. This data results in a reportable range of 7.0 to 1200.0 U/L.

## PRECISION STUDIES (CLSI EP5)<sup>(8)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Activity (U/L)	Total SD (U/L)	Total CV %	Activity (U/L)	Within Run SD (U/L)	Within Run CV %
30.6	1.1	3.6	30.4	0.6	1.9
139.1	5.1	3.7	135.8	0.6	0.5

Within run precision data was collected on two concentrations of control sera each run twenty times in a single assay.

## ACCURACY (CLSI EP9)<sup>(8)</sup>

The performance of this method (y) was compared with the performance of a similar GGT method (x) on a Roche/Hitachi® analyzer. Forty (40) patient serum samples ranging from 13.5 to 924.0 U/L gave a correlation coefficient of 1.0000. Linear regression analysis gave the following equation:

$$\text{This method} = 0.981 \text{ (reference method)} + 0.07 \text{ U/L.}$$

The performance of this method (plasma) was compared with the performance of this method (serum) on a Roche/Hitachi® 911 analyzer. Forty (40) patient serum and plasma samples ranging from 11 to 239 U/L gave a correlation coefficient of 0.9995. Linear regression analysis gave the following equation:

$$\text{This method (plasma)} = 0.991 \text{ (serum)} + 0.4 \text{ U/L.}$$

All trademarks, brands, product names and trade names are the property of their respective companies.

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## REFERENCES / BIBLIOGRAFÍA

1. Burtis, C.A., Ashwood, E.R. (Eds). *Tietz Textbook of Clinical Chemistry*, W.B. Saunders Co., Toronto, pp 848-849 (1994).
2. Szewczuk, A. and Baranowski, T., *Purification and Properties of  $\gamma$ -glutamyl Transpeptidase from Beef Kidney*, Biochem. Z. 338, 317-329 (1963).
3. Orłowski, M. and Meister, A., *Biochem. Biophys. Acta* 73:679 (1963).
4. Szasz, G., *Clin. Chem.* 15: 124 (1976).
5. Persijn, J.P. and Van der Slik, W., *A New Method for the Determination of  $\gamma$ -glutamyltransferase in Serum*, *J. Clin. Chem. Clin. Biochem.*, 14, 421-427 (1976).
6. Shaw, L.M., Stromme, J.H. et.al., *IFCC Methods for the Measurement of Catalytic Concentration of Enzymes*, Part 4, *IFCC Method for  $\gamma$ -Glutamyltransferase*, *J. Clin. Chem. Clin. Biochem.* 21, 533 (1983).
7. Kaplan, L. and Pesce, A. (Eds) , *Clinical Chemistry: Theory, Analysis, Correlation*, Mosby-Year Book Inc, St Louis, Missouri, p 1072 (1996).
8. *CLSI Method Evaluation Protocols*, Clinical and Laboratory Standards Institute, Wayne, PA.
9. Young, D.S., *Effects of Drugs on Clinical Laboratory Tests*, AACC Press, 3rd Edition, Washington DC, (1990).
10. Tietz, N.W., *Clinical Guide to Laboratory Tests*, 3<sup>rd</sup> Edition, W.B. Saunders Company, Philadelphia, Pennsylvania, p. 286 (1995).

IN33410-14  
February 2, 2022

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

### The Americas

SEKISUI Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada  
Phone: 800-565-0265  
Fax: 902-628-6504  
Email: [questions@sekisuidiagnostics.com](mailto:questions@sekisuidiagnostics.com)  
[techservices@sekisuidiagnostics.com](mailto:techservices@sekisuidiagnostics.com)

### International

SEKISUI Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK  
Email: [info@sekisuidiagnostics.com](mailto:info@sekisuidiagnostics.com)

[sekisuidiagnostics.com](http://sekisuidiagnostics.com)