EN

N-geneous® LDL-ST Cholesterol Calibrator

CATALOGUE NUMBER: 80-5666-02 **SIZE:** 3 x 1 mL

Note: Changes are highlighted.

INTENDED USE

For calibration of the N-geneous® LDL-ST Cholesterol assay.

TEST SUMMARY AND PRINCIPLE

See the N-geneous® LDL-ST Cholesterol Reagent package insert.

REAGENTS

Composition

The N-geneous® LDL-ST Cholesterol Calibrator contains human serum, containing lipoproteins from the various lipoprotein classes including low-density lipoproteins and sodium azide as preservative (0.1%).

Note: The value of the N-geneous® LDL-ST Cholesterol Calibrator was assigned by procedures traceable to the CDC LDL cholesterol reference method.²

WARNINGS AND PRECAUTIONS FOR USE

Dispose of contents/container in accordance with local/regional /national/international regulations.

IVD

For *in vitro* diagnostic use

R_x ONLY

Do not pipette by mouth.

Do not use the calibrator after the expiration date printed on the label.

Warning: Human source material. Treat as potentially infectious. Each plasma donor unit used in the preparation of this product has been tested by an FDA approved method and found non-reactive for HBsAg, HCV, HIV 1&2 and HIV-1 Antigen. Because no known test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good laboratory practices using appropriate precautions.¹

Caution: Contains sodium azide, which may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush drain with large volume of water to prevent buildup.

See Safety Data Sheet for additional information.

PRODUCT PREPARATION, STORAGE AND STABILITY

Allow the vial of calibrator to equilibrate to room temperature before reconstitution. Reconstitute each vial with 1.0 mL of distilled or deionized water (or equivalent) using a Class A volumetric pipette. Replace the cap/stopper and allow to stand for 5 minutes at room temperature. Dissolve the contents of the vial by swirling gently to avoid the formation of foam.

DO NOT SHAKE.

Unopened calibrator is stable until the expiration date printed on the label when stored at 2-8°C.

Once reconstituted, the calibrator is stable up to 2 weeks at 2-8°C.

Reconstituted stability of the calibrator may be extended by aliquotting and freezing the reconstituted calibrator at -80°C for 1 month.

Freeze only once.

PRODUCT DETERIORATION

The following may indicate deterioration:

Inability to recover control values.

Presence of turbidity.

DISPOSAL

Product must be disposed of in accordance with all Federal, Provincial, State and local regulations.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

Description	Configuration	Catalog Number
N-geneous® LDL-ST Calibrator	3 x 1 mL	80-5666-02

MATERIALS REQUIRED (BUT NOT PROVIDED)

Description	Configuration	Catalog Number
N-geneous® LDL-ST Cholesterol Reagent 1	1 x 30 mL	80-5727-02
N-geneous® LDL-ST Cholesterol Reagent 2	1 x 10 mL	80-5729-02
N-geneous® LDL-ST Cholesterol Reagent 1	1 x 250 mL	80-5665-01
N-geneous® LDL-ST Cholesterol Reagent 2	1 x 80 mL	80-5714-00

1. Class A volumetric pipettes.
2. Distilled, deionized, Type II water or equivalent.

TEST CONDITION

See the N-geneous® LDL-ST Reagent package insert.

CALIBRATION

The value of the N-geneous® LDL-ST Cholesterol Calibrator, found on the vial label, was assigned by procedures traceable to the CDC LDL cholesterol reference method.² Refer to the instrument manufacturer's recommendation for calibration frequency. If assistance is required, please contact SEKISUI Diagnostics Technical Services at 800-565-0265.

Quality Control values should be within the expected range.

All trademarks, brands, product names and trade names are the property of their respective companies.

FR

N-geneous® LDL-ST Cholesterol Calibrator

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IVD

In vitro diagnostic medical device
Appareil médical de diagnostic *in vitro*.

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- Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," (1988).
- National Reference System for Cholesterol. CRMLN LDL Cholesterol Protocol, May 2004.

Refer to the N-geneous® LDL-ST Cholesterol Reagent package insert for additional references.

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N-geneous® is a registered trademark of SEKISUI Diagnostics, LLC.
©2022 SEKISUI Diagnostics, LLC - All rights reserved.

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The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

IN80566602-5
February 25, 2022

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada
Phone: 800-565-0265
Fax: 902-628-6504
Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com

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N-geneous® est une marque de commerce déposée de SEKISUI Diagnostics, LLC.
©2022 SEKISUI Diagnostics, LLC - Tous droits réservés.