EN

Lipase Color

CATALOGUE NUMBER: 905-B **SIZE:** Reagent 5 x 30 mL
Solvent 1 x 200 mL
Activator 1 x 60 mL
Calibrator 2 x 3 mL

Note: Changes are highlighted.

INTENDED USE

For the quantitative measurement of pancreatic lipase activity in serum or plasma.

TEST SUMMARY

Pancreatic lipase in serum or plasma is closely associated with pancreatic diseases. The activity of this enzyme is measured as an important marker for diagnosing pancreatic diseases and the associated monitoring of therapeutic effects. Pancreatic lipase measurement has been reported using titrimetric, turbidimetric, fluorometric, and colorimetric methodologies.¹⁻⁵

The Pancreatic Lipase Color test is a colorimetric, kinetic assay that uses a clear substrate solution of 1,2-diglyceride that is a “natural” substrate. The assay is highly sensitive and specific for pancreatic lipase, using colipase and deoxycholate as activators.⁴

The assay shows excellent reproducibility and stability. Furthermore, the simplicity of this procedure makes it readily adaptable for use on automated analyzers.

TEST PRINCIPLE

Serum lipase acts on a natural substrate, 1,2-diglyceride, to liberate 2-monoglyceride. This is hydrolyzed by monoglyceride lipase (a highly specific enzyme for monoglyceride) into glycerol and free fatty acid. Glycerol kinase acts on glycerol to form glycerol-3-phosphate, which is in turn acted on by glycerol-3-phosphate oxidase to generate hydrogen peroxide. Peroxidase converts the hydrogen peroxide, 4-Aminoantipyrine and TOOS (N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine) into a quinine dye. The rate of formation of the dye, measured as an increase in absorbance at 550 nm, is proportional to the lipase concentration in the sample.⁴

REAGENTS

Composition

Component	Ingredients	Concentration
Lipase Color Reagent*	1,2-Diglyceride (Egg)	0.077%
	Monoglyceride Lipase (<i>micro-organism</i>)	≤0.88 U/mL
	Glycerol Kinase (<i>S. canus</i>)	≤1.34 U/mL
	Glycerol-3-Phosphate Oxidase (<i>Streptococcus</i> sp.)	≤40 U/mL
	TOOS	0.07%
	ATP (Bacterial)	0.66 mM
	Peroxidase (Horseradish)	≤1.34 U/mL
	Colipase (Porcine)	≤40 U/mL
	Buffer	pH 6.8
	Human Serum Albumin (Human)	0.27%
	Ascorbate Oxidase (Cucumber)	≤2.66 U/mL
	Stabilizers	
	Lipase Color Solvent*	Cholic Acid (Ox or Sheep)
Buffer		pH 6.8
Sodium Azide		0.05%

Composition (continued)

Component	Ingredients	Concentration
Lipase Color Activator	Deoxycholate (Ox or Sheep)	36 mM
	4-Aminoantipyrine	0.12%
	Buffer	pH 8.7
	Sodium Azide	0.05%

The Lipase Color Calibrator contains human pancreatic lipase, serum albumin (bovine), and preservative.

***NOTE:** The concentrations shown for the Lipase Color Reagent and Lipase Color Solvent are the concentrations after reconstitution.

WARNINGS AND PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

R_{ONLY}

Lipase Color Reagent

Danger

Contains: 4-nonylphenol, branched, ethoxylated (CAS No) 127087-87-0; Lipase cofactor; (CAS No) 55126-92-6; Oxidase, glycerol phosphate (CAS No) 9046-28-0

Hazard statements

H334 - May cause an allergy or asthma symptoms or breathing difficulties if inhaled.

H411 - Toxic to aquatic life with long lasting effects.

Precautionary statements

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray.

P273 - Avoid release to the environment.

P284 - [In case of inadequate ventilation] wear respiratory protection.

P304 + P341 - If inhaled: If breathing is difficult, remove person to fresh air and keep comfortable for breathing.

P342 + P311 - If experiencing respiratory symptoms: Call a poison center or doctor.

P391 - Collect spillage.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Do not pipette by mouth.

Do not use the calibrator after the expiration date printed on the label.

Warning: Human source material. Treat as potentially infectious. Each donor unit used in the preparation of this product has been tested by an FDA approved method and found non-reactive for HBsAg, HCV, HIV 1/2 and HIV-1 antigen. Because no known test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good laboratory practices using appropriate precautions.⁶

The enzymes in triglyceride and cholesterol reagents may contaminate the Lipase Color Assay. To avoid contamination, ensure probes, cuvettes or tubes of automated analyzers are thoroughly washed between triglyceride or cholesterol assays and use of the Lipase Color Assay.

Caution: Contains sodium azide, which may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush drain with a large volume of water to prevent buildup. Dispose of in accordance with local, state, and federal regulations.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Lipase Color Reagent: Reconstitute the Reagent vial with the 30 mL Solvent, using a Class A volumetric pipette. Allow to stand a minimum of ten minutes at room temperature. Mix gently by inversion before use.

Lipase Color Solvent: Supplied ready to use.

Lipase Color Activator: Supplied ready to use.

Lipase Calibrator: Reconstitute the calibrator vial with 3.0 mL of distilled, deionized, Type II water or equivalent using a Class A volumetric pipette. Replace the cap/stopper and mix gently by inversion. Allow to stand for a minimum of 10 minutes at room temperature. Mix gently by inversion before use.

DO NOT SHAKE

Unopened Lipase Color Reagent, Solvent, Activator and Calibrator are stable until the expiration date shown on the label when stored at 2-8°C.

DO NOT FREEZE

Once reconstituted, the Lipase Color Reagent is stable for 28 days at 2-8°C.

DO NOT FREEZE

Once opened the Lipase Color Solvent and Activator are stable until the expiration date on the label when stored capped at 2-8°C.

DO NOT FREEZE

Once reconstituted, the Lipase Calibrator is stable up to 28 days at 2-8°C. Reconstituted stability of the calibrator may be extended by freezing the reconstituted calibrator at -20°C for up to 4 months. Freezing in the original vial is recommended.

FREEZE ONLY ONCE

REAGENT DETERIORATION

The following would indicate deterioration:

Presence of extreme turbidity.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State and local regulations.

SPECIMEN

Serum, EDTA-treated plasma or lithium and sodium heparinized plasma are the recommended sample types. Specimens should be collected as per the National Committee for Clinical Laboratory Standards Guideline H4-A3.⁸

Serum: Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection.⁸ (Within 3 hours)

Plasma: Specimens may be collected in EDTA or lithium and sodium heparin. Centrifuge and remove the plasma as soon as possible after collection.⁸ (Within 3 hours)

If the assay is not performed immediately, the serum and plasma must be refrigerated or frozen until use. If not analyzed promptly, specimens may be stored at 2-8°C for up to 3 weeks. If specimens need to be stored for more than 3 weeks, they may be preserved at -20°C or below for up to 3 months. Samples may be frozen once. Refer to NCCLS Document H18-A for further instructions on specimen collection, handling and storage.

ANALYTICAL SPECIFICITY

No significant interference was detected in the Lipase Color assay up to and including the concentrations stated below:

Substance Tested	Concentration with no significant ($\pm 10\%$) interference
Unconj. Bilirubin	20 mg/dL
Conj. Bilirubin	25 mg/dL
Hemoglobin	2000 mg/dL
Triglyceride	1000 mg/dL
Liposyn	1%
Glycerol	250 mg/dL
Ascorbic Acid	50 mg/dL

Samples containing interfering substances greater than the levels listed exhibited bias $>10\%$. Dilutions are not recommended to reduce interferences.

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

Samples containing the following should not be used:
N-acetylcysteine (NAC).

1. Refer to the work of Young et al⁹ for a review of the effects of drugs on clinical laboratory tests.

2. The enzymes in triglyceride and cholesterol reagents may contaminate the Lipase Color Assay. To avoid contamination, ensure probes, cuvettes or tubes of automated analyzers are thoroughly washed between triglyceride or cholesterol assays and use of the Lipase Color Assay.

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

The Lipase Color Reagent, Solvent, Activator and Calibrator are required for the measurement of Lipase. The Lipase Color reagents are packaged and sold together. A configuration of each of the following items will be included in the package you receive.

Description	Configuration	Catalog Number
Lipase Color	Reagent	5 x 30 mL
	Solvent	1 x 200 mL
	Activator	1 x 60 mL
	Calibrator	2 x 3 mL
		905-B

MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Spectrophotometer or other instrument capable of reading at 550 nm.
2. Class A volumetric pipettes.
3. Distilled or deionized water meeting specifications equivalent to USP purified water.
4. Lipase control sera or quality control material (See "Quality Control").

TEST CONDITION

All analyzer applications should be validated in accordance with CLIA recommendations. For assistance with applications on automated analyzers, please contact SEKISUI Diagnostics Technical Services at 800-565-0265.

CALIBRATION

The Lipase Color Calibrator is required for calibration. Refer to the instrument manual for analyzer specific calibrator procedures and for guidance in determining calibration frequency. The calibrator value can be found on the calibrator vial label.

QUALITY CONTROL

Reliability of test results should be routinely monitored with quality control materials or serum pools that reasonably represent performance on patient specimens. Controls or serum pools should be run with each assay to ensure that the reagents are functioning properly and that correct procedures have been followed. Quality control materials are intended for use only as monitors of accuracy and precision. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range confirm procedures were performed correctly and follow normal troubleshooting measures. If the problem persists call SEKISUI Diagnostics Technical Services at 800-565-0265.

Quality control requirements should be established in accordance with local, state and/or Federal regulations or accreditation requirements.

REFERENCE INTERVALS

A normal range study was performed using the Lipase Color assay. A serum range of 21-67 U/L (138 healthy donors) was obtained on the Roche/Hitachi 911 system. Ranges were calculated as recommended by NCCLS guideline C28-A.¹⁰ These results were obtained using a specific lot of reagent. It is recommended that each laboratory establish the normal range for its patient population.

PERFORMANCE CHARACTERISTICS

RESULTS

Unit Definition: One unit (U) is defined as the amount of enzyme activity which liberates 1 μ mole of 2-monoglyceride from 1,2-diglyceride per minute at 37°C. Lipase activity is expressed in U/L.

Patient results may be reported in U/L. To convert from conventional units to S.I. units, multiply the conventional units by 1.67×10^{-8} Katal/U. $U/L \times (1.67 \times 10^{-8} \text{ Katal/U}) = \text{Katal/L Lipase}$.

REPORTABLE RANGE

Limit of Detection

Limit of detection of the Lipase Color assay, quantified as 2 S.D. plus the mean of twenty replicate measurements of saline, is 2 U/L on a Roche/Hitachi 911 analyzer.

Linearité

The Lipase Color method is linear to 750 U/L lipase on the Roche/Hitachi 911 analyzer. If the result is greater than 750 U/L, dilute with saline, multiply the result by the dilution factor to obtain the lipase activity of the sample.

PRECISION STUDIES

Within-run precision of the Lipase Color assay was determined using three frozen spiked human lipase pools. Each run consisted of 20 replicate samples. Within-run precision studies produced the following results on the Roche/Hitachi 911 Analyzer.

Serum Pool	Low	Mid	High
n	20	20	20
Mean (U/L)	33	118	269
S.D. (U/L)	0.8	1.5	2.1
CV (%)	2.4	1.2	0.8

Between run precision of the Lipase Color assay was determined using three frozen spiked human lipase pools tested in duplicate, once or twice per day, for 15 days.

Serum Pool	Low	Mid	High
n	20	20	20
Mean (U/L)	34	120	275
S.D. (U/L)	1.5	2.7	6.3
CV (%)	4.4	2.3	2.3

ACCURACY

Accuracy of the Lipase Color method (y), performed using the Roche/Hitachi 911 Analyzer was verified by comparison to a similar method (x) on the Beckman CX analyzer. Thirty-eight patient serum samples from 50 – 572 U/L gave a correlation coefficient of 0.9988. Deming regression analysis gave the following equation:

$$\text{This method} = 1.045 (\text{reference method}) - 0.7 \text{ U/L}$$

The word SEKURE and the Sekure logo are trademarks of SEKISUI Diagnostics, LLC.

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IVD

In vitro diagnostic medical device
Appareil médical de diagnostic *in vitro*

Use by Date
YYYY-MM-DD or YYYY-MM
Date limite d'utilisation
AAAA-MM-JJ ou AAAA-MM

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1. Imamura S., Misaki H., "A sensitive method for assay of lipase activity by coupling with β -oxidation enzymes of fatty acid." Selected Topics in Clinical Enzymology; 2:73 (1984).
2. Imamura S., et al, "A sensitive method for assay of Lipase activity using 1,2-diglyceride as substrate and coupling with β -oxidation enzymes of fatty acid as an indicator reaction." Collection of summaries of lectures in the 126th general meeting of Kinki branch, analytical section, Japan Society of Clinical Chemistry; p.11-31 (1986).
3. Hayashi C., et al, "Assay methods for human lipases." Clinical Examination, Instrument and Reagent, 2:225 (1986).
4. Imamura S., et al, "An Enzymatic Method Using 1,2-Diglyceride for Pancreatic Lipase Test in Serum." Clin. Chem. 1989; 35 (6): 1126.
5. Tietz NW, "Clinical Guide To Laboratory Tests", 3rd ed., Philadelphia, PA; WB Saunders Co.; 364 (1995).
6. Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 1999.
7. Clinical and Laboratory Standards Institute, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline - Fourth Edition, GP40-A4-AMD
8. National Committee for Clinical Laboratory Standards, "Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture", Approved Standard, Third Edition, NCCLS publication H4-A3, Villanova, PA (1991).
9. Young DS, Effects of Drugs on Clinical Laboratory Tests, 4th Edition, AACC Press, Washington, DC; 3-398 to 3-400 (1995).
10. National Committee for Clinical Laboratory Standards, "How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline", NCCLS Document C28-A. Vol. 15, No.4, June 1995.

©2024 SEKISUI Diagnostics, LLC - All rights reserved.

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techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK

Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com

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©2024 SEKISUI Diagnostics, LLC - Tous droits réservés.
SEKISUI Diagnostics P.E.I. Inc.