EN

IRON-SL ASSAY

CATALOGUE NUMBER: 157-30 SIZE: R1: 3 x 100 mL
 R2: 1 x 75 mL
 Calibrator: 1 x 10 mL

Note: Changes are highlighted.

INTENDED USE


For the *in vitro* quantitative measurement of iron in serum.

TEST SUMMARY

Iron measurements are used in the diagnosis and treatment of diseases and conditions such as iron deficiency anemia, hemochromatosis (characterized by a progressive increase in iron stores leading to organ impairment), chronic inflammatory disorders, hepatitis, and lead poisoning⁽¹⁾.

Various photometric methods have been used to measure serum iron. In 1970 Stookey⁽²⁾ reported the synthesis of 3-(2-pyridyl)-5,6-bis(4-phenylsulfonic acid)-1,2,4-triazine, monosodium salt (Ferrozine[®]) which complexed with ferrous iron to form a tris ferrozine/iron, Fe(Fz)₃ complex. An additional ferrozine type compound called 5,5'-(3-(2-pyridyl)-1,2,4-triazine-5,6 diyl)-bis-2-furansulfonic acid, disodium salt (Ferene[®]) has become available^(3,4,5) and is used in this reagent. Ferene[®] is a superior iron chelating agent forming a tris complex with ferrous iron with a maximum absorption at 593 nm and a molar absorptivity of 35,500. The compound has a 27% higher molar absorption than ferrozine, absorbs at a longer wavelength and has the other advantages of ferrozine; namely, its solubility and stability over the pH range 4-9.

TEST PRINCIPLE

In an acidic medium transferrin bound iron dissociates into ferric ions which are reduced to ferrous ions in the presence of ascorbic acid. The ferrous iron reacts with the chromogen Ferene[®] to form a blue chromophore which absorbs at 595 nm. The absorbance is directly proportional to the serum iron concentration.

REAGENTS

Iron-SL Acid Dissociating Reagent (R1): A buffer solution (pH 4.5 at 25°C) containing a surfactant, preservatives, and stabilizers.

Iron-SL Color Reagent (R2): A solution containing 5.4 mmol/L ferene and stabilizers.

Iron Calibrator (for manual assay only): 1 x 10 mL of a solution containing 195 µg/dL (35 µmol/L) iron (included with 157-10 and 157-30 only).

WARNINGS AND PRECAUTIONS FOR USE

IVD

For *in vitro* diagnostic use

R₁ ONLY

IRON-SL Acid Dissociating Reagent (R1)

Warning

Contains: Acetic acid (CAS No) 64-19-7, thiourea, thiocarbamide (CAS No) 62-56-6

Hazard statements

H315 - Causes skin irritation
 H319 - Causes serious eye irritation
 H351 - Suspected of causing cancer
 H361 - Suspected of damaging fertility or the unborn child

Precautionary statements

P201 - Obtain special instructions before use.
 P202 - Do not handle until all safety precautions have been read and understood.
 P264 - Wash hands, forearms and face thoroughly after handling.
 P280 - Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352 - If on skin: Wash with plenty of water.
 P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P308+P313 - If exposed or concerned: Get medical advice/attention.
 P321 - Specific treatment (see supplemental first aid instruction on this label).
 P332+P313 - If skin irritation occurs: Get medical advice/attention.
 P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

P405 - Store locked up.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

IRON-SL Color Reagent (R2)

Warning

Contains: thiourea, thiocarbamide (CAS No) 62-56-6

Hazard statements

H351 - Suspected of causing cancer
 H361 - Suspected of damaging fertility or the unborn child

Precautionary statements

P201 - Obtain special instructions before use.
 P202 - Do not handle until all safety precautions have been read and understood.
 P280 - Wear protective gloves/protective clothing/eye protection/face protection.
 P308+P313 - If exposed or concerned: Get medical advice/attention.
 P405 - Store locked up.
 P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Iron Calibrator

Contains: Formaldehyde (50-00-0). May produce an allergic reaction.

Avoid contact with skin and eyes.

See Safety Data Sheet for additional information.

REAGENT PREPARATION, STORAGE AND STABILITY

Reagents and calibrator are ready to use.

The reagents and calibrator included are stable until the expiry date stated on the labels at 2-8°C. The R2 Color Reagent should be protected from light.

Stability claims are based on real time studies.

REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations. Reagent 1 contains Triton detergent. Dispose of waste and unused Reagent 1 by incineration.

SPECIMEN

Freshly drawn, clear, unhemolysed serum from fasting patients is the specimen of choice. Avoid anticoagulants. Serum should be separated from the cells promptly to minimize hemolysis. It is recommended that specimen collection be carried out in accordance with CLSI document H4.⁽¹¹⁾

GLASSWARE PREPARATION

All glassware and equipment used in an iron assay must be free of contaminating iron. Glassware may be prepared by soaking overnight in 1 N HCl or sulfuric acid-dichromate cleaning solution. If stronger concentrations of HCl are used, the time necessary for decontamination may be decreased. The glassware should be rinsed with deionized water before it is used.

SAMPLE STORAGE

Serum iron is reported to be stable for 4 days at 18-26°C or 7 days at 2-8°C.⁽⁶⁾

ANALYTICAL SPECIFICITY (CLSI EP7)⁽¹¹⁾

Cross Contamination studies have not been performed on automated instruments. Certain reagent / instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Copper is the only cation of the trace metals normally present in serum capable of forming a colored complex with ferene. Copper interference with ferene is similar to that encountered with ferrozine and studied by Duffy and Gaudin⁽⁷⁾. Ninety-five percent of the copper interference is eliminated by chelation of free copper.

Sera obtained from patients receiving heparin for therapy (for example, Hemodialysis) may recover falsely high results with Iron-SL reagent.⁽⁹⁾

Interferences from icterus, lipemia and hemolysis were evaluated for this iron method on a Roche/Hitachi[®] analyzer using a significance criterion of >10% variance from control. Interference data was collected in serum.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
73 µg/dL	13.1 µmol/L	Hemoglobin	100 mg/dL	15.5 µmol/L
75 µg/dL	13.4 µmol/L	Bilirubin	40 mg/dL	684 µmol/L
69 µg/dL	12.4 µmol/L	Intralipid	1000 mg/dL	3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

Samples containing elevated levels of Immunoglobulin M (IgM) or samples from patients with Waldenstrom's Macroglobulinemia may produce unreliable results.

Samples containing the following should not be used: deferoxamine.

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.⁽⁸⁾

ANALYTICAL PROCEDURE

MATERIALS PROVIDED

SEKISUI Diagnostics' Iron-SL reagents.

MATERIALS REQUIRED

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelengths as per instrument application.
2. Calibration material.
3. Quality Control materials.

TEST CONDITIONS

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:15:3 and a wavelength reading of 600/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration using an automated system is dependent on the system and the parameters used. A serum based calibrator is required for automated systems.

QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required. The results should fall within the acceptable range as established by the laboratory. **NOT INTENDED FOR USE WITH BOVINE BASED CONTROL SERUM.**⁽¹⁰⁾

CALCULATIONS

The analyzer automatically calculates the iron concentration of each sample.

TEST LIMITATIONS

A sample with an iron concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

REFERENCE INTERVALS⁽¹¹⁾

Male: 65-170 µg/dL (11.6-30.4 µmol/L)
Female: 50-170 µg/dL (8.9-30.4 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located. Serum iron concentrations show a diurnal variation with peak values seen in the early morning.

PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® analyzer unless otherwise stated.

RESULTS

Iron concentration is reported as µg/dL (µmol/L).

REPORTABLE RANGE (CLSI EP6)⁽¹²⁾

The linearity of the procedure described is 1000 µg/dL (179 µmol/L). The lower limit of detection of the procedure described is 8 µg/dL (2 µmol/L). This data results in a reportable range of 8-1000 µg/dL (2-179 µmol/L).

PRECISION STUDIES (CLSI EP5)⁽¹²⁾

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration		Total SD		Total CV %	Within Run SD		Within Run CV %
µg/dL	µmol/L	µg/dL	µmol/L		µg/dL	µmol/L	
84	15	1.4	0.25	1.6	0.7	0.13	0.8
277	50	3.0	0.54	1.1	1.0	0.18	0.3

ACCURACY (CLSI EP9)⁽¹²⁾

The performance of this method (y) was compared with the performance of a similar commercially available method (x) on a Roche/Hitachi® analyzer. Forty-five patient serum samples ranging from 9-288 µg/dL (2-52 µmol/L) gave a correlation coefficient of 0.9995. Linear regression analysis gave the following equation:
This method = 1.00 (reference method) - 5 µg/dL (1 µmol/L).

The information presented above is based on results from SEKISUI Diagnostics studies and is current at the date of publication.

TRADEMARKS

All other trademarks are the property of their respective owners.

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REFERENCES / BIBLIOGRAFÍA

- Burtis, Carl A., Ashwood, Edward R. (Ed), Tietz Textbook of Clinical Chemistry, second edition, W.B. Saunders Company, Philadelphia, London, Toronto, Montreal, Sydney, Tokyo, p. 2062, 1994.
- Stookey, L.L., Ferrozine - A New Spectrophotometric Reagent for Iron, Anal. Chem. 42, 779 (1970).
- Artiss, J.D., Vinogradov, S., Zak, B., Spectrophotometric Study of Several Sensitive Reagents for Serum Iron, Clin. Biochem. 14, 311-315 (1981).
- Higgins, T., Novel Chromogen for Serum Iron Determinations, Clin. Chem. 27, 1619 (1981).
- Artiss, J.D., Strandbergh, D.R., Zak, B., Study of Continuous Flow Automation for Serum Iron on Comparing Several Sensitive Reagents. Microchemical Journal, 28, 275-284 (1983).
- Weissman, N., Pileggi, V.J, In Clinical Chemistry Principles and Technics, 2a ed., R.J. Henry, D.C. Cannon, J.W. Winkelman, Editors, Harper & Row, Hagerstown (MD), 1974, pg. 684, 685, 695.
- Duffy, J.R., Gaudin, J., Copper Interference in the Determination of Iron in Serum Using Ferrozine, Clin. Biochem. 10, 122-123 (1977).
- Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, third edition, Washington, 1990.
- Bunting, Peter S., Aggarwal, Mahesh, Interference from Renal Dialysis Patients' Specimens in a Direct Method for Serum Iron, Clin. Chem. 29/6, 1106-1108 (1983).
- Barbour, H.M., Logan S., Interaction between bovine control materials and the ferene iron method, Annals of Clin. Biochem. 1993, 221-222 (1993).
- Burtis, Carl A., Ashwood, Edward R. (Ed), Tietz Textbook of Clinical Chemistry, second edition, W.B. Saunders Company, Philadelphia, London, Toronto, Montreal, Sydney, Tokyo, p. 2195, 1994.
- CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.

The word SEKURE and the Sekure logo are registered trademarks of SEKISUI Diagnostics, LLC.

IN15710-21
September 26, 2022

The Americas

SEKISUI Diagnostics P.E.I. Inc.
70 Watts Avenue
Charlottetown, PE C1E 2B9
Canada

Phone: 800-565-0265
Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com
techservices@sekisuidiagnostics.com

International

SEKISUI Diagnostics (UK) Limited
Liphook Way
Allington, Maidstone
KENT, ME16 0LQ, UK
Email: info@sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

sekisuidiagnostics.com