EN

## MICROALBUMIN ASSAY

**CATALOGUE NUMBER:** 252-20 **SIZE:** R1: 4 x 20 mL; R2: 2 x 10 mL

### INTENDED USE

For the IN VITRO quantitative measurement of low concentrations of albumin in urine.

### TEST SUMMARY

Low concentrations of protein are normally excreted in the urine of healthy individuals. The uriniferous tubules and glomeruli filter out most of these excreted mucoproteins. Albumin, a protein of molecular weight of 50,000, is not easily filtered out and small amounts are excreted into the urine. Increased excretion of albumin (microalbuminuria) is an early indicator of glomerular disease<sup>(1,2)</sup>.

Microalbuminuria is characterized by increased urinary excretion of albumin in the absence of overt nephropathy<sup>(3,4)</sup>. Microalbumin is recognized as a strong predictor of impending nephropathy in Type I Diabetics and its mortality risk in diabetic patients<sup>(5)</sup>. Early detection of microalbuminuria may be beneficial for treatment programs for diabetics because renal damage may be reversible if diabetes is well controlled at this stage.

Many of the methods traditionally used for measuring albumin lack the sensitivity and precision required for measuring microalbumin. Sekisui Diagnostics' Microalbumin method uses an immunoturbidimetric format which provides the sensitivity required for accurate measurement of urinary microalbumin.

### TEST PRINCIPLE

When a sample is mixed with anti-human albumin goat antiserum, agglutination is caused by the antigen-antibody reaction. The turbidity is measured at 340 nm and 700 nm and albumin in the sample is quantitatively determined.

### REAGENTS

**Microalbumin Buffer Reagent (R1):** A solution containing 100 mM tris (hydroxymethyl) aminomethane (pH 7.6 at 25°C) and a preservative.

**Microalbumin Antibody Reagent (R2):** A solution containing 20% anti-human albumin, goat anti serum, 100 mM tris (hydroxymethyl) aminomethane, and a preservative.

### WARNINGS AND PRECAUTIONS

**IVD**

For *In Vitro* Diagnostic Use.

**Rx ONLY**

Avoid contact with skin and eyes. See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

The reagents are provided in a ready to use format. The reagents included are stable until the expiry date stated on the labels when stored at 2-8°C.

Opened reagents can be used for one month if stored at 2-8°C. Recap and return reagents to 2-8°C promptly after use.

### REAGENT DETERIORATION

The reagent solutions should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations.

### SPECIMEN

The specimen should be a fresh or 24 hour urine.

### SAMPLE STORAGE

Urine specimens should be stored at 2-8°C. The specimens may be stored refrigerated up to a week or frozen at -20°C. USE PLASTIC TUBES FOR STORING THE SAMPLES. DO NOT USE GLASS.

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(6)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/ instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interferences from icterus, lipemia and hemolysis were evaluated for this method on a Roche/Hitachi® 717 analyzer using a significance criterion of >10% from control.

Lipemia produces significant interference with this method; therefore, lipemic samples should not be used.

Concentration of Analyte		Substance Tested	Concentration of Interferent Where Interference is Insignificant	
Conventional Units	SI Units			
17.1 mg/L	1.7 mg/dL	Hemoglobin	400 mg/dL	62 µmol/L
15.6 mg/L	1.6 mg/dL	Bilirubin	4 mg/dL	68.4 µmol/L

The information presented above is based on results from Sekisui Diagnostics' studies and is current at the date of publication.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(7)</sup>

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

Sekisui Diagnostics' Microalbumin reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at appropriate wavelength as per instrument application.
2. Calibration material.
3. Quality Control materials.

#### TEST CONDITION

For data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:40 and wavelength readings of (primary/secondary) 340/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact Sekisui Diagnostics Technical Services at (800)565-0265. Outside Canada and the U.S., please contact your local distributor.

#### CALIBRATION

For albumin concentrations less than 100 mg/L a two point calibration curve can be made using a saline blank (0 mg/dL) and an albumin standard. When a sample with a high albumin concentration (> 100 mg/L) is assayed, it is recommended that a multipoint calibration curve be made. It is recommended that Sekisui Diagnostics' Microalbumin Multi-Calibrator Set, catalogue number SE-252 be used in conjunction with this assay. The frequency of calibration on automated systems is dependent on the system and the parameters used.

#### QUALITY CONTROL

A normal and abnormal concentration control serum should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

#### CALCULATIONS

The analyzer calculates the albumin concentration of each sample.

#### TEST LIMITATIONS

A sample with an albumin concentration exceeding the linearity limit should be diluted one part sample with four parts isotonic saline including 0.05% Tween 20 and reassayed, incorporating the dilution factor of five into the calculation of the value.

## REFERENCE INTERVALS<sup>(8)</sup>

American Diabetes Association Definition of Moderately increased urine albumin.

Urine	24 h Collection (mg/24 h)	Timed Collection (µg/min)	Spot collection (µg/mg creatinine)
Normal	< 30	< 20	< 30
Moderately increased	30-299	20-199	30-299
Clinical Albuminuria	≥ 300	≥ 200	≥ 300

These values are suggested guidelines. It is recommended that each laboratory establish the normal range for the area in which it is located.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 717 analyzer unless otherwise stated.

## RESULTS

Albumin concentration is expressed as mg/L (mg/dL).

## REPORTABLE RANGE (CLSI EP6)<sup>(6)</sup>

The linearity of the procedure described is 300 mg/L (30.0 mg/dL). The lower limit of detection of the procedure described is 5 mg/L (0.5 mg/dL). This data results in a reportable range of 5-300 mg/L (0.5-30.0 mg/dL).

## PRECISION STUDIES (CLSI EP5)<sup>(6)</sup>

Total precision data was collected on two concentrations of control urine in five runs over three days.

Concentration		Total SD		Total CV (%)	Concentration		Within Run SD		Within Run CV (%)
mg/L	mg/dL	mg/L	mg/dL		mg/L	mg/dL	mg/L	mg/dL	
7.9	0.79	0.5	0.05	6.0	5.8	0.58	0.3	0.03	4.8
54.9	5.49	2.0	0.20	3.6	45.0	4.50	1.8	0.18	3.9

Within run precision data was collected on two concentrations of control urine each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(6)</sup>

The performance of this method (y) was compared with the performance of a similar method (x) on a Roche/Hitachi® analyzer. Sixty-four patient urine samples ranging from 0-182 mg/L (0.0-18.2 mg/dL) were tested and gave a correlation coefficient of 0.9981. Linear regression analysis gave the following equation:

$$\text{This method} = 0.949 (\text{reference method}) + 0.4 \text{ mg/L (0.04 mg/dL)}$$

## TRADEMARKS

All trademarks, brands, product names and trade names are the property of their respective companies.

Manufactured by:

**SEKISUI**  
DIAGNOSTICS

### The Americas

Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada

Phone: 800-565-0265

Fax: 902-628-6504

Email: questions@sekisuidiagnostics.com

peidiagnosticttechnical@sekisuidiagnostics.com

www.sekisuidiagnostics.com

### International

Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK

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The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.
MDSS GmbH