## N-geneous® LDL Cholesterol Calibrator For In Vitro Diagnostics Use

CATALOGUE NUMBER: 80-4610-02 SIZE: 3 x 1 mL

**Note:** Changes are highlighted.

### INTENDED USE

For calibration of the N-geneous® LDL Cholesterol and LDL Direct Liquid Select Cholesterol assays in serum or plasma.

### TEST SUMMARY

See the N-geneous® LDL Cholesterol and LDL Direct Liquid Select Cholesterol Reagent package inserts.

### REAGENTS

#### Composition

The N-geneous® LDL Cholesterol Calibrator is a preparation of lyophilized human serum containing lipoproteins from the various lipoprotein classes including low-density lipoproteins.

**NOTE:** The LDL cholesterol value is traceable to the Reference Method (Ultracentrifugation) for the determination of LDL Cholesterol.

### WARNINGS AND PRECAUTIONS FOR USE

Dispose of contents/container in accordance with local/regional /national/international regulations.

- IVD**  
For *in vitro* diagnostic use.
- FOR ONLY**
- Do not pipette by mouth.
- WARNING:** Human source material. Treat as potentially infectious. Each plasma donor unit used in the preparation of this product has been tested by an FDA approved method and found non-reactive for HBsAg, HCV, HIV 1&2 and HIV-1 Antigen. Because no known test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus (HIV) or other infectious agents are absent, all human-based products should be handled in accordance with good laboratory practices using appropriate precautions.
- Do not use the calibrator after the expiration date printed on the box.
- Contains 0.1% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush drain with large volume of water to prevent buildup.

See Safety Data Sheet for additional information.

### PRODUCT PREPARATION, STORAGE AND STABILITY

Lyophilized serum calibrator should be reconstituted by adding 1.0 mL of distilled or deionized water. Close the vial and let stand for 5 minutes. Dissolve the contents of the vial by swirling gently, avoiding the formation of foam. **Do not shake.**

Unopened calibrator is stable at 2-8°C until the expiration date printed on the vial.

After reconstitution, calibrator is stable up to 2 weeks at 2-8°C.

The reconstituted calibrator may be aliquoted and stored at -80°C for 1 month.

### PRODUCT DETERIORATION

Presence of extreme turbidity or growth may indicate deterioration.

### DISPOSAL

Product must be disposed of in accordance with all Federal, Provincial, State and local regulations.

### ANALYTICAL PROCEDURE

### MATERIALS PROVIDED

Description	Configuration	Catalog Number
N-geneous® LDL Cholesterol Calibrator	3 x 1 mL	80-4610-02

### MATERIALS REQUIRED (BUT NOT PROVIDED)

- N-geneous® LDL Cholesterol Reagent.

N-geneous® LDL Cholesterol Reagent.	Catalog Number
Reagent 1, 250 mL fill	80-4598-00
Reagent 2, 80 mL fill	80-4601-00
Reagent 1, 2 Liter fill	80-4604-00
Reagent 2, 2 Liter fill	80-4607-00

LDL Direct Liquid Select Cholesterol Reagent.	Catalog Number
LDL Direct Liquid Select Cholesterol	7120

- Volumetric pipette (1.0 mL).
- Distilled or deionized water meeting specifications equivalent to USP purified water.

### CALIBRATION

The calibrator should be run along side patient samples in accordance with the instruction outlined in Sekisui Diagnostics N-geneous® LDL Cholesterol Reagent and LDL Direct Liquid Select Cholesterol Reagent package inserts. The values of the N-geneous® LDL Cholesterol Calibrator, found on the vial label, were assigned by procedures traceable to the CDC LDL cholesterol reference method.<sup>9</sup> Calibration materials have concentrations around the medical decision level. Refer to the instrument manufacturer's recommendation for calibration frequency. For additional assistance please contact Sekisui Diagnostics Technical Services at 800-565-0265.

### REFERENCE INTERVALS

The following cutpoints for patient classification are used for the prevention and management of coronary heart disease.

### LDL Cholesterol

<130 mg/dL  
(<3.36 mmol/L)

130-159 mg/dL  
(3.36-4.11 mmol/L)

≥160 mg/dL  
(≥4.14 mmol/L)

### Classification

Desirable

Borderline High Risk

High Risk

All trademarks, brands, product names and trade names are the property of their respective companies.

**SEKISUI**  
DIAGNOSTICS

#### The Americas

Sekisui Diagnostics P.E.I. Inc.  
70 Watts Avenue  
Charlottetown, PE C1E 2B9  
Canada  
Phone: 800-565-0265  
Fax: 902-628-6504  
Email: questions@sekisuidiagnostics.com  
peidiagnosticttechnical@sekisuidiagnostics.com

#### International

Sekisui Diagnostics (UK) Limited  
Liphook Way  
Allington, Maidstone  
KENT, ME16 0LQ, UK  
Email: info@sekisuidiagnostics.com

www.sekisuidiagnostics.com

FR

## N-geneous® LDL Cholesterol Calibrator

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*In vitro* diagnostic medical device  
Appareil médical de diagnostic *in vitro*.

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1. Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia, *Hospital Practice*, 23:Suppl. 1,4 (1988).
2. Crouse, J.R. et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, *Lipid Res.*, 26;566 (1985).
3. Kannel, W.B., et al., Cholesterol in the prediction of atherosclerotic disease. New perspectives base on the Framingham study, *Ann Intern. Med.*, 90-85-91, (1979).
4. Kannel, W.B., et al., Optimal resources for primary prevention of atherosclerotic disease, *Circulation*, 70:157A – 205A, (1984).
5. Schaefer, E.J., et al., Pathogenesis and management of lipoprotein disorders., *N. Engl. J. Med.*, 312(20) 1300-1310 (1985).
6. The Lipid research Clinics, coronary primary prevention trial results., *JAMA*, 251(3) 351-374.
7. Consensus Development conference: Lowering Blood cholesterol to prevent coronary heart disease, *JAMA*, 253(14), 2080-2086 (1985).
8. Report of the National Cholesterol Education Program Expert panel on Detection, Evaluation and Treatment of High Blood Cholesterol in adults, *Arch. Intern. Med.*, 148:36-69 (1988).
9. National Reference System for Cholesterol. CRMLN LDL Cholesterol Protocol, May 2004.

N-geneous® is a registered trademark of Sekisui Diagnostics, LLC.  
©2011 Sekisui Diagnostics, LLC - All rights reserved.

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The word SEKURE and the Sekure logo are trademarks of Sekisui Diagnostics, LLC.

IN80461002-3  
May 7, 2020

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N-geneous® est une marque de commerce déposée de Sekisui Diagnostics, LLC.
©2011 Sekisui Diagnostics, LLC - Tous droits réservés.