

OSOM[®] RSV/Adeno Test 197E

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

EN **OSOM® RSV/Adeno Test** Catalog number 197E

**FOR EXPORT USE ONLY. NOT FOR SALE IN THE U.S.
FOR LABORATORY AND PROFESSIONAL USE ONLY.**

INTENDED USE

The OSOM® RSV/Adeno test is a rapid chromatographic immunoassay for the qualitative detection of Respiratory Syncytial Virus (RSV) and/or Adenovirus antigens directly from nasal swabs or nasal suction fluid in patients suspected of having a viral respiratory infection. Throat swabs are also an acceptable sample type for the detection of Adenovirus antigens. This test is intended for use as an aid in the diagnosis of RSV and/or Adenovirus infections in symptomatic patients. It is recommended that negative test results be confirmed by viral culture. Negative results should not be used as the sole basis for treatment.

SUMMARY AND EXPLANATION OF TEST

RSV is an RNA virus that is responsible for outbreaks of respiratory tract infections. RSV infections can occur throughout the year, but typically peak during the winter months.¹ It is said that approximately 50% and nearly 100% of children suffer from infection with these viruses by 1 and 2 years of age, respectively. RSV viruses not only cause upper respiratory tract infections but also bronchiolitis of the lower respiratory tract, which often becomes severe in infants and toddlers with underlying diseases.² Children who were born premature, or who have pre-existing lung, heart or immune dysfunctions have the greatest risk of developing RSV associated infections. Diagnosis of RSV is difficult because the initial symptoms can be similar to those caused by other infectious agents. Considering that the

RSV virus is highly contagious, accurate diagnosis and prompt treatment of patients can have a positive effect on public health.

Adenovirus causes respiratory tract infections, such as pharyngitis/tonsillitis, pharyngoconjunctival fever, and pneumonia, but also diarrhea, epidemic keratoconjunctivitis, and other diseases. Pharyngoconjunctival fever is epidemic in summer, while other adenovirus infections are nearly perennial. The incidence of adenovirus respiratory tract infections is higher in children. These viruses account for up to 10% of such infections in children, and the infections sometimes become severe.^{3,4,5}

Accurate diagnosis can also help reduce the inappropriate use of antibiotics and gives the physician the opportunity to prescribe an appropriate therapy. The OSOM RSV/Adeno Test can provide rapid detection of RSV and/or Adenovirus antigens from symptomatic patients.

PRINCIPLE OF THE TEST

The OSOM RSV/Adeno Test consists of a test stick that specifically detects RSV and/or Adenovirus antigens. The test procedure requires the solubilization of the antigen from a nasal swab or wash/aspirate sample, by mixing with Extraction Buffer. For adenovirus, throat swabs are also acceptable. The test stick is then placed in the sample mixture, which migrates along the membrane surface.

If RSV antigens are present in the sample, it will form a complex with gold colloidal-labeled mouse monoclonal IgG antibodies to RSV. The resulting complex will then be bound by another mouse anti-RSV antibody coated on the nitrocellulose membrane to form a red/purple line in the presence of RSV antigen.

If adenovirus antigens are present in the sample, it will form a complex with gold colloidal-labeled mouse monoclonal IgG antibodies to adenovirus. The resulting complex will then be bound by another mouse anti-adenovirus antibody coated on the nitrocellulose membrane to form a red/purple line in the presence of adenovirus antigen.

A black control line must appear in the control region of the stick for results to be valid. The appearance of red/purple lines of any intensity in the test line region indicates a positive result

REAGENTS AND MATERIALS PROVIDED

- 25 Test Sticks
- 25 Test Tubes
- 25 Sterile Nasal Swabs
- 1 Extraction Buffer vial
 - 12 mL – phosphate buffered salt solution (pH 7.6), and 0.09% sodium azide (as a preservative)
- 1 Extraction buffer dropper top
- 1 Directional Insert
- 1 Workstation

MATERIALS REQUIRED BUT NOT PROVIDED

- A timer or watch
- Throat swabs (for adenovirus detection only) Part # 7703
- Nasal suction collection container

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- Only use throat swabs to detect adenovirus. Detection of RSV in throat swabs is not an intended use of this product.
- Follow your clinical and/or laboratory safety guidelines in the collection, handling, storage and disposal of patient specimens and all items exposed to patient specimens.⁶
- Swabs, tubes and test sticks are for single use only.
- The extraction buffer contains a solution with a preservative (0.09% sodium azide). If solution comes in contact with the skin or eyes, flush with ample volumes of water.
- The Sample Buffer contains 0.1-1.0% Triton detergent. Dispose of waste and unused Sample Buffer by incineration.
- Do not interchange or mix components from different kit lots.

STORAGE CONDITIONS


- Store Test Sticks and Extraction Buffer tightly capped at room temperature (15°-30°C/59°-86°F).
- Do not freeze any of the test kit components.
- Do not use test sticks and reagents after expiration date.
- Recap the desiccated container immediately after removing a test stick.

SPECIMEN COLLECTION AND PREPARATION

Methods of sample collection


1) Nasal swabs

Imagine a line connecting the nasal cavity and the external auditory meatus (EAM), and insert a swab for the collection of nasal samples (included in kit) along it. When the swab reaches the base of the nasal cavity and stops going forward, collect a sample from the nasal membrane by rubbing it several times.


2) Nasal suction fluid

Attach a suction catheter with mucus trap to the suction apparatus. Insert one end of the catheter into the nasal cavity, and suck out the nasal fluid. Collect a sample by immersing a swab included with kit into the nasal fluid collected in the trap.

3) Throat swabs (Adenovirus only)

While the mouth is open wide, depress the tongue with a tongue depressor. Collect a sample carefully by strongly rubbing the palate, tonsil, and posterior wall of the pharynx with a separately sold swab (for the collection of throat samples). Caution should be exercised to prevent contact between the swab and the inside of the buccal cavity/tongue/teeth.

PRECAUTIONS

- Use only the nasal swab or nasal suction method for detection of RSV antigens.
- Throat swab specimens may be used for the detection of adenovirus antigens. These are not provided with the kit but can be ordered separately. (Part # 7703)
- The swabs provided with the kit are intended for use as nasal swabs only. Do not use these swabs to collect throat samples.
- Test the sample as soon as possible after collecting the specimen. If samples cannot be processed immediately, specimens may be held at 2°-8°C (36°-46°F) for up to 24 hours.
- To transport patient nasal swabs, place swab in a clean, dry container such as a plastic or glass tube.
- **A separate swab must be collected if culture is desired.**
- The test performance depends on the quality of the sample obtained, and the handling and transport of the sample. Training in specimen collection is recommended because of the importance of specimen quality. False negative results may occur from inadequate specimen collection and/or handling.
- If a sample is too viscous, the pretreated sample may not develop a line on the test stick. In this case, classify the result of the test as invalid or collect another sample and repeat the test with a different test stick.

QUALITY CONTROL (QC)

The OSOM RSV/Adeno test provides two types of controls: Procedural internal controls to aid in determining test validity and external controls.


Internal Procedural Controls

Several controls are incorporated into each test stick for routine quality checks.

1. The appearance of the black control line in the results window is an internal procedural control:
 - **Test System:** The appearance of the black control line assures that adequate test volume was present and that adequate capillary migration of the extracted sample has occurred. It also verifies proper assembly of the Test Stick.

- **Operator:** The appearance of the control line indicates that adequate test volume was present for capillary flow to occur. If the control line does not appear at the end of read time, the test is invalid.

2. The clearing of the background in the results area may also be documented as an internal procedural control. It also serves as an additional capillary flow control. At the read time, the background should appear white to light pink and not interfere with the reading of the test. If the background color does not clear and interferes with the test result, the test is invalid.

TEST PROCEDURE, NASAL OR THROAT SWAB

STEP 1: ADD EXTRACTION BUFFER

Using the supplied dropper top, add 0.3 mL of Extraction buffer to each test tube. Fill the dropper to the line indicated on the barrel of the dropper top and expel entire contents into tube.

Note: Add Extraction buffer to the test tube before putting in the specimen swab to prevent contaminating the Extraction Buffer vial.

STEP 2: MIX SWAB IN BUFFER

Put the specimen swab into the tube. Vigorously mix the solution by rotating the swab forcefully against the side of the tube at least ten times (while submerged). Best results are obtained when the specimen is vigorously mixed in the solution.

STEP 3: SQUEEZE LIQUID FROM SWAB


Squeeze out as much liquid as possible from the swab by pinching the side of the flexible test tube as the swab is removed. Discard the swab in a suitable biohazardous waste container.

STEP 4: ADD TEST STICK

Remove a Test Stick from the canister. Recap the canister immediately. Place the test stick (arrows pointing downward) into the tube with the extraction buffer solution. Set a timer for 10 minutes.

STEP 5: READ RESULTS

At 10 minutes read results by removing stick from the tube and holding against a white background. Some positive results may be seen and reported earlier. See interpretation of results section. Test is invalid beyond the stated read time. Discard used test tubes and Test Sticks in suitable biohazardous waste container.

TEST PROCEDURE, NASAL SUCTION FLUID


STEP 1: ADD EXTRACTION BUFFER

Using the supplied dropper top, add 0.3 mL of Extraction buffer to each test tube. Fill the dropper to the line indicated on the barrel of the dropper top and expel entire contents into tube.

Note: Add Extraction buffer to the test tube before putting in the specimen swab to prevent contaminating the Extraction Buffer vial.

STEP 2: COLLECT SAMPLE

To obtain sample from suction fluid, immerse a swab, included with kit, into the nasal fluid collected in the trap.

STEP 3: MIX SWAB IN BUFFER

Put the specimen swab into the tube. Vigorously mix the solution by rotating the swab forcefully against the side of the tube at least ten times (while submerged). Best results are obtained when the specimen is vigorously mixed in the solution.

STEP 4: SQUEEZE LIQUID FROM SWAB


Squeeze out as much liquid as possible from the swab by pinching the side of the flexible test tube as the swab is removed. Discard the swab in a suitable biohazardous waste container.

STEP 5: ADD TEST STICK

Remove a Test Stick from the canister. Recap the canister immediately. Place the test stick (arrows pointing downward) into the tube with the extraction buffer solution. Set a timer for 10 minutes.

STEP 6: READ RESULTS

At 10 minutes read the results by removing stick from the tube and holding against a white background. Some positive results may be seen and reported earlier. See interpretation of results section. Test is invalid beyond the stated read time. Discard used test tubes and Test Sticks in suitable biohazardous waste container.

INTERPRETATION OF TEST RESULTS

The appearance of a black Control Line, with or without a red/purple Test Line(s), indicates a valid result. A black or red/purple line that appears uneven in color shading is still considered a valid line. In cases of moderate or high positive specimens, some color behind the Test Line may be seen. As long as the Test Line and the Control Line are visible, the results are valid.

RSV positive

One black line in the control position and one pink/purple line in the RSV test line position.

Adenovirus positive

One black line in the control position and one pink/purple line in the Adeno test line position.

Both RSV and Adenovirus positive

One black line in the control position and a pink/purple line at both the RSV and the Adeno test line positions.


Note that the red and black lines can be any shade of that color. The lines can be lighter or darker than the lines in the picture. Any visible red/purple line should be considered positive.

Negative

A black Control Line but no red/purple Test Line is a presumptive negative result. A negative result means that no RSV or adenovirus antigen was present or the level of the antigen in the sample is below the detection limit of the assay.

Invalid

If no black Control Line appears or background color makes reading the black Control Line impossible, the result is invalid. If this occurs, repeat the test using a new Test Stick.

REPORTING RESULTS

- Report negative test results when only the black control line is observed which indicates that RSV viral antigen and adenoviral antigen is not detected. Infection due to RSV or adenovirus cannot be ruled out since the antigens may be present in the specimen below the detection limit of the test. Negative tests are presumptive and should be confirmed by culture.
- Report positive test results when both the black control line and any shade of red line are observed which indicates the presence of RSV and/ or adenovirus antigen. This result does not rule out co-infections with other pathogens or identify any specific viral subtype.
- If result is considered invalid, repeat the test using a new sample and a new test dipstick.

LIMITATIONS

- The OSOM RSV/Adeno Test is for the qualitative detection of RSV viral and/or adenovirus antigens. The test performance depends on antigen load and may not correlate with viral culture performed on the same specimen. Negative test results are not intended to rule out other viral infections.
- Sensitivity can differ with various strains of RSV and/or adenovirus due to difference in antigen expression.
- This test detects both viable and non-viable RSV and/or adenovirus, and may yield a positive result in the absence of living organisms.
- Only use throat swabs to detect adenovirus. Detection of RSV in throat swabs is not an intended use of this product.
- The test performance depends on the quality of the sample obtained as well as the handling and transport of the sample. Negative results can occur from inadequate specimen collection and/or handling.
- As with all diagnostic assays, the results obtained with this test kit yield data that must be used only as an adjunct to other information available to the physician.
- High levels of blood on specimen swabs might cause an intense red background on the test strip that could interfere with the test interpretation. Avoid samples that have been heavily contaminated with whole blood.
- Positive and negative predictive values of these diagnostic assays are highly dependent on prevalence or current level of RSV and/or adenovirus activity. During peak RSV/Adenovirus activity in a season, positive predictive values are higher, with false positives less likely; and negative predictive values are lower, with false negatives more likely. Conversely, during low RSV/Adenovirus activity (e.g., off-season or beginning of a season), negative predictive values are higher and positive predictive values lower, with false positive test results more likely.

EXPECTED RESULTS

RSV/Adenovirus viruses can cause epidemics which typically occur during the winter months.

PERFORMANCE CHARACTERISTICS

Comparison of OSOM RSV/Adeno Test to In Vitro Diagnostic Agent (Immunochromatography):

NASAL SWAB OR NASAL SUCTION FLUID FOR RSV

	In Vitro Diagnostic Agent (Immunochromatography)			total
		+	-	
OSOM RSV/Adeno Test	+	178	3	181
	-	2	216	218
total		180	219	399

Clinical sensitivity: 98.9% (178/180)

Clinical specificity: 98.6% (216/219)

THROAT SWABS FOR ADENOVIRUS

	In Vitro Diagnostic Agent (Immunochromatography)			total
		+	-	
OSOM RSV/Adeno Test	+	79	5	84
	-	7	190	197
total		86	195	281

Clinical sensitivity: 91.9% (79/86)

Clinical specificity: 97.4% (190/195)

Comparison of OSOM RSV/Adeno Test to PCR:

NASAL SWAB OR NASAL SUCTION FLUID FOR RSV

	PCR			total
		+	-	
OSOM RSV/Adeno Test	+	179	2	181
	-	19	199	218
total		198	201	399

Clinical sensitivity: 90.4% (179/198)

Clinical specificity: 99.0% (199/201)

THROAT SWABS FOR ADENOVIRUS

	PCR			total
		+	-	
OSOM RSV/Adeno Test	+	75	9	84
	-	7	190	197
total		82	199	281

Clinical sensitivity: 91.5% (75/82)

Clinical specificity: 95.5% (190/199)

NASAL SWAB OR NASAL SUCTION FLUID FOR ADENOVIRUS

	PCR			total
		+	-	
OSOM RSV/Adeno Test	+	96	8	104
	-	17	500	517
total		113	508	621

Clinical sensitivity: 85.0% (96/113)
Clinical specificity: 98.4% (500/508)

Assay Reproducibility
Accuracy

- (1) The test was performed using a negative control for accuracy assessment as the sample: the test is negative.
- (2) The test was performed using a strongly positive RSV control for accuracy assessment (RS virus titer: 5.5×10^5 - 6.5×10^5 TCID₅₀/mL*) and a moderately positive RSV control for accuracy assessment (RS virus titer: 5.5×10^4 - 6.5×10^4 TCID₅₀/mL*) as the samples: the test is only positive for RS virus.
- (3) The test was performed using a strongly positive adenovirus control for accuracy assessment (adenovirus titer: 7.5×10^5 - 8.5×10^5 TCID₅₀/mL*) and a moderately positive adenovirus control for accuracy assessment (adenovirus titer: 7.5×10^4 - 8.5×10^4 TCID₅₀/mL*) as the samples: the test is only positive for adenovirus.

Within-run reproducibility

- (1) A negative control was used for within-run reproducibility assessment as the sample and repeat the test 5 times: the test is negative every time.
- (2) A strongly positive RSV control was used for within-run reproducibility assessment (RS virus titer: 5.5×10^5 - 6.5×10^5 TCID₅₀/mL*) as the sample and repeat the test 5 times: the test is only positive for RS virus every time.
- (3) A strongly positive adenovirus control was used for within-run reproducibility assessment (adenovirus titer: 7.5×10^5 - 8.5×10^5 TCID₅₀/mL*) as the sample and repeat the test 5 times: the test is only positive for adenovirus every time.

((1)-(3): tested by the in-house methods of Sekisui Medical)

*TCID₅₀/mL(50% Tissue Culture Infective Dose)
TCID₅₀/mL is defined as the virus titer in a sample measured by the TCID₅₀ method, as mentioned below. Prepare serial 10ⁿ dilutions of a sample, and add a constant amount of each dilution to 6 wells of cultured cells. Define the dilution (10ⁿ) at which a cytopathic effect (CPE) is noted in 3 out of 6 wells as the virus titer of the sample at 10ⁿ TCID₅₀/mL.

Analytical Sensitivity

RS virus antigen: 2.5×10^4 TCID₅₀/mL
Adenovirus antigen: 3.5×10^4 TCID₅₀/mL
(Tested by in-house methods of Sekisui Medical.)

Analytical Specificity and Cross-reactivity

The OSOM RSV/Adeno Test was evaluated with 48 bacterial and viral isolates. Bacterial isolates were tested at a concentration of approximately $>1 \times 10^7$ cfu/mL. All of the bacteria listed gave negative responses. Viral isolates were tested at approximately $>1 \times 10^6$ TCID₅₀/mL. All of the viruses listed produced negative responses.

Bacterial Panel

Acinetobacter calcoaceticus
Bordetella pertussis
Candida albicans
Corynebacterium diphtheriae
Enterococcus faecalis
Enterococcus gallinarum
Escherichia coli
Haemophilus influenza

Neisseria meningitidis
Proteus mirabilis
Proteus vulgaris
Pseudomonas aeruginosa
Serratia marcescens
Staphylococcus epidermidis
Streptococcus Group A
Streptococcus Group B

Klebsiella pneumoniae
Moraxella catarrhalis
Mycobacterium avium

Streptococcus mutans
Streptococcus pneumoniae
Torulopsis glabrata (*Candida glabrata*)

Viral Panel

Echovirus Type 3	Influenza virus Type A/New Caledonia/20/99 (H1N1)
Echovirus Type 6	Influenza virus Type A/Hiroshima/52/2005 (H3N2)
Echovirus Type 9	Influenza virus Type B/Shanghai/361/2002 (Yamagata lineage)
Echovirus Type 11	Influenza virus Type B/Malaysia/2506/2004 (Victoria lineage)
Echovirus Type 14	Influenza virus Type C/JJ/50
Echovirus Type 18	Influenza virus Type C/Yamagata/3/96
Echovirus Type 25	Parainfluenza virus Type 1
Echovirus Type 30	Parainfluenza virus Type 2
Coxsackievirus Type B1	Parainfluenza virus Type 3
Coxsackievirus Type B2	Mumps virus
Coxsackievirus Type B3	Rubella virus
Coxsackievirus Type B4	Measles virus
Coxsackievirus Type B5	
Herpes simplex virus Type 1	

RSV Panel testing

This kit produces a positive reaction with the following RS virus strains.

Respiratory Syncytial Virus Type A/Long
Respiratory Syncytial Virus Type A/A-2
Respiratory Syncytial Virus Type B/Washington/18537/62
Respiratory Syncytial Virus Type B/WV/14617/85[B-1 wild type]
Respiratory Syncytial Virus Type B/9320

Adenovirus Panel testing

This kit produces a positive reaction with the following adenovirus strains.

Adenovirus Type 1	Adenovirus Type 8
Adenovirus Type 2	Adenovirus Type 11
Adenovirus Type 3	Adenovirus Type 18
Adenovirus Type 4	Adenovirus Type 19
Adenovirus Type 5	Adenovirus Type 37
Adenovirus Type 6	Adenovirus Type 41
Adenovirus Type 7	

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1. Thompson, William W. PhD, Shay, David K. MD, MPH, Eric Weintraub MPH et al. Mortality Associated with Influenza and Respiratory Syncytial Virus in the United States. JAMA, January 8, 2003 (289:2) 179-186.
2. Weisman, Leonard E. MD. Populations at risk for developing respiratory syncytial virus and risk factors for respiratory syncytial virus severity: infants with predisposing conditions. Journal of Pediatric Infectious Diseases, 2003 (22:2) S33-39.
3. Rocholl, Christian, Gerber, Kris, Judy Daly et al. Adenoviral Infections in Children: The Impact of Rapid Diagnosis. Pediatrics, 2004 (113) e51-56.
4. Abd-Jamil, Juriana, Teoh, Boon-Teong, Eddy H Hassan et al. Molecular identification of adenovirus causing respiratory tract infection in pediatric patients at the University of Malaya Medical Center. BMC Pediatrics, 2010 (10:46).
5. Moon, Rachel Y. Adenovirus Infections. Pediatrics in Review. 1999 (20) 56.
6. CDC, Biosafety in Microbiological and Biomedical Laboratories, 2nd Ed., HHS Publication No. 8808395, 4-6, 1988.

DEFINITIONS OF SYMBOLS / SYMBOLFORKLARING / DÉFINITIONS DES SYMBOLES / DEFINIZIONE DEI SIMBOLI / DEFINITIONEN DER SYMBOLE / SYMBOLER / SYMBOLDEFINITIONER / OBJAŚNIENIE SYMBOLI

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In Vitro Diagnostic Medical Device

Medicinsk udstyr til *in vitro*-diagnostik
Dispositif médical de diagnostic *in vitro*
Dispositivo medico-diagnostico *in vitro*
In-vitro-Diagnostikum
Til *in vitro* diagnostisk bruk
Medicintekniska produkter för *in vitro* diagnostik
Wyrób do diagnostyki *In Vitro*

Contains sufficient for <n> tests

Indeholder tilstrækkeligt til "n" test
Contenu suffisant pour "n" tests
Contenuto sufficiente per "n" saggi
Inhalt ausreichend für <n> Prüfungen
Antall tester
Räcker till "n" antal tester
Zawiera ilość wystarczającą do wykonania <n> badań

Authorised representative in the EC

Repræsentant i det Europæiske Fællesskab
Mandataire dans la Communauté européenne
Mandatario nella Comunità Europea
Bevollmächtigter in der Europäischen Gemeinschaft
Autorisert representant
Auktoriserad representant i Europeiska Gemenskapen
Autoryzowany przedstawiciel w Unii Europejskiej

osom[®]

is a registered U.S. trademark of Sekisui Diagnostics, LLC.
er et registreret varemærke i USA tilhørende Sekisui Diagnostics, LLC
est une marque de commerce déposée aux États-Unis de Sekisui Diagnostics, LLC
è un marchio registrato negli USA di Sekisui Diagnostics, LLC
ist eine in den USA eingetragene Marke von Sekisui Diagnostics, LLC
er et registrert varemerke for Sekisui Diagnostics, LLC
är ett registrerat varumärke som tillhör Sekisui Diagnostics, LLC
to zastrzeżony amerykański znak towarowy firmy Sekisui Diagnostics, LLC

SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121
USA

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MDSS GmbH