EN

RPR REAGENT

INTENDED USE

For the quantitative measurement of syphilitic anti-lipid antibodies concentration in serum or plasma on the Hitachi 917 clinical analyzer, the SK500/Bolis 50i clinical analyzer and multiple other automated clinical analyzers.

SUMMARY

Syphilis is a chronic infection caused by *Treponema pallidum* which can be transmitted congenitally or by sexual contact. It is characterized by episodes of active disease interrupted by periods of latency¹.

There are two types of serological test for syphilis, non-treponemal and treponemal. The most widely used non-treponemal antibody tests for syphilis are the rapid plasma reagin (RPR) and venereal disease research laboratory (VDRL) tests, both measuring antibodies against a cardiolipin-lecithin-cholesterol antigen complex. Treponemal tests measure antibodies to native (Nichols strain) or recombinant *T. pallidum* antigens.¹ This automated assay is based on a latex immunological agglutination test principle.

PRINCIPLE

The polystyrene latex, coated with lipid antigens (cardiolipin and lecithin), is exposed to the test sample under certain conditions, to induce formation of anti-lipid antibody-latex aggregate. The elevation of turbidity due to the formation of this aggregate (the change in turbidity) relative to the pre-exposure level is measured, to determine the anti-lipid antibodies titer in the test sample.

REAGENTS

Composition

Component	Ingredients	Concentration
Reagent 1	Bovine serum albumin Buffer (pH 7.1 – 7.5) Sodium Azide	1% <0.1%
Reagent 2	Latex particles coated with cardiolipin, lecithin and cholesterol antigens Sodium Azide Buffer (pH 7.1 – 7.5)	≤ 4.0 mg/mL <0.1%

Precautions and Warnings

1. For In Vitro Diagnostic Use.
2. Do not use the reagents beyond the expiration date printed on the label.
3. **Warning :** All specimens used in the test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.²
4. RPR Reagents must be used with the RPR Calibrator Set.
5. **Caution :** Avoid freezing reagents.
6. **Caution :** Reagent 1 and 2 contains <0.1% sodium azide as an antimicrobial agent. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.
7. Disposal of all waste material should be in accordance with local guidelines.

Preparation

Reagent 1: Liquid, ready to use

Reagent 2: Liquid, ready to use

Invert to mix daily before use. Avoid the formation of foam.

Storage and Stability

Unopened reagent is stable until the expiration date shown on the label when stored at 2 - 8°C.

Once opened, the reagent is stable up to 4 weeks at 2 - 8°C.

DO NOT FREEZE.

Onboard stability

Reagents are stable open on the Hitachi 917 analyzers for 28 days at 2 - 8°C.

Indications of Deterioration

Presence of turbidity in reagent 1 or microbial growth in either reagents may indicate deterioration.

Inability to recover control values.

SPECIMEN COLLECTION AND PREPARATION

Serum and sodium heparin plasma are the recommended collection media. Use standard sample collection and preparation methods.³

If not analyzed promptly, specimens may be stored as shown below:
Stability for serum⁸

7 days at 2 - 8°C
1 day at 15 - 25°C
4 weeks at (-15) – (-25)°C

Stability for plasma⁸

7 days at 2 - 8°C
1 day at 15 - 25°C
7 days at (-15) – (-25)°C


If the samples have been frozen, centrifuge at 15000xg for 10 minutes before measurement. Samples may be frozen and thawed once.⁴

PROCEDURE


Assay

Below are two examples of the RPR assay procedure for the Hitachi 7180 automated clinical analyzer and the SK500/Bolis 50i automated analyzer. All analyzer applications should be validated. Sekisui Diagnostics has applications for several automated analyzers available upon request. The performance of applications not validated by Sekisui must be defined and validated by the user.

(Hitachi 7180)

(SK500/ Bolis 50i)

Materials Provided

RPR Reagents 1 and 2 are required for the measurement of syphilitic anti-lipid antibodies. The RPR reagents are packaged and sold as a kit.

Description	Configuration	Catalog Number
RPR Reagent 1	1 x 60mL	486616
RPR Reagent 2	1 x 20mL	

Materials Required but not Provided

Description	Configuration	Catalog Number
RPR Calibrator Set	5 levels x 1mL	486623
RPR Control Set	Positive 2 x 1mL Negative 2 x 1mL	486630

• Analyzer capable of running two-reagent chemistries.

Calibration

Only the RPR Calibrators should be used to calibrate the RPR assay. The assigned values of the RPR calibrators are traceable to the WHO International standard (the International Standard for Syphilitic human serum (1st International standard preparation) established in 1958).

Calibration frequency should be determined by the user.

Quality Control values should be within the expected ranges.

Quality Control

Reliability of test results should be monitored routinely with quality control materials or serum pools that reasonably represent performance with patient specimens. Controls or serum pools should be used to monitor that the reagents are functioning properly and that correct procedures are being followed. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range, follow normal troubleshooting procedures. If assistance is required, please contact your local distributor.

Quality control requirements should be established in accordance with local, state and/or federal regulations, or accreditation requirements.

RESULTS

Results are expressed in RPR units (R.U.). Please note that R.U. is based on the WHO IU (International Standard for Syphilitic Human Serum). 1R.U. = 0.4IU. 1R.U. is equivalent to a 1-fold RPR card test titer.

To convert from R.U. to IU, divide the R.U. units by 2.5.

Limitations/Interfering Substances

Criterion: Recovery within \pm 20% of initial value

All studies were conducted on the Hitachi 7180 automated clinical analyzer.

Lipemia: No significant lipemic interference was observed up to 1% (Intralipid). If you suspect the sample is lipemic centrifuge at 15000xg for 10 minutes before measurement.

Chyle (formazin turbidity unit): No significant interference was observed up to 620 formazin degree. Samples exceeding 620 formazin degree should be centrifuged at 15000xg for 10 minutes before measurement.

Hemoglobin concentration of up to 488 mg/dL (75.7 μ mol/L) did not interfere in samples with syphilitic anti-lipid antibodies levels of 3.2 R.U.

Conjugated Bilirubin concentration of up to 21 mg/dL (359.1 μ mol/L) did not interfere in samples with syphilitic anti-lipid antibodies levels of 3.2 R.U.

Unconjugated Bilirubin concentration of up to 19.7 mg/dL (336.9 μ mol/L) did not interfere in samples with syphilitic anti-lipid antibodies levels of 3.2 R.U.

Rheumatoid factor was tested up to 450 IU/mL with no significant interference observed in samples with syphilitic anti-lipid antibodies levels of 3.8 R.U.

Serum and plasma samples from patients in the early stage of antibody production due to compromised immune function contain a small amount of antibody and may test negative.

A non-specific immune response may occur in serum and plasma samples from patients with autoimmune diseases. The test result should be evaluated based on other test results and clinical symptoms.

Serum and plasma samples from patients receiving blood products containing immunoglobulin may test positive. Evaluate the test result carefully

Due to the large range of specimen concentrations possible, sample probe washing must be adequate to prevent sample carryover.

Expected Values

A measurement of 1 R.U. or higher indicates that the sample is antibody positive.

A positive antibody test should be followed by subsequent tests and should be evaluated along with other test results and clinical symptoms. A final diagnosis of syphilis should be made by a physician. Results that do not match the clinical symptoms should be re-tested.

Each laboratory should confirm the reference interval for the patient population it serves.


SPECIFIC PERFORMANCE CHARACTERISTICS

Method Comparison

Comparative performance studies were conducted using the RPR Reagent on the Roche Hitachi 917 analyzer and a commercially available RPR card test. 83 serum samples, with RPR concentrations between 0.0 and 8.6 R.U. were tested.

The regression analysis is provided below:

Sekisui Diagnostics RPR vs. RPR Card Test (n = 83)	
Slope	1.01
Intercept (R.U.)	0.35
Correlation Coefficient (r)	0.909

Precision

Within-run precision of the RPR reagent was determined using 3 levels of positive serum samples and 1 negative serum sample according to internal protocol. Samples were measured in triplicate, 10 times using 3 lots of reagents on the Hitachi 7180 clinical analyzer.

Total precision of the RPR reagent was determined using 1 negative control and 1 positive control samples according to an internal protocol, 3 times measurements per day for 21 days on the Hitachi 7180 clinical analyzer.

Within-Run Precision

Serum pool / Control	Mean Recovery (R.U.)	Standard Deviation (R.U.)	%CV
Human Serum (Syphilis Negative)	0.0	0.0	—
Human Serum (Syphilis Positive)	1.7	0.1	4.3
Human Serum (Syphilis Positive)	2.9	0.1	1.6
Human Serum (Syphilis Positive)	6.9	0.1	1.6

Between- Run Precision

RPR Control	Mean Recovery (R.U.)	Standard Deviation (R.U.)	%CV
RPR Control (negative)	0.0	0.05	—
RPR Control (positive)	2.1	0.11	5.1

Limit of Detection (LoD)

The limit of detection is the actual concentration at which an observed test result is 2 SD above that of the lowest calibrator (calibrator 1). Limit of detection was established using calibrator 1 and 10 measurements with RPR Reagent on the Roche Hitachi 7180 clinical analyzer

The limit of detection is 0.2 R.U.

Specificity

For qualitative use of the assay, false negative results (<1 R.U.) (by a high dose hook effect) were not observed at analyte concentrations up to 100.5 R.U. syphilitic anti-lipid antibodies on the Hitachi 7180 clinical analyzer.

For the quantitative application, specimens above 4.0 R.U. should be re-run after dilution (1:10) because of a possible high-dose hook effect.

A negative result that does not match the clinical signs may occur with patients with hyperglobulinemia. In such cases the sample may be diluted with 0.9% NaCl solution and re-run to permit correct measurement.

Linearity

The RPR method is linear from 0.2 to 8.0 R.U. on the Hitachi 7180 automated clinical analyzer.

Specimens above 8.0 R.U. may be diluted with 0.9% NaCl solution. Multiply the result by the dilution factor to obtain the RPR concentration for the sample.

Other Performance Studies

Cross Reactivity (Analytical specificity)

Samples containing potentially interfering substances were analysed for cross-reactivity. Samples tested with the RPR assay were:

- from patients with collagenosis
- from patients undergoing dialysis
- from pregnant women

The following results were obtained.

Specimen	Number	RPR reactive	Rate
Collagenosis	28	1	4%
Pregnant	26	0	0%
Dialysis	50	3	6%

Clinical sensitivity^{5,6}

A total of 187 selected confirmed syphilis-positive samples in various stages of the disease were tested. The sensitivity with these samples was 99.5%.

		RPR	
Syphilis- positive samples	Number	Positive	Negative
187	186	1	

Clinical specificity⁷

A total of 2639 syphilis-negative samples were tested. The specificity with these samples was 99.5%.

		RPR	
Syphilis- negative samples	Number	Positive	Negative
	2639	12	2627

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1. Longo, DL; Fauci, AS; Kasper DL, Hauser SL, Jameson JL and Loscalzo J. Eds Harrison's Online-Harrison's Principles of Internal Medicine. 18e. Part 8. Infectious Diseases. McGraw-Hill
2. Richardson JH and Barkley WE, eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC: 1984.
3. Clinical and Laboratory Standards Institute, Procedures for the Handling and Processing of Blood Specimens: Approved Guideline. CLSI Document H18-A, Wayne, PA: 1990.
4. Shibasaki M et al. An Automated Measurement of Anti-Treponema Antibody Titer by MEDIACE TPLA, a Latex Agglutination Test using Hitachi 7170 Automatic Analyzer. The Journal of Clinical Laboratory Instruments and Reagents 1996; 19 (4):635–639.
5. Osato K et al. Clinical Evaluation of Latex Agglutination Test Kits for Detecting Anti-syphilitic Lipoidal Antibodies and Anti-treponemal Antibodies. Japanese Journal of Sexually Transmitted Diseases 2002; 13 (1): 124–130.
6. Kawai K et al. The possibility of assessing the stage of infection by using Mediace TPLA and RPR. The Journal of Clinical Laboratory Instruments and Reagents 2003; 26 (4): 301–304.
7. Kinjo T et al. Laboratory -based evaluation of Latex Agglutination Turbidimetric Assay by Mediace RPR on P Module of Hitachi Auto analyzer 7600 to Quantitatively Determine Serum RPR Antibody. Japanese Journal of Clinical Laboratory Automation (JJCLA) 2005; 30 (3): 257–262.
8. Data on file at SEKISUI.

Definitions for Symbols/ Definition Der Symbole/ Définition des symboles/ Definizioni dei simboli/ Definiciones de los símbolos/ Definice symbolů/ Определение символов

REF

Catalog number
Bestellnummer
Référence
Número de catálogo
Número de catálogo
Katalogové číslo
Номер по каталогу

Temperature limitation
Temperaturbegrenzung
Limites de température
Limite di temperatura
Limite temperature
Teplotní omezení
Температурное ограничение

Use by
Zu verwenden bis
Utiliser avant
Usare entro il
Fecha de caducidad
Datum spotřeby
Использовать до

LOT

Batch code
Chargennummer
Code de lot
Codice del lotto
Código del lote
Kód šarže
Код партии

EC REP

Authorized Representative in the European Community
Autorisierte Vertretung in der Europäischen Gemeinschaft
Représentant agréé dans la Communauté européenne
Rappresentante autorizzato nell'Unione Europea
Representante autorizado en la Comunidad Europea
Autorizovaný zástupce pro Evropské společenství
Уполномоченный представитель в Европейском сообществе

SEKISUI MEDICAL CO., LTD.
1-3, Nihonbashi 2-chome, Chuo-ku
Tokyo, 103-0027, Japan

EC REP

Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41
30175 Hannover, Germany

SEKURE™
CHEMISTRY

Sep2022 KI486616.05

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.

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Das Wort SEKURE und das SEKURE-Logo sind Markenzeichen von Sekisui Diagnostics, LLC.
Le terme SEKURE et le logo Sekure sont des marques déposées de Sekisui Diagnostics, LLC.
La parola SEKURE e il logo Sekure sono marchi commerciali registrati di Sekisui Diagnostics, LLC.
La palabra SEKURE y el logotipo de SEKURE son marcas comerciales registradas de Sekisui Diagnostics, LLC.
Slovo SEKURE a logo Sekure jsou registrované ochranné známky společnosti Sekisui Diagnostics, LLC.