# EN TPLA REAGENT

## INTENDED USE

For the quantitative determination of anti-*Treponema pallidum* antibodies in human serum or plasma on the Hitachi 9000 automated clinical analyzer, the SK500/Biolis 50i automated clinical analyzer and multiple other analyzers.

## SUMMARY

Syphilis is a chronic infection caused by *Treponema pallidum* which can be transmitted congenitally or by sexual contact. It is characterized by episodes of active disease interrupted by periods of latency<sup>1</sup>.

There are two types of serological test for syphilis, non-treponemal and treponemal. The most widely used non-treponemal antibody tests for syphilis are the rapid plasma reagin (RPR) and venereal disease research laboratory (VDRL) tests, both measuring antibodies against a cardiolipin-lecithin-cholesterol antigen complex. Treponemal tests measure antibodies to native (Nichols strain) or recombinant *T. pallidum* antigens.<sup>1</sup> This automated assay is based on a latex immunological agglutination test principle.

## PRINCIPLE

The polystyrene latex, coated with antigen components derived from *Treponema pallidum* (Nichols strain), is exposed to the test sample under certain conditions, to induce formation of anti-*Treponema* antibody-latex aggregate. The elevation of turbidity due to the formation of this aggregate (the change in turbidity) relative to the pre-exposure level is measured, to determine the anti-*Treponema* antibody titer in the test sample.

## REAGENTS

### Composition

Component	Ingredients	Concentration
Reagent 1	Bovine serum albumin Sodium Azide Buffer (pH 7.2-7.4)	10% <0.1%
Reagent 2	Latex particles coated with <i>Treponema pallidum</i> derived-antigen Bovine serum albumin Sodium Azide Buffer (pH 7.2-7.5)	≤3.5mg/mL 1% <0.1%

### Precautions and Warnings

- For In Vitro Diagnostic Use.
- Do not use the reagents beyond the expiration date printed on the label.
- Warning:** All specimens used in the test should be considered potentially infectious. Universal precautions as they apply to your facility should be used for handling and disposal of materials during and after testing.<sup>2</sup>
- TPLA Reagents must be used with the TPLA Calibrator Set.
- Caution:** Avoid freezing reagents.
- Caution:** Reagents 1 and 2 contain <0.1% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive metal azide buildup. Flush with copious amounts of water when discarding material.
- Caution:** Do not pool reagents within a kit or between reagent kits.
- Disposal of all waste material should be in accordance with local guidelines.

### Preparation

Reagent 1: Liquid, ready to use.  
Reagent 2: Liquid, ready to use.  
Invert to mix before use. Avoid the formation of foam.

### Storage and Stability

Unopened reagent is stable until the expiration date shown on the label when stored at 2 - 8°C.  
Once opened, the reagent is stable up to 4 weeks at 2 - 8°C.  
DO NOT FREEZE.

### Onboard stability

Reagents are stable open on the Hitachi 9000 analyzer for 4 weeks at 2 - 8°C.

### Indications of Deterioration

Presence of turbidity in reagent 1 or microbial growth in either reagents may indicate deterioration.  
Inability to recover control values.

## SPECIMEN COLLECTION AND PREPARATION

Serum and lithium heparin or EDTA-2K plasma are the recommended collection media. Use standard sample collection and preparation methods.<sup>3</sup>

If not analyzed promptly, serum or plasma specimens may be stored at 2 - 8°C for 1 week or at 15 - 25°C for 1 day.<sup>4</sup> If specimens need to be stored for more than 1 week, they may be preserved at -20°C or below for up to 4 weeks.<sup>4</sup>


If the samples have been frozen, centrifuge at 15000xg for 10 minutes before measurement. Samples may be frozen and thawed once.<sup>4</sup>

## PROCEDURE


### Assay

Below are two examples of the TPLA assay procedure for the Hitachi 9000 automated clinical analyzer and SK500/Biolis 50i automated analyzer. All analyzer applications should be validated. Sekisui Diagnostics has applications for several automated analyzers available upon request. The performance of applications not validated by Sekisui must be defined and validated by the user.

(Hitachi 9000)

(SK500/ Biolis 50i)

### Materials Provided

TPLA Reagents 1 and 2 are required for the measurement of anti-*Treponema pallidum* antibodies. The TPLA reagents are packaged and sold as a kit.

Description	Configuration	Catalog Number
TPLA Reagent 1	1 x 60mL	486647
TPLA Reagent 2	1 x 10mL	

### Materials Required but not Provided

Description	Configuration	Catalog Number
TPLA Calibrator Set	5 levels x 2mL	515132
TPLA Control Set	Level A 1 x 3mL	515149
	Level B 1 x 3mL	

- Analyzer capable of running two-reagent chemistries.

### Calibration

Only the TPLA Calibrator set should be used to calibrate the TPLA assay. The assigned values of the TPLA Calibrators are traceable to an in-house standard.

Calibration frequency should be determined by the user.

Quality Control values should be within the expected ranges.

### Quality Control

Reliability of test results should be monitored routinely with quality control materials or serum pools that reasonably represent performance with patient specimens. Controls or serum pools should be used to monitor that the reagents are functioning properly and that correct procedures are being followed. An acceptable range for each lot of control material should be established by the laboratory. If control values are not within the expected range, follow normal troubleshooting procedures. If assistance is required, please contact your local distributor.

Quality control requirements should be established in accordance with local, state and/or federal regulations, or accreditation requirements.

## RESULTS

Results are expressed in titer units (T.U.). T.U. are based on the *Treponema pallidum* haemagglutination (TPHA) assay. 1280 T.U. is

equivalent to a titer of 1:1280 TPHA.

To convert from T.U. (Titer Units) to mIU multiply T.U. units by 2.

### Limitations/Interfering Substances

Criterion: Recovery within  $\pm 20\%$  of initial value. All studies were conducted on the Hitachi 9000 automated clinical analyzer.

**Lipemia:** No significant lipemic interference was observed up to 0.25% (intralipid) in samples with syphilitic anti-lipid antibodies of approximately 80 T.U. If you suspect the sample is lipemic centrifuged at 15000xg for 10 minutes before measurement.

**Chyle (formazin turbidity unit):** No significant interference was observed up to 1450 formazin degree in samples with syphilitic anti-lipid antibodies of approximately 80 T.U. Samples exceeding 1450 formazin degree should be centrifuged at 15000xg for 10 minutes before measurement.

Hemoglobin concentration of up to 490 mg/dL (76.0  $\mu\text{mol/L}$ ) did not interfere in samples with syphilitic anti-lipid antibodies levels of 80 T.U.

Conjugated Bilirubin concentration of up to 18.5 mg/dL (316.4  $\mu\text{mol/L}$ ) did not interfere in samples with syphilitic anti-lipid antibodies levels of 80 T.U.

Unconjugated Bilirubin concentration of up to 18.5 mg/dL (316.4  $\mu\text{mol/L}$ ) did not interfere in samples with syphilitic anti-lipid antibodies levels of 80 T.U. Rheumatoid factor was tested to 500 IU/mL and did not affect performance on the Roche Hitachi H9000 clinical analyzer.

Serum and plasma samples from patients in the early stage of antibody production due to compromised immune function contain a small amount of antibody and may test negative.

A non-specific immune response may occur in serum samples from patients with autoimmune diseases. The test result should be evaluated based on other test results and clinical symptoms.

Serum samples from patients receiving blood products containing immunoglobulin may test positive. Evaluate the test result carefully.

Due to the large range of specimen concentrations possible, sample probe washing must be adequate to prevent sample carryover.

### Expected Values

A measurement of 10 T.U. or higher indicates that the sample is antibody positive.

A positive antibody test should be followed by subsequent tests and should be evaluated along with other test results and clinical symptoms. A final diagnosis of syphilis should be made by a physician. Results that do not match the clinical symptoms should be re-tested.

Each laboratory should confirm the reference interval for the patient population it serves.


## SPECIFIC PERFORMANCE CHARACTERISTICS

### Method Comparison

Comparative performance studies were conducted using the TPLA Reagent on the Roche Hitachi 7050 analyzer and a commercially available *Treponema pallidum* hemagglutination (TPHA) method. 171 serum samples, with positive TPLA sample concentrations between 27 and 585 TU.

The regression analysis is provided below:

Sekisui Diagnostics TPLA vs. an existing TPHA test method (n = 171)	
Slope	1.06
Intercept (T.U.)	10.3
Correlation Coefficient (r)	0.828

### Precision

Within-run precision of the TPLA reagent was determined using 3 levels of positive serum samples and 1 negative serum sample according to an internal protocol. Samples were measured in triplicate, 10 times using 3 lots of reagents on the Hitachi H9000 automated clinical analyzer.

Between-Run precision of the TPLA reagent was determined using 1 negative serum sample and 2 positive control samples according to an internal protocol, 2 measurements per day for 21 days on the Hitachi H9000 automated clinical analyzer.

### Within-Run Precision

Samples	Mean (T.U.)	SD (T.U.)	CV (%)
Human serum 1	0.0	0.0	—
Human serum 2	28.8	1.1	3.7
Human serum 3	77.4	2.0	2.5
Human serum 4	190.9	8.4	4.4

### Between-Run Precision

Sample/Control	Mean (T.U.)	SD (T.U.)	CV (%)
Human serum	0.0	0.0	—
TPLA Control Level A	26.5	1.4	5.2
TPLA Control Level B	80.3	3.6	4.5

A positive control measured 42 times gave a CV of  $\leq 15\%$ .

### Limit of Detection (LoD)

The limit of detection is the actual concentration at which an observed test result is 2.0 SD above that of the lowest calibrator (calibrator 1). Limit of detection was established using the level 1 calibrator and 10 measurements with TPLA Reagent on the Roche Hitachi H9000 clinical analyzer.

The limit of detection has been defined as: 4.6 T.U.

### Specificity

A high dose hook effect was not observed at analyte concentrations up to 362 T.U. on the Hitachi H9000 automated clinical analyzer

A negative result that does not match the clinical signs may occur with patients with hyperglobulinemia. In such cases the sample may be diluted with 0.9% NaCl solution and rerun.

### Linearity

The TPLA method is linear from 5 T.U. to 250 T.U. on the Roche Hitachi H9000 clinical analyzer.

Specimens above 250 T.U. may be diluted with 0.9% NaCl solution. Multiply the result by the dilution factor to obtain the TPLA concentration for the sample.

### Other Performance Studies

#### Cross Reactivity (Analytical specificity)

Samples containing potentially interfering substances were analysed for cross-reactivity. Samples tested with the Sekisui Diagnostics TPLA assay were:

- from patients with collagenosis and patients undergoing dialysis
- from pregnant women

Specimen	Number	TPLA Reactive	Rate
Collagenosis patients	28	0	0%
Pregnant women	26	0	0%
Dialysis patients	50	0	0%

No false-positive results were found.

#### Clinical sensitivity<sup>5,6,7</sup>

A total of 268 confirmed syphilis-positive samples in various stages of the disease were tested with the Sekisui Diagnostics TPLA assay. The sensitivity was found to be 100%

		TPLA	
Syphilis- positive samples	Number	Positive	Negative
	268	268	0

#### Clinical specificity<sup>4,5,6</sup>

A total of 3427 syphilis-negative samples were tested with the Sekisui Diagnostics TPLA assay. The specificity with these samples was 99.6%.

		TPLA	
Syphilis-negative samples	Number	Positive	Negative
	3427	13	3414

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
1. Longo, DL, Fauci, AS, Kasper DL, Hauser SL, Jameson JL and Loscalzo J. Eds Harrison's Online-Harrison's Principles of Internal Medicine. 18e. Part 8. Infectious Diseases. McGraw-Hill
2. Richardson JH and Barkley WE, eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, HHS Publication No. (CDC) 84-8395, Washington, DC: 1984.
3. Clinical and Laboratory Standards Institute, Procedures for the Handling and Processing of Blood Specimens: Approved Guideline. CLSI Document H18-A, Wayne, PA: 1990.
4. Shibazaki M et al. An Automated Measurement of Anti-Treponema Antibody Titer by MEDIACE TPLA, a Latex Agglutination Test using Hitachi 7170 Automatic Analyzer. The Journal of Clinical Laboratory Instruments and Reagents 1996; 19 (4):635-639.
5. Osato K et al. Clinical Evaluation of Latex Agglutination Test Kits for Detecting Anti-syphilitic Lipoidal Antibodies and Anti-treponemal Antibodies. Japanese Journal of Sexually Transmitted Diseases 2002; 13 (1):124-130.
6. Osato K et al. Clinical Evaluation of Automated Latex Agglutination Test Kits (TPLA) for Syphilis Diagnosis. The Journal of Clinical Laboratory Instruments and Reagents 1991; 14 (4):739-743.
7. Kataniwa Y et al. Clinical Evaluation of Latex Reagent Samedia TPLA for Diagnosis of Syphilis. The Journal of Clinical Laboratory Instruments and Reagents 1991; 14 (4):735-738.
8. Data on file at SEKISUI.

## Definitions for Symbols/ Definition Der Symbole/ Définition des symboles/ Definizioni dei simboli/ Definiciones de los símbolos/ Definicje symboló/ Определение символов

REF	IVD
Catalog number Bestellnummer Référence Numero di catalogo Número de catálogo Katalogové číslo Номер по каталогу	For in vitro diagnostic use Für die In-Vitro-Diagnostik bestimmt Pour diagnostic in vitro Per uso diagnostico in vitro Para uso en diagnóstico in vitro K diagnostickému použití in vitro Только для in vitro диагностики

Temperature limitation  
Temperaturbegrenzung  
Limites de température  
Limite di temperatura  
Limite temperatura  
Teplotní omezení  
Температурное ограничение

Manufactured by  
Hergestellt von  
Fabriqué par  
Prodotto da  
Fabricante  
Výrobce  
Изготовлено

Use by  
Zu verwenden bis  
Utiliser avant  
Usare entro il  
Fecha de caducidad  
Datum spotřeby  
Использовать до

Consult instructions for use  
Gebrauchsanweisung beachten  
Consulter le mode d'emploi  
Consultare le istruzioni per l'uso  
Consulte las instrucciones de uso  
Prostudujte si návod k použití  
Обратитесь к инструкции по применению

Batch code  
Chargenkennung  
Code de lot  
Codice del lotto  
Código del lote  
Kód šarže  
Код партии

Caution, consult accompanying documents  
Achtung, Begleitunterlagen beachten  
Attention, consulter la documentation jointe  
Attenzione, consultare i documenti di accompagnamento  
Precaución, consulte los documentos anexos  
Upozornění, prostudujte si příloženou dokumentaci  
Внимание! Обратитесь к сопровождающей документации

Authorized Representative in the European Community  
Autorisierte Vertretung in der Europäischen Gemeinschaft  
Représentant agréé dans la Communauté européenne  
Rappresentante autorizzato nell'Unione Europea  
Representante autorizado en la Comunidad Europea  
Autorizovaný zástupce pro Evropské společenství  
Уполномоченный представитель в Европейском сообществе

Biological risks  
Biologisches Risiko  
Risques biologiques  
Rischi biologici  
Riesgos biológicos  
Biologická rizika  
Биологический риск

  
SEKISUI MEDICAL CO., LTD.  
1-3 , Nihonbashi 2 -chome, Chuo-ku  
Tokyo, 103-0027, Japan

  
Medical Device Safety Service GmbH (MDSS)  
Schiffgraben 41  
30175 Hannover, Germany

  
April 2019 KI486647. 05

The word SEKURE and the Sekure logo are registered trademarks of Sekisui Diagnostics, LLC.  
Das Wort SEKURE und das SEKURE-Logo sind Marken-zeichen von Sekisui Diagnostics, LLC.  
Le terme SEKURE et le logo Sekure sont des marques déposées de Sekisui Diagnostics, LLC.  
La parola SEKURE e il logo Sekure sono marchi commerciali registrati di Sekisui Diagnostics, LLC.  
La palabra SEKURE y el logotipo de SEKURE son marcas comerciales registradas de Sekisui Diagnostics, LLC.  
Slovo SEKURE a logo Sekure jsou registrované ochranné známky společnosti Sekisui Diagnostics, LLC.  
Слово SEKURE и логотип Sekure являются зарегистрированными товарными знаками компании Sekisui Diagnostics, LLC.