## UNSATURATED IRON BINDING CAPACITY (UIBC) ASSAY

**CATALOGUE NUMBER:** 153-10    **SIZE:** R1: 1 x 100 mL, R2: 1 x 25 mL  
 153-30                  R1: 3 x 100 mL, R2: 1 x 75 mL  
 153-50                  R1: 1 x 500 mL, R2: 1 x 300 mL  
 153-90                  R1: 1 x 1000 mL, R2: 1 x 240 mL

**Note:** Changes are highlighted.

### INTENDED USE


For the *in vitro* quantitative measurement of unsaturated iron binding capacity (UIBC) in serum.

### TEST SUMMARY

The measurement of unsaturated iron binding capacity (UIBC) in combination with serum iron is a useful diagnostic tool in the determination of various iron disorders. The combined value of UIBC and serum iron gives a value for the total iron binding capacity (TIBC). This represents the maximum concentration of iron that serum proteins can bind. Serum UIBC levels vary in disorders of iron metabolism, where iron binding capacities are often increased in iron deficiency and decreased in chronic inflammatory disorders or malignancies.<sup>(1)</sup>

In 1970, Stokey<sup>(2)</sup> reported the synthesis of 3-(2-pyridyl)-5,6-bis (4-phenylsulfonic acid)-1,2,4-triazine, monosodium salt (Ferrozine®) which complexed with ferrous iron to form a tris ferrozine/iron, Fe(Fz)<sub>3</sub> complex. The major advantages of ferrozine are the high molar absorptivity of the ferrozine/iron complex (28,000), its water solubility, and stability over the pH range of 4-9. This assay uses a ferroin type compound called 5,5'(3-[2-pyridyl]-1,2,4-triazine-5,6 diyl)-bis-2-furansulfonic acid, disodium salt (Ferene®).<sup>(3,4,5)</sup> This reagent is a superior iron chelating agent forming a Ferene® complex with ferrous iron with a maximum absorbance at 593 nm and a molar absorptivity of 35,500. The compound has a 27% higher molar absorption than ferrozine, absorbs at a longer wavelength, and has the other advantages of ferrozine, namely its solubility and stability.

### TEST PRINCIPLE

A known ferrous ion concentration incubated with serum binds specifically with transferrin at unsaturated iron binding sites.

Remaining unbound ferrous ions are measured with the ferene reaction. The difference between the amount of unbound iron and the total amount added to the serum is equivalent to the quantity bound to transferrin. This is the UIBC of the sample.

### REAGENTS

Binding Buffer Reagent (R1): A buffer (pH 8.7 at 25°C) containing 12.8 μmol/L Ferrous Ammonium Sulfate, surfactants, stabilizers, and a preservative.

Color Reagent (R2): A solution containing 240 mmol/L Ascorbic Acid, 5.4 mmol/L Ferene, and a stabilizer.

### WARNINGS & PRECAUTIONS FOR USE

For *in vitro* diagnostic use

### UIBC-SL Binding Buffer Reagent (R1)

#### Warning

Contains: thiourea; thiocarbamide (CAS No) 62-56-6, sodium hydroxide; caustic soda (CAS No) 1310-73-2

#### Hazard statements

H315 - Causes skin irritation

H319 - Causes serious eye irritation

H351 - Suspected of causing cancer

H361 - Suspected of damaging fertility or the unborn child

H412 - Harmful to aquatic life with long lasting effects

#### Precautionary statements

P201 - Obtain special instructions before use.

P202 - Do not handle until all safety precautions have been read and understood.

P264 - Wash hands, forearms and face thoroughly after handling.

P273 - Avoid release to the environment.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 - If on skin: Wash with plenty of water.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 - If exposed or concerned: Get medical advice/attention.

P332+P313 - If skin irritation occurs: Get medical advice/attention.

P337+P313 - If eye irritation persists: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

P405 - Store locked up.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

### UIBC-SL Color Reagent (R2)

#### Warning

Contains: thiourea; thiocarbamide (CAS No) 62-56-6

#### Hazard statements

H351 - Suspected of causing cancer

H361 - Suspected of damaging fertility or the unborn child

#### Precautionary statements

P201 - Obtain special instructions before use.

P202 - Do not handle until all safety precautions have been read and understood.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P308+P313 - If exposed or concerned: Get medical advice/attention.

P405 - Store locked up.

P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

Avoid contact with skin and eyes.

Harmful if swallowed.

May cause sensitization by skin contact.

See Safety Data Sheet for additional information.

### REAGENT PREPARATION, STORAGE AND STABILITY

Reagents are ready for use.

Supplied reagent is stable until expiry date when stored at 2-8°C and protected from light. Stability claims are based on real time studies.

### REAGENT DETERIORATION

The reagent solution should be clear. Turbidity would indicate deterioration.

### DISPOSAL

Reagents must be disposed of in accordance with all Federal, Provincial, State, and local regulations. Reagent 1 contains Triton detergent. Dispose of waste and unused Reagent 1 by incineration.

### SPECIMEN

Fresh, clear, unhemolysed serum. The sample should be drawn in the morning following a 12 hour fast. Blood collecting glassware should be free of iron contamination.<sup>(1)</sup>

### SAMPLE STORAGE

If a sample is to be stored for longer than 8 hours, refrigeration at 2-8°C is recommended.<sup>(6)</sup>

### ANALYTICAL SPECIFICITY (CLSI EP7)<sup>(7)</sup>

Cross contamination studies have not been performed on automated instruments. Certain reagent/instrument combinations used in sequence with this assay may interfere with reagent performance and test results. The existence of, or effects of, any potential cross contamination issues are unknown.

Interference from icterus, lipemia, and hemoglobin was evaluated for this method on a Roche/Hitachi® 911 analyzer using a significance criterion of >10% variance from control.

Concentration of Analyte	Substance Tested	Concentration of Interferent Where Interference is Insignificant
μg/dL	μmol/L	
323	57.6	Hemoglobin      400 mg/dL      62 μmol/L
236	42.2	Bilirubin      40 mg/dL      684 μmol/L
310	55.5	Intralipid      1000 mg/dL      3000 mg/dL (33.9 mmol/L) Simulated Triglycerides

Copper is the only cation of the trace metals usually present in serum to form a colored complex with ferene. Copper interference with ferene is similar to that encountered with ferrozine and studied by Duffy and Gaudin.<sup>(8)</sup> Ninety five percent of the copper interference is eliminated by chelation of free copper.

Samples containing the following should not be used: deferoxamine.

A summary of the influence of drugs on clinical laboratory tests may be found by consulting Young, D.S.<sup>(9)</sup>

The information presented above is based on results from SEKISUI Diagnostics' studies and is current at the date of publication.

### ANALYTICAL PROCEDURE

#### MATERIALS PROVIDED

SEKISUI Diagnostics' UIBC reagent.

#### MATERIALS REQUIRED (BUT NOT PROVIDED)

1. Automated analyzer capable of accurately measuring absorbance at the appropriate wavelength as per the instrument application.
2. Calibration material.
3. Quality Control materials.

#### TEST CONDITION

For the data presented in this insert, studies using this reagent were performed on an automated analyzer using an endpoint test mode, with a sample to reagent ratio of 1:18.5, and wavelength readings of (primary/secondary) 600/700 nm. For assistance with applications on automated analyzers within Canada and the U.S., please contact SEKISUI Diagnostics Technical Services at 1-800-565-0265. Outside Canada and the U.S., please contact your local distributor.

### CALIBRATION

Calibration material should be used to calibrate the procedure. The frequency of calibration is dependent on the system and the parameters used.

## QUALITY CONTROL

A normal and abnormal concentration control should be analyzed as required in accordance with local, state and federal guidelines. The results should fall within the acceptable range as established by the laboratory.

## CALCULATIONS

The analyzer automatically calculates the UIBC concentration of each sample.

## TEST LIMITATIONS

A sample with a UIBC concentration exceeding the linearity limit should be diluted with 0.9% saline and reassayed incorporating the dilution factor in the calculation of the value.

## REFERENCE INTERVALS

UIBC reference intervals were determined from an apparently healthy population of 203 people following CLSI C28<sup>(7)</sup> guidelines. Serum samples were analyzed on automated instrumentation at 37°C.

126-382 µg/dL (22.6-68.4 µmol/L)

These values are suggested guidelines. It is recommended that each laboratory establish its own expected range.

## PERFORMANCE CHARACTERISTICS

Data presented was collected on a Roche/Hitachi® 704 analyzer unless otherwise stated.

## RESULTS

UIBC concentration is reported as µg/dL (µmol/L).

## REPORTABLE RANGE (CLSI EP6)<sup>(7)</sup>

The linearity of the procedure described is 600 µg/dL (107.4 µmol/L). The lower limit of detection of the procedure described is 17 µg/dL (3.0 µmol/L). This data results in a reportable range of 17 to 600 µg/dL (3.0 to 107.4 µmol/L).

## PRECISION STUDIES (CLSI EP5)<sup>(7)</sup>

Total precision data was collected on two concentrations of control sera in 40 runs conducted over 20 days.

Concentration	Total SD		Total CV (%)	Concentration		Within SD		Within Run CV (%)
	µg/dL	µmol/L		µg/dL	µmol/L	µg/dL	µmol/L	
126	22.6	7.2	1.29	5.7	129	23.0	6.4	1.15
209	37.4	8.6	1.54	4.1	203	36.3	5.5	0.98

Within run precision data was collected on two concentrations of control sera each run 20 times in a single assay.

## ACCURACY (CLSI EP9)<sup>(7)</sup>

The performance of this method (y) was compared with the performance of a similar UIBC method (x) on a Roche/Hitachi® 704 analyzer. Forty patient serum samples ranging from 94-437 µg/dL (16.8-78.2 µmol/L) were tested and gave a correlation coefficient of 0.9937. Linear regression analysis gave the following equation:

This method = 1.065 (reference method) + 2.0 µg/dL (0.4 µmol/L).

## TRADEMARKS

The word SEKURE and the Sekure logo are trademarks of SEKISUI Diagnostics, LLC.

All other trademarks are the property of their respective owners.

## **FR**

# **UNSATURATED IRON BINDING CAPACITY (UIBC) ASSAY**

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1. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., p. 1578 (1986).
2. Stokey, L.L., Ferrozine-A New Spectrophotometric Reagent for Iron, Anal. Chem. 42, 779 (1970).
3. Artiss, J.D., Vinogradov, S., Zak, B., Spectrophotometric Study of Several Sensitive Reagents for Serum Iron, Clin. Biochem. 14, 31-315 (1981).
4. Higgins, T., Novel Chromogen for Iron Determinations, Clin. Chem. 27, 1619 (1981).
5. Artiss, J.D., Strandbergh, D.R., Zak, B., Study of Continuous Flow Automation for Serum Iron on Comparing Several Sensitive Reagents, Microchemical Journal, 28, 275-284 (1983).
6. World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations. Geneva: World Health Organization; 2002:39.
7. CLSI Method Evaluation Protocols, Clinical and Laboratory Standards Institute, Wayne, PA.
8. Duffy, J.R., Gaudin, J., Copper Interference in the Determination of Iron in Serum using Ferrozine, Clin. Biochem. 10, 122-123 (1977).
9. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, AACC Press, Third Edition, Washington, 1990.


MDSS GmbH  
Schiffgraben 41,  
30175 Hannover, Germany  
**Phone:** (+49)-511-6262 8630  
**Fax:** (+49)-511-6262 8633

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