

(Diagnostic Reagent Grade)


ASAHI KASEI ENZYMES T-101

CHOLESTEROL OXIDASE [CON II -FD]

Lyophilized type

from *Rhodococcus* sp.

(Cholesterol: oxygen oxidoreductase, EC 1.1.3.6)



Preparation and Specification

Appearance : Yellowish amorphous powder, lyophilized
 Specific activity : More than 15 U/mg solid

Properties

| | | |
|-----------------------|--|-----------------------|
| Substrate specificity | : See Table 1 | |
| Molecular weight | : 61.8 KDa (SDS-PAGE) | |
| Isoelectric point | : pH 4.5 | |
| Michaelis constant | : Cholesterol $6.0 \times 10^{-5}\text{M}$ | |
| Optimum pH | : 7.0-7.5 | Figure 1 |
| pH stability | : 5.7-7.8 (65°C, 10 min) | Figure 2 |
| Optimum temperature | : 50°C (Tris-HCl buffer) | Figure 3 |
| Thermal stability | : Stable at 65°C and below (pH 7.0, 10 min) | Figure 4 and Figure 5 |
| Effects of detergents | : See Table 2 | |

Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of **total cholesterol**, **HDL-C**, and **LDL-C** when coupled with cholesterol esterase (T-18 and T-98).





Table 1. Substrate specificity

| Substrate (1mM) | Relative activity (%) |
|--------------------------|-----------------------|
| Cholesterol | 100 |
| β -Cholesterol | 93 |
| Pregnenolone | 98 |
| Dehydro-iso-androsterone | 10 |
| β -Sitosterol | 94 |
| Stigmasterol | 66 |
| Androsterol | 2 |
| Teststerone | 1 |
| Cholic acid | 3 |

Table 2. Effect of detergents on CON II -FD activity


| Detergents (0.1%) | Relative activity (%) |
|-------------------|-----------------------|
| Triton X-100 | 100 |
| Emulgen 810 | 101 |
| Emulgen 911 | 113 |
| Emulgen 709 | 107 |
| Emulgen 109P | 118 |
| Adekamol B-797 | 103 |
| Adekamol SO-120 | 100 |
| Adekamol 720N | 114 |
| RHEODOL 460 | 63 |
| SM 1080 | 122 |

Fig.1 pH Optimum




○ : MES buffer
 ▲ : PIPES buffer
 □ : Tris-HCl buffer
 ● : Glycine-NaOH buffer

Fig.2 pH Stability




65°C, 10 min.
 ○ : MES buffer
 ▲ : PIPES buffer
 □ : Tris-HCl buffer
 ● : Glycine-NaOH buffer

Fig.3 Optimum Temperature




pH 7.0
 100 mM Tris-HCl buffer

Fig.4 Thermal Stability



pH 7.0, 10 min.
 100 mM Tris-HCl buffer

Fig.5 Thermal Stability




100mM, 10min.
 ○ : MES buffer(pH5.95)
 □ : Tris-HCl buffer(pH7.63)
 ● : Glycine-NaOH buffer(pH8.95)

Assay

Principle

The assay is based on the increase in absorbance at 240 nm as Δ^4 -cholesten-3-one is produced in the following reactions:



Unit definition

One unit is defined as the amount of enzyme which liberates 1 μ mole of Δ^4 -cholesten-3-one per minute at 37°C under the conditions specified in the assay procedure.

Reagents

- Substrate solution (6 mM cholesterol solution)
 Dissolve 232 mg of cholesterol with isopropanol to make a total of 100 ml.

2. Enzyme dilution buffer

0.1 M KH_2PO_4 - Na_2HPO_4 buffer pH 7.0 containing
0.05% (W/V) Triton X-100

※ Prepare the enzyme dilution buffer two days before
use and keep it in the refrigerator until use.

3. Reagents

Cholesterol : NACALAI TESQUE, INC. Special grade
#08721

Triton X-100 : The Dow Chemical Company

■ Enzyme solution

Accurately weigh about 20 mg of the sample and add enzyme dilution buffer to make a total of 20 ml. After 1-1.5 hour incubation at room temperature, dilute it with enzyme dilution buffer to adjust the concentration to within 0.1-0.2 U/ml.

■ Procedure

- Pipette accurately 3.0ml of enzyme dilution buffer and 50 μl of enzyme solution and preincubate at 37°C.
※ In the case of a test blank, add 50 μl of enzyme dilution buffer in place of enzyme solution.
- After 5 min, add 50 μl of substrate solution and mix to start the reaction at 37°C.
- After starting the reaction, measure the rate of increase per minute in absorbance at 240nm. The rate must be measured within the linear portion of the absorbance curve.

Absorbance sample : As/min
blank : Ab/min

$0.010 \text{ Abs/min} \leq \Delta A/\text{min} = (\text{As/min} - \text{Ab/min})$
 $\leq 0.060 \text{ Abs/min}$

CON II-FD 活性測定法 (Japanese)

I. 試薬液

- 基質溶液 (6mM コレステロール溶液)
コレステロール 232mg をイソプロパノールに溶解して全容 100ml とする。
- 酵素溶解希釈用液
0.05% (W/V) トリトン X-100 を含む
0.1M KH_2PO_4 - Na_2HPO_4 緩衝液 pH7.0
※ 酵素溶解希釈用液は使用する2日前に調製し、
使用まで冷蔵保存する。
- 試薬
コレステロール : ナカライテスク社製 特級 #08721
トリトン X-100 : Dow Chemical 社製

II. 酵素試料液

検品約 20mg を精密に量り、酵素溶解希釈用液で溶解して全容 20ml とする。室温にて 1~1.5 時間放置し、その液を酵素溶解希釈用液で約 0.1-0.2U/ml 濃度となるように適宜希釈する。

III. 測定操作法

- 小試験管に酵素溶解希釈用液 3.0ml と酵素試料液 50 μl を正確に加え 37°C で予備加温する。
※ 盲検は酵素試料液の代りに酵素溶解希釈用液 50 μl を加える。

■ Calculation

$$\text{Activity (U/mg)} = \frac{\Delta A/\text{min}}{12.2} \times \frac{3.10}{0.05} \times \frac{1}{X}$$

12.2 : millimolar extinction coefficient of Δ^4 -Cholesten-3-one at 240 nm ($\text{cm}^2 / \mu\text{mole}$)

3.10 : final volume (ml)

0.05 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution (mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended. Enzyme activity will be retained for at least one year under this condition.

References

- Richmond, W. (1973) Clin. Chem., **19**, 1350.
- Flegg, H. M. (1973) Ann. Clin. Biochem., **10**, 79.
- Alain, C. C. et. al. (1973) Clin. Chem., **20**, 470.
- Tarbutton, P. N. and Gunter, C. R. (1974) Clin. Chem., **20**, 724.
- Nomoto, S. (1976) Rinsho Kensa, **20**, 688.
- Kameno, K., Nakano, N. and Baba, S. (1976) Jap. J. Clin. Path., **24**, 650.

IV. 計算

$$\text{活性 (U/mg)} = \frac{\Delta A/\text{min}}{12.2} \times \frac{3.10}{0.05} \times \frac{1}{X}$$

12.2: Δ^4 -コレステレン-3-オンの 240nm におけるミリモル分子吸光係数 ($\text{cm}^2 / \mu\text{mole}$)

3.10: 反応総液量 (ml)

0.05: 反応に供した酵素試料液量 (ml)

X : 酵素試料液の検品濃度 (mg/ml)