ADP-DEPENDENT HEXOKINASE [ADP-HKT II]

from *Thermococcus litoralis*
(ADP-specific glucokinase, EC 2.7.1.147)

Preparation and Specification

Appearance : White amorphous powder, lyophilized
Specific activity : More than 25 U/mg solid

Properties

Substrate specificity	: See table 1	
Molecular weight	: 50 kDa (gel filtration) 50 kDa (SDS-PAGE)	
Isoelectric point	: pH 4.1	
Michaelis constants	: Glucose 0.4 mM (at 37°C) ADP 0.057 mM (at 37°C)	
Optimum pH	: 7-7.5	Figure 1
pH stability	: 6.5-10.5 (37°C, 1 hr)	Figure 2
Optimum temperature	: 100°C	Figure 3
Thermal stability	: Stable at 95°C and below	Figure 4
Activators	: Mg ²⁺	
Effect of chemicals	: See Table 2	

Applications for Diagnostic Test

This enzyme is useful for enzymatic determination of 1, 5-Anhydroglucitol (1, 5-AG) when coupled with 1, 5-Anhydroglucitol-6-phosphate dehydrogenase (T-95).


Table 1. Substrate specificity

Substrate (20mM)	Relative activity (%)
D-Glucose	100
D-Mannose	13
D-Galactose	9
D-Fructose	2
2-deoxy-D-Glucose	4
D-Glucosamine	67
D-1,5-Anhydroglucitol	166
D-Sorbitol	0
D-Myoinocitol	0
D-Mannitol	0
Sucrose	0

Table 2. Effect of chemicals on ADP-HKTII activity

Chemical	Concentration	Relative activity (%)
None	-	100
NaF	1mM	107
NaN ₃	1mM	107
EDTA	10mM	0
Cholic acid	0.1%	97
Tween-20	0.1%	103
Tween-40	0.1%	104
Tween-60	0.1%	104
Tween-80	0.1%	103


Fig.1 pH Optimum

○ : Acetate buffer
 □ : Citrate buffer
 △ : Phosphate buffer
 ● : Tris-HCl buffer
 ■ : Glycine-NaOH buffer


Fig.2 pH Stability

37°C, 60 min.
 ○ : Acetate buffer
 □ : Citrate buffer
 △ : Phosphate buffer
 ● : Tris-HCl buffer
 ■ : Glycine-NaOH buffer


Fig.3 Optimum Temperature

pH 7.5
 50 mM Tris-HCl buffer



Fig.4 Thermal Stability

pH 7.5, 10min.
 50 mM Tris-HCl buffer

Assay

Principle

The assay is based on the increase in absorbance at 340 nm as the formation of NADPH proceeds in the following reaction.

Unit definition

One unit is defined as the amount of enzyme which converts 1 μ mole of D-glucose to D-Glucose-6-phosphate per minute at 37°C under the conditions specified in the assay procedure.

Reagents

- Reaction mixture

0.1M Tris-HCl Buffer pH7.5	1.50 ml
0.2M Glucose solution	0.30 ml

- | | |
|---------------------------------|---------|
| 20mM ADP solution | 0.30 ml |
| 20mM MgCl ₂ solution | 0.30 ml |
| 10mM NADP solution | 0.30 ml |
| 50U/ml G6PDH solution | 0.30 ml |
- Enzyme dilution buffer

10mM Tris-HCl Buffer pH7.5 (25°C)

 - Reagents

Tris (hydroxymethyl) aminomethane:Sigma Chemical Co. #T-1503
Glucose:Wako Pure Chemical Industries, Ltd. #041-00595
ADP (Adenosine diphosphate·2Na):Oriental yeast Co.,Ltd.
MgCl ₂ ·6H ₂ O:Wako Pure Chemical Industries, Ltd. #131-00162
NADP (Nicotinamide adenine dinucleotide oxidized form): Wako Pure Chemical Industries, Ltd. #308-50463
G6PDH II (Glucose-6-phosphate dehydrogenase) : Asahi Kasei Pharma Corporation #T-51

Enzyme solution

Accurately weigh about 20 mg of the sample and add enzyme dilution buffer to make a total of 20 ml. Dilute it with enzyme dilution buffer to adjust the concentration as required.

■ Procedure

- Pipette accurately 3.0 ml of reaction mixture into a small test tube and preincubate at 37°C.
- After 5 min. add accurately 50 μl of enzyme solution and mix to start the reaction at 37°C.
※ In the case of a test blank, add 50 μl of enzyme dilution buffer in place of enzyme solution.
- After starting the reaction, measure the rate of increase per minutes in absorbance at 340 nm. The rate must be measured within the linear portion of the absorbance curve. (Ex. Linear range from 2 min. to 6 min.)

Absorbance sample : As/min.

blank : Ab/min.

$$\Delta A/\text{min.} = (As/\text{min.} - Ab/\text{min.}) \leq 0.120 \text{ Abs/min.}$$

■ Calculation

$$\begin{aligned} \text{Activity (U/mg of powder)} &= \frac{\Delta A/\text{min.}}{6.3} \times \frac{3.05}{0.05} \times \frac{1}{X} \\ &= \Delta A/\text{min.} \times 9.68 \div X \end{aligned}$$

6.3 : millimolar extinction coefficient of NADPH at 340nm
($\text{cm}^2/\mu\text{mol}$)

3.05 : final volume (ml)

0.05 : volume of enzyme solution (ml)

X : concentration of the sample in enzyme solution
(mg/ml)

Storage

Storage at -20°C in the presence of a desiccant is recommended.

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