

## Acetaminophen 506 / 03R11-20

**Note:** Unless otherwise indicated, data was collected on a Roche/Hitachi® 717 analyzer.

<b>510 (K) Number:</b>	K081938
<b>Manufactured by</b>	Sekisui Diagnostics
<b>Test Principle:</b>	The enzyme, acyl amidohydrolase, cleaves the amide bond of the acetaminophen molecule, leaving p-aminophenol and acetate. The p-aminophenol is reacted with 2,5-dimethylphenol in the presence of manganese ions to form a colored compound, 4-(4-aminophenyl)-2,5-dimethylcyclohexadiene-1-one. The increased absorbance at 605 (660 <sup>4</sup> ) nm due to the formation of 4-(4-aminophenyl)-2,5-dimethylcyclohexadiene-1-one is directly proportional to the concentration of acetaminophen in the sample.
<b>Methodology</b>	Enzymatic/Colorimetric
<b>Sample Types:</b>	Fresh, clear, unhemolysed serum or lithium heparinized plasma. EDTA is not suitable for use.
<b>Fill Requirements:</b>	Use a minimum volume of 20ml of R2 reagent at a time, using only 20 ml wedges. When adding additional reagent to the analyzer use a new wedge. <sup>4</sup>
<b>On Board Stability</b>	8 days (192 hours) <sup>4</sup>
<b>Calibration Stability</b>	24 hours <sup>4</sup>
<b>Precision</b>	<b>Within-run:</b> ≤ 1.5% <b>Total Precision:</b> ≤ 2.9%
<b>Accuracy</b>	<b>Serum<sup>1</sup></b> Slope: 1.064 Intercept: 1.1 µg/mL (7.0 µmol/L) Correlation Coefficient: 0.9998  <b>Plasma<sup>3</sup></b> Slope: 0.999 Intercept -0.3 µg/mL (-2.2 µmol/L) Correlation Coefficient: 0.9999
<b>Linearity</b>	0.6 – 377.5 ug/mL (4 – 2500 µmol/L)
<b>No Significant Interference to levels indicated</b>  (See insert/IFU for complete listing)	<ul style="list-style-type: none"> <li>▪ N-Acetylcysteine: 1500 mg/L (9.2 mmol/L)</li> <li>▪ Hemoglobin: 200 mg/dL (31 µmol/L)</li> <li>▪ Conjugated Bilirubin: 2 mg/dL (23.7 µmol/L)</li> <li>▪ Unconjugated Bilirubin: 2 mg/dL (34.2 µmol/L)</li> <li>▪ Ascorbic Acid: 3000 µg/dL (170 µmol/L)</li> <li>▪ Intralipid: 200 mg/dL (600 mg/dL Simulated Triglyceride)</li> </ul>

<sup>1</sup> SERUM: The performance of this method (y) was compared with the performance with a similar acetaminophen method (x) on a Roche/Hitachi® 717 analyzer.

<sup>2</sup> PLASMA: The performance of this method with plasma (y) was compared to the performance of this method with serum (x) on an Advia® 1650 analyzer.

<sup>3</sup> PLASMA: The performance of this method with plasma (y) was compared to the performance of this method with serum (x) on a Roche/Hitachi® 717 analyzer

<sup>4</sup> Testing completed on Architect c8000 system