

CLSI EP9 Method Comparison

X Method: ACE506 (Hit 717)

Y Method: ACE506 c8000

Regression Analysis

	Deming	Regular
Slope:	0.996 (0.993 to 0.998)	0.995 (0.993 to 0.998)
Intercept:	-3.8 (-6.8 to -0.9)	-3.6 (-6.6 to -0.6)
Std Err Est:	12.5	12.5

95% Confidence Intervals are shown in parentheses

Preliminary Data Examination

Within-method outlier analysis	PASS	Acceptability Limits: X: 40 or 5.2%, Y: 13 or 1.7% Number of outliers: 0 of 101 pairs
Between-method outlier analysis	PASS	Acceptability Limits: 45 or 5.8% Number of outliers: 0 of 101 pairs
Visual check for linear relationship	_____	SubRange Bounds: None
Test for adequate number of results	PASS	N = 101 pairs
Test for adequate range of results	PASS	R = 0.9999
Visual check for uniform scatter	FAIL	Computed precision ratio of 6.09 is outside of normal limits

User's Specifications

Medically Allowable Error	24.0 umol/L or 30.0%	SubRange Bounds: None Scatter Plot Bounds: None
---------------------------	----------------------	--

Experiment Description

	X Method	Y Method
Date:	26 Jun 2015	26 Jun 2015
Rep SD:	10.202	3.2889
Result Ranges:	8 to 2481	8 to 2482
Units:	umol/L	umol/L
Comment:		
Analyst:	Lorna	Lorna

Accepted by: _____

Signature



_____ Date

CLSI EP9 Method Comparison



X Method: ACE506 (Hit 717)

Y Method: ACE506 c8000

Individual Results Plotted



Means of Replicate Results Plotted



CLSI EP9 Method Comparison

X Method: ACE506 (Hit 717)

Y Method: ACE506 c8000

Experimental Results

Spec ID	Y (Test)		X (Comp)		Bias (Y-X)	Calc'd		Spec ID	Y (Test)		X (Comp)		Bias (Y-X)	Calc'd	
	Result	Diff	Result	Diff		Y	Resid		Result	Diff	Result	Diff		Y	Resid
S00001	23		20		3	16.1	6.9		1355	9	1389	-19	-34	1379.2	-24.2
	22	1	23	-3	-1	19.1	2.9	S00029	1200		1220		-20	1210.9	-10.9
S00002	307		302		5	296.9	10.1		1197	3	1227	-7	-30	1217.9	-20.9
	306	1	302	0	4	296.9	9.1	S00030	496		506		-10	500.0	-4.0
S00003	1325		1331		-6	1321.4	3.6		497	-1	510	-4	-13	504.0	-7.0
	1328	-3	1346	-15	-18	1336.4	-8.4	S00031	760		770		-10	762.8	-2.8
S00004	1085		1085		0	1076.5	8.5		759	1	768	2	-9	760.9	-1.9
	1089	-4	1089	-4	0	1080.5	8.5	S00032	1743		1749		-6	1737.6	5.4
S00005	503		502		1	496.0	7.0		1743	0	1779	-30	-36	1767.5	-24.5
	503	0	512	-10	-9	506.0	-3.0	S00033	2332		2329		3	2315.1	16.9
S00006	762		745		17	738.0	24.0		2338	-6	2355	-26	-17	2341.0	-3.0
	753	9	752	-7	1	744.9	8.1	S00034	2043		2083		-40	2070.2	-27.2
S00007	1722		1705		17	1693.8	28.2		2038	5	2077	6	-39	2064.2	-26.2
	1721	1	1733	-28	-12	1721.7	-0.7	S00035	1581		1612		-31	1601.2	-20.2
S00008	2310		2302		8	2288.2	21.8		1587	-6	1617	-5	-30	1606.2	-19.2
	2304	6	2305	-3	-1	2291.2	12.8	S00036	425		432		-7	426.3	-1.3
S00009	1994		2010		-16	1997.5	-3.5		425	0	429	3	-4	423.3	1.7
	1996	-2	1998	12	-2	1985.6	10.4	S00037	91		96		-5	91.7	-0.7
S00010	1575		1564		11	1553.4	21.6		91	0	94	2	-3	89.8	1.2
	1570	5	1577	-13	-7	1566.4	3.6	S00038	1181		1191		-10	1182.0	-1.0
S00011	416		413		3	407.4	8.6		1182	-1	1200	-9	-18	1191.0	-9.0
	415	1	411	2	4	405.4	9.6	S00039	1424		1452		-28	1441.9	-17.9
S00012	69		68		1	63.9	5.1		1426	-2	1448	4	-22	1437.9	-11.9
	68	1	67	1	1	62.9	5.1	S00040	596		607		-11	600.5	-4.5
S00013	1135		1152		-17	1143.2	-8.2		599	-3	616	-9	-17	609.5	-10.5
	1153	-18	1160	-8	-7	1151.2	1.8	S00041	406		406		0	400.4	5.6
S00014	1405		1412		-7	1402.1	2.9		406	0	404	2	2	398.4	7.6
	1402	3	1409	3	-7	1399.1	2.9	S00042	2404		2400		4	2385.8	18.2
S00015	598		599		-1	592.6	5.4		2410	-6	2400	0	10	2385.8	24.2
	600	-2	602	-3	-2	595.6	4.4	S00043	174		177		-3	172.4	1.6
S00016	415		414		1	408.4	6.6		175	-1	176	1	-1	171.4	3.6
	413	2	413	1	0	407.4	5.6	S00044	324		326		-2	320.8	3.2
S00017	2392		2375		17	2360.9	31.1		324	0	328	-2	-4	322.7	1.3
	2384	8	2403	-28	-19	2388.8	-4.8	S00045	X 2481		2499				
S00018	173		172		1	167.4	5.6		X 2488		2511				
	174	-1	178	-6	-4	173.4	0.6	S00046	1252		1247		5	1237.8	14.2
S00019	318		317		1	311.8	6.2		1256	-4	1251	-4	5	1241.8	14.2
	316	2	321	-4	-5	315.8	0.2	S00047	378		378		0	372.5	5.5
S00020	X 2523		2546						377	1	380	-2	-3	374.5	2.5
	X 2504		2531					S00048	664		665		-1	658.3	5.7
S00021	1252		1259		-7	1249.7	2.3		665	-1	680	-15	-15	673.2	-8.2
	1246	6	1260	-1	-14	1250.7	-4.7	S00049	2115		2135		-20	2122.0	-7.0
S00022	373		379		-6	373.5	-0.5		2112	3	2128	7	-16	2115.0	-3.0
	374	-1	381	-2	-7	375.5	-1.5	S00050	749		754		-5	746.9	2.1
S00023	658		661		-3	654.3	3.7		749	0	756	-2	-7	748.9	0.1
	659	-1	673	-12	-14	666.3	-7.3	S00051	8		8		0	4.1	3.9
S00024	2087		2105		-18	2092.1	-5.1		8	0	8	0	0	4.1	3.9
	2082	5	2121	-16	-39	2108.0	-26.0	S00052	272		273		-1	268.0	4.0
S00025	738		753		-15	745.9	-7.9		274	-2	282	-9	-8	276.9	-2.9
	737	1	749	4	-12	741.9	-4.9	S00053	1288		1290		-2	1280.6	7.4
S00026	27		29		-2	25.0	2.0		1287	1	1301	-11	-14	1291.6	-4.6
	27	0	31	-2	-4	27.0	0.0	S00054	1079		1087		-8	1078.5	0.5
S00027	297		301		-4	295.9	1.1		1078	1	1091	-4	-13	1082.5	-4.5
	296	1	303	-2	-7	297.9	-1.9	S00055	479		486		-7	480.1	-1.1
S00028	1364		1370		-6	1360.3	3.7		476	3	485	1	-9	479.1	-3.1

Calculated values and residuals from Deming regression. Values marked with an "X" were excluded from the calculations. Outliers "O" were also excluded.

CLSI EP9 Method Comparison

X Method: ACE506 (Hit 717)

Y Method: ACE506 c8000

Experimental Results

Spec ID	Y (Test)		X (Comp)		Bias (Y-X)	Calc'd		Spec ID	Y (Test)		X (Comp)		Bias (Y-X)	Calc'd	
	Result	Diff	Result	Diff		Y	Resid		Result	Diff	Result	Diff		Y	Resid
S00056	736		748		-12	740.9	-4.9		769	2	766	12	3	758.9	10.1
	732	4	754	-6	-22	746.9	-14.9	S00082	1809		1816		-7	1804.3	4.7
S00057	1767		1735		32	1723.7	43.3		1814	-5	1810	6	4	1798.4	15.6
	1763	4	1767	-32	-4	1755.6	7.4	S00083	2387		2368		19	2354.0	33.0
S00058	2277		2289		-12	2275.3	1.7		2375	12	2409	-41	-34	2394.8	-19.8
	2268	9	2310	-21	-42	2296.2	-28.2	S00084	2070		2073		-3	2060.2	9.8
S00059	2167		2154		13	2140.9	26.1		2070	0	2080	-7	-10	2067.2	2.8
	2158	9	2206	-52	-48	2192.7	-34.7	S00085	1635		1659		-24	1648.0	-13.0
S00060	1528		1541		-13	1530.5	-2.5		1635	0	1645	14	-10	1634.1	0.9
	1529	-1	1539	2	-10	1528.5	0.5	S00086	442		445		-3	439.2	2.8
S00061	403		402		1	396.4	6.6		437	5	443	2	-6	437.3	-0.3
	402	1	407	-5	-5	401.4	0.6	S00087	95		94		1	89.8	5.2
S00062	52		47		5	43.0	9.0		94	1	98	-4	-4	93.7	0.3
	51	1	50	-3	1	45.9	5.1	S00088	1225		1226		-1	1216.9	8.1
S00063	1097		1103		-6	1094.4	2.6		1220	5	1226	0	-6	1216.9	3.1
	1100	-3	1125	-22	-25	1116.3	-16.3	S00089	1455		1466		-11	1455.8	-0.8
S00064	1347		1362		-15	1352.3	-5.3		1447	8	1472	-6	-25	1461.8	-14.8
	1348	-1	1360	2	-12	1350.3	-2.3	S00090	624		640		-16	633.4	-9.4
S00065	572		581		-9	574.7	-2.7		626	-2	632	8	-6	625.4	0.6
	571	1	582	-1	-11	575.7	-4.7	S00091	434		435		-1	429.3	4.7
S00066	393		392		1	386.5	6.5		428	6	438	-3	-10	432.3	-4.3
	391	2	397	-5	-6	391.5	-0.5	S00092	2482		2481		1	2466.5	15.5
S00067	2355		2320		35	2306.2	48.8		2476	6	2465	16	11	2450.5	25.5
	2340	15	2339	-19	1	2325.1	14.9	S00093	195		192		3	187.3	7.7
S00068	162		162		0	157.5	4.5		195	0	197	-5	-2	192.3	2.7
	163	-1	169	-7	-6	164.4	-1.4	S00094	347		349		-2	343.7	3.3
S00069	305		308		-3	302.8	2.2		347	0	355	-6	-8	349.6	-2.6
	304	1	310	-2	-6	304.8	-0.8	S00095	X 2466		2517				
S00070	2446		2475		-29	2460.5	-14.5		X 2472		2497				
	2436	10	2454	21	-18	2439.6	-3.6	S00096	1332		1330		2	1320.4	11.6
S00071	1219		1237		-18	1227.8	-8.8		1329	3	1332	-2	-3	1322.4	6.6
	1216	3	1257	-20	-41	1247.7	-31.7	S00097	396		399		-3	393.4	2.6
S00072	356		365		-9	359.6	-3.6		389	7	402	-3	-13	396.4	-7.4
	356	0	361	4	-5	355.6	0.4	S00098	679		686		-7	679.2	-0.2
S00073	647		660		-13	653.3	-6.3		684	-5	697	-11	-13	690.2	-6.2
	646	1	669	-9	-23	662.3	-16.3	S00099	2179		2211		-32	2197.6	-18.6
S00074	2062		2098		-36	2085.1	-23.1		2178	1	2189	22	-11	2175.7	2.3
	2057	5	2063	35	-6	2050.3	6.7	S00100	769		777		-8	769.8	-0.8
S00075	733		739		-6	732.0	1.0		764	5	783	-6	-19	775.8	-11.8
	734	-1	757	-18	-23	749.9	-15.9	S00101	X 2537		2555				
S00076	42		40		2	36.0	6.0		X 2525		2531				
	41	1	42	-2	-1	38.0	3.0	S00102	2123		2115		8	2102.1	20.9
S00077	316		319		-3	313.8	2.2		2124	-1	2151	-36	-27	2137.9	-13.9
	314	2	324	-5	-10	318.8	-4.8	S00103	2243		2246		-3	2232.5	10.5
S00078	1404		1402		2	1392.1	11.9		2231	12	2258	-12	-27	2244.4	-13.4
	1399	5	1434	-32	-35	1424.0	-25.0	S00104	2108		2141		-33	2127.9	-19.9
S00079	1239		1248		-9	1238.8	0.2		2108	0	2129	12	-21	2116.0	-8.0
	1241	-2	1284	-36	-43	1274.6	-33.6	S00105	2180		2192		-12	2178.7	1.3
S00080	515		532		-17	525.9	-10.9		2185	-5	2216	-24	-31	2202.6	-17.6
	515	0	529	3	-14	522.9	-7.9								
S00081	771		778		-7	770.8	0.2								

Calculated values and residuals from Deming regression. Values marked with an "X" were excluded from the calculations. Outliers "O" were also excluded.