



ARK™ Zonisamide Calibrator











This ARK Diagnostics, Inc. package insert for the ARK Zonisamide Calibrator must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of the assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Customer Service


 48089 Fremont Blvd
 Fremont, CA 94538 USA
 Tel: 1-877-869-2320
 Fax: 1-510-270-6298
 customersupport@ark-tdm.com
 www.ark-tdm.com

 
 Emergo Europe
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands


Key to Symbols Used

	Batch code	 YYYY-MM-DD	Use by/Expiration date
	Catalog Number		Manufacturer
	Authorized Representative		CE Mark
	In Vitro Diagnostic Medical Device		Temperature limitation
	Consult Instructions for Use		Calibrator
Rx Only	For Prescription Use Only		

1 Name

ARK™ Zonisamide Calibrator

2 Intended Use

ARK Zonisamide Calibrator is intended for use in calibration of the ARK Zonisamide Assay,  5022-0001-00.

3 Content

ARK Zonisamide Calibrator is comprised of a synthetic protein matrix with the following concentrations of zonisamide:

REF	Product Description	Quantity/Volume
5022-0002-00	ARK Zonisamide Calibrators Zonisamide, buffer, bovine serum albumin, and preservatives	Dropper vials
	A	0.0 µg/mL
	B	2.5 µg/mL
	C	7.5 µg/mL
	D	20.0 µg/mL
	E	40.0 µg/mL
	F	80.0 µg/mL

4 Standardization

There is no internationally recognized standard for zonisamide. ARK Zonisamide Calibrators are prepared by gravimetric dilution of high purity zonisamide into a synthetic proteinaceous matrix free of zonisamide.

5 Warnings and Precautions

- For *In Vitro Diagnostic* Use.
- Do not mix calibrators from different lot numbers.
- Use each lot as a set.

6 Instructions For Use

- For a complete summary and explanation of the zonisamide assay, refer to the package insert for the ARK Zonisamide Assay.
- Calibrators are ready to use. Mix each level by gentle inversion before dispensing.
- Squeeze sufficient volume (~40µL/drop) into individual sample cups for each level. Consult instrument-specific sample volume requirements. Return caps to their original containers and keep tight.
- Store at 2-8°C. Use prior to the expiration date.

7 Procedure

Calibration

Perform a full calibration (6- point) procedure; test calibrators in duplicate. Verify the calibration curve with at least two levels of quality controls according to the established laboratory quality assurance plan.

When to Re-Calibrate

- Whenever a new lot number of reagents is used
- Whenever indicated by quality control results
- Whenever required by standard laboratory protocols

8 Limitations of Procedure

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, calibrators, controls, storage of product as directed, and good laboratory technique.

9 Trademarks

ARK™ is a trademark of ARK Diagnostics, Inc. Other brand or product names are trademarks of their respective holders.