EN	1
BG	7
DA	13
NL	18
FI	24
FR	29
DE	35
IT	41
LT	47
PL	52
PT	58
RO	64
RU	69
ES	75
SV	81
TR	86

REFERENCES / СПРАВКИ / LITTERATUR / BIBLIOGRAFIE / VIITTAUKSET / RÉFÉRENCES / REFERENZEN / BIBLIOGRAFIA / LITERATÜRA / PIŚMIENNICTWO / BIBLIOGRAFIA / REFERİNTE / ССЫЛКИ / REFERENCIAS / REFERENSER / REFERANSLAR.....	93
SYMBOLS / СИМВОЛИ / SYMBOLER / SYMBOLEN / MERKIT / SYMBOLES / SYMBOLE / SIMBOLI / ŽENKLAI / SYMBOLE / SÍMBOLOS / SIMBOLURI / УСЛОВНЫЕ ОБОЗНАЧЕНИЯ / SÍMBOLOS / SYMBOLER / SEMBOLLER	94

EN

OSOM® Strep A Test

Catalog number 141E or 141E-20

CLIA Complexity: Waived

The user should read the instructions for use before performing the test.

INTENDED USE

The OSOM Strep A Test is a manual, rapid chromatographic assay intended for the qualitative detection of Group A Streptococcal antigen from throat swabs from patients presenting with pharyngitis as an aid in the diagnosis of Group A Streptococcal infections. For near-patient, laboratory professional, and healthcare professional *in vitro* diagnostic use only. Additionally, the OSOM Strep A test is intended for use for confirmation of presumptive Group A Streptococcal colonies recovered from culture, for laboratory professional *in vitro* diagnostic use only.

SUMMARY AND EXPLANATION OF TEST

Group A *Streptococcus* is one of the most important causes of acute upper respiratory tract infection. Early diagnosis and treatment of Group A Streptococcal pharyngitis has been shown to reduce the severity of symptoms of patients suspected of Strep A infection and further complications such as rheumatic fever and glomerulonephritis.¹ Conventional identification procedures for Group A *Streptococcus* from throat swabs involve the isolation and subsequent identification of viable pathogens by techniques that require 24 to 48 hours or longer.² The OSOM Strep A Test is a manual (aid to diagnostic) test that detects either viable or nonviable organisms directly from a throat swab, providing results within 5 minutes.

PRINCIPLES OF TEST

The OSOM Strep A Test uses color immunochromatographic dipstick technology with rabbit antibodies coated on the nitrocellulose membrane. In the test procedure, a throat swab is subjected to a chemical extraction of a carbohydrate antigen unique to Group A *Streptococcus*. The Test Stick is then placed in the extraction mixture and the mixture migrates along the membrane. If Group A *Streptococcus* is present in the sample, it will form a complex with the anti-Group A *Streptococcus* antibody conjugated color particles. The complex will then be bound by the anti-Group A *Streptococcus* capture antibody and a visible blue Test Line will appear to indicate a positive result.

KIT CONTENTS AND STORAGE

50/20 Test Sticks

50/20 Test Tubes

50/20 Sterile Swabs

1 Strep A Reagent 1 (Strep A REAG 1) (2 M Sodium Nitrite)

1 Strep A Reagent 2 (Strep A REAG 2) (0.3 M Acetic Acid)

1 Positive Control (Nonviable Group A Streptococci, 0.1% Sodium Azide)

1 Negative Control (Nonviable Group C Streptococci, 0.1% Sodium Azide)

1 Instructions for Use

Note: Two extra test sticks have been included in the kit for external QC testing. In addition, extra components (swabs, tubes) have been provided for your convenience.

Store Test Sticks and reagents tightly capped at 15-30°C (59-86°F).

Do not use Test Sticks or reagents after expiration date.

Components are not intended for reuse.

MATERIALS REQUIRED BUT NOT PROVIDED

A timer or watch

WARNINGS AND PRECAUTIONS

Component (s)	Pictogram	Signal Word	Hazardous Ingredients
Strep A REAG 1		Warning	sodium nitrite (CAS No) 7632-00-0
Hazard statements	H302 - Harmful if swallowed.		
Precautionary statements	P264 - Wash hands, forearms and face thoroughly after handling. P270 - Do not eat, drink or smoke when using this product. P301+P312 - IF SWALLOWED: Call a POISON CENTRE or doctor if you feel unwell. P330 - Rinse mouth. P501 - Dispose of contents and container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.		
Component (s)	Pictogram	Signal Word	Hazardous Ingredients
Strep A REAG 2		Warning	Acetic acid (CAS No) 64-19-7
Hazard statements	H315 - Causes skin irritation. H319 - Causes serious eye irritation.		
Precautionary statements	P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection. P302+P352 - IF ON SKIN: Wash with plenty of water. P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313 - If skin irritation occurs: Get medical advice/attention. P337+P313 - If eye irritation persists: Get medical advice/attention. P362+P364 - Take off contaminated clothing and wash it before reuse.		

For *in vitro* diagnostic use.

Follow your laboratory safety guidelines in the collection, handling, storage and disposal of controls, patient specimens and all items exposed to patient specimens.³

Dispose of all used reagents and any other contaminated disposable materials following procedures for infectious or potentially infectious waste. It is the responsibility of each laboratory to handle solid and liquid waste according to their nature and degree of hazard and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

The Positive and Negative Controls contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azide. For sites permitted to dispose of material down a sink: large quantities of water must be used to flush discarded control material down a sink.

Test Stick Container should remain capped when not in use. Test Sticks are sensitive to prolonged exposure to moisture.

Do not interchange or mix components from different lots.

Do not use the sterile swab if packaging has been breached.

SPECIMEN COLLECTION AND PREPARATION

- Collect specimens with a sterile swab from the tonsils and/or the back of the throat² taking care to avoid the teeth, gums, tongue or cheek surfaces.
- Do not use swabs with cotton tips, wooden shafts or calcium alginate swabs.
- Do not use a collection system that contains charcoal or semisolid transport media.
- If your lab requires a culture result as well as the OSOM Strep A Test result, streak the culture plate with the swab before starting the OSOM Strep A Test procedure as the extraction reagents will cause the specimen to become nonviable.
- Process the swab as soon as possible after collecting the specimen. If you do not perform the OSOM Strep A test immediately, store the swabs either at room temperature or refrigerated for up to 72 hours. The swabs and the test kit must be at room temperature prior to running the test.
- Sample Transport:
 - Because the performance characteristics of this product were established with the sterile rayon swabs supplied with the kit, we recommend using these swabs to assure optimal performance. You may purchase the kit swabs in a double swab/dry tube format as an accessory (SEKISUI Diagnostics Part #7784).
 - Because the test does not require live organisms for processing, a rayon transport swab containing Stuart's or Amies media may also be used; however, swabs from other suppliers have not been validated.

CULTURE CONFIRMATION

The OSOM Strep A Test can also be used to confirm the identification of Group A *Streptococcus* on blood agar plates. The plates must be less than 72 hours old. Lightly touch 1-3 suspect colonies (showing characteristic beta hemolysis) using a sterile swab. Do not sweep the plate. Follow the instructions in the TEST PROCEDURE section to test the swab.

QUALITY CONTROL

Internal Procedural Controls

The OSOM Strep A Test provides three levels of procedural controls with each test run.

- The color of the liquid changes from pink to light yellow as you add Extraction Reagent 2 to Extraction Reagent 1. This is an internal extraction reagent control. The color change means that you mixed the extraction reagents properly. The color change also means that the reagents are functioning properly.
- The red Control Line is an internal positive procedural control. The Test Stick must absorb the proper amount of sample and the Test Stick must be working properly for the red Control Line to appear. For the Test Stick to be working properly, the capillary flow must occur.
- A clear background is an internal background negative procedural control. If no interfering substances are in the specimen and the Test Stick is working properly, the background in the Control Line area will clear. A discernible result will be seen.

If the red Control Line does not appear, the test may be invalid. If the background does not clear and interferes with the test result, the test may be invalid.

External Quality Control Testing

Each kit contains Positive and Negative Control material. The Controls are for external quality control testing. Use the Controls to test that the extraction reagents and the Test Sticks are working. Also use the Controls to test that you are able to correctly perform the test procedure. If you choose, you may use Group A and non Group A *Streptococcus* ATCC reference strains as controls. Some commercial controls may contain interfering additives. Therefore it is recommended that you do not use other commercial controls with the OSOM Strep A Test.

Quality Control requirements should be established in accordance with local regulations or accreditation requirements. Minimally, SEKISUI Diagnostics recommends that positive and negative external controls be run with each new lot and with each new untrained operator.

QC Testing Procedure:

- Dispense 3 drops Reagent 1 and 3 drops Reagent 2 into Test Tube.
- Vigorously mix the control contents. Add 1 free falling drop of Control from dropper bottle.
- Place a clean swab into the Tube.
- Continue as you would for a patient sample, as instructed in the PROCEDURE section.

LIMITATIONS

- The OSOM Strep A Test has been categorized as CLIA waived only for the application of qualitative detection of Group A Streptococcal Antigen from throat swabs. The application for the confirmation of presumptive Group A Streptococcal colonies recovered from culture is not waived.
- The results obtained with this kit yield data that must be used only as an adjunct to other information available to the physician. The OSOM Strep A Test is a qualitative test for the detection of Group A Streptococcal antigen. This test does not differentiate between viable and nonviable Group A Streptococci.
- The OSOM Strep A Test should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
- This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A *Streptococcus* (GAS).^{1,2}
- A negative result may be obtained if the specimen is inadequate or antigen concentration is below the sensitivity of the test.
- The American Academy of Pediatrics states "Several rapid diagnostic tests for GAS pharyngitis are available. The specificities of these tests generally are very high, but the reported sensitivities vary considerably. As with throat cultures, the accuracy of these tests is most dependent on the quality of the throat swab specimen, which must contain pharyngeal and tonsillar secretions, and on the experience of the person who is performing the test. Therefore, when a patient suspected of having GAS pharyngitis has a negative rapid streptococcal test a throat culture should be obtained to ensure that the patient does not have GAS infection." It also states; "Cultures that are negative for GAS infection after 24 hours should be incubated for a second day to optimize isolation of GAS."⁴

EXPECTED RESULTS

Approximately 19% of all upper respiratory tract infections are caused by Group A Streptococci.⁵ Streptococcal pharyngitis displays a seasonal variation and is most prevalent during winter and early spring. The highest incidence of this disease is found in crowded populations such as in school-age children.⁶

PERFORMANCE CHARACTERISTICS

In a multi-center evaluation, a total of 639 throat swabs were collected from patients presenting with pharyngitis. Each swab was inoculated to a sheep blood agar plate, then tested by the OSOM Strep A Test. Plates were incubated for 18-24 hours at 35-37°C at 5-10% CO₂, with a Bacitracin disk. Presumptive GAS colonies were confirmed with commercially available Strep A testing kits.

Of the 639 total specimens, 464 were found to be negative by culture and 454 were also negative by the OSOM Strep A Test, for a specificity of 97.8%. Of the 175 specimens found to be positive by culture, 168 were also positive by the OSOM Strep A Test for a sensitivity of 96.0%. The 95% confidence intervals were calculated to be 96.6-99.0% for specificity and 94.4-97.6% for sensitivity. Overall agreement between culture and the OSOM Strep A Test was 97.3% (622/639). The results are summarized below:

Culture Classification	OSOM/Culture	% Correct
Negative (Specificity)	454/464	97.8%
1+ (\leq 10 colonies)	3/6	50.0%
2+ (11-50 colonies)	9/13	69.2%
3+ ($>$ 50 colonies)	44/44	100%
4+ (predominant growth)	112/112	100%
Total Positive (Sensitivity)	168/175	96.0%
Total (Overall Agreement)	622/639	97.3%

In addition, the OSOM Strep A Test was used to confirm the identification of Group A *Streptococcus* on blood agar plates. As a culture confirmation test, the Strep A Test was 100% sensitive (62/62) and 100% specific (39/39).

TEST PROCEDURE


Absorbent End

Result Window

Handle End

- Just before testing, add 3 drops Reagent 1 (pink) and 3 drops Reagent 2 to the Test Tube (the solution should turn light yellow).
- Immediately put the swab into the Tube.
- Vigorously mix the solution by rotating the swab forcefully against the side of the Tube at least ten (10) times. Best results are obtained when the specimen is vigorously extracted in the solution.
- Let stand for 1 minute.
- Express as much liquid as possible from the swab by squeezing the sides of the tube as the swab is withdrawn.
- Discard the swab.
- Remove Test Stick(s) from the container; recap container immediately.
- Place the Absorbent End of the Test Stick into the extracted sample.
- **Read results at 5 minutes. Positive results may be read as soon as the red Control Line appears.**
- Results are invalid after the stated read time. The use of a timer is recommended.


INTERPRETATION OF TEST RESULTS

Notes


A blue or red line which appears uneven in color density is considered a valid result. In cases of moderate or high positive specimens, some blue color behind the Test Line may be seen; as long as the Test Line and Control Line are visible, the results are valid. Ensure adequate lighting conditions for viewing results.

Positive

A blue Test Line and a red Control Line is a positive result for the detection of Group A Streptococcus antigen. Note that the blue line can be any shade of blue and can be lighter or darker than the line in the picture.

Negative

A red Control Line but no blue Test Line is a presumptive negative result.

Invalid

If no red Control Line appears or background color makes reading the red Control Line impossible, the result is invalid. If this occurs, repeat the test on a new Test Stick.

ANALYTICAL PERFORMANCE

Analytical Sensitivity

The limit of detection (LoD) for the OSOM Strep A Test was established via dilution experiments that were performed with three reagent lots and two strains of *S. pyogenes* (ATCC strain 19615 and ATCC strain 12363). The *S. pyogenes* were introduced onto swabs, which were then extracted and tested. The LoD was identified as the most dilute concentration tested that produced positive results in both replicates for each of the three lots. The results of the LoD experiment are summarized in the table below.

S. pyogenes strain	Limit of Detection (LoD)
ATCC strain 19615	8 x 10 ⁴ organisms/test
ATCC strain 12363	4 x 10 ⁴ organisms/test

Cross-reactivity

The following organisms tested at levels of approximately 1x10⁸ organisms/test were all found to be negative when tested with the OSOM Strep A Test.

Group B <i>Streptococcus</i>	<i>Enterococcus faecalis</i>	<i>Pseudomonas aeruginosa</i>
Group C <i>Streptococcus</i>	<i>Staphylococcus aureus</i>	<i>Bordetella pertussis</i>
Group D <i>Streptococcus</i>	<i>Staphylococcus epidermidis</i>	<i>Neisseria meningitidis</i>
Group F <i>Streptococcus</i>	<i>Corynebacterium diphtheriae</i>	<i>Neisseria gonorrhoeae</i>
Group G <i>Streptococcus</i>	<i>Serratia marcescens</i>	<i>Neisseria sicca</i>
<i>Streptococcus pneumoniae</i>	<i>Candida albicans</i>	<i>Neisseria subflava</i>
<i>Streptococcus sanguis</i>	<i>Klebsiella pneumoniae</i>	<i>Branhamella catarrhalis</i>
<i>Streptococcus mutans</i>	<i>Escherichia coli</i>	<i>Hemophilus influenza</i>
<i>Staphylococcus aureas Cowan strain</i>		

POL Studies

An evaluation of the OSOM Strep A Test was conducted at three physicians offices where testing was performed by personnel with diverse educational backgrounds. Each site tested the randomly coded panel consisting of negative (6), low positive (3) and moderate positive (3) specimens for three days. The results obtained had >99% agreement (107/108) with the expected results.

REORDER

No. 141E - OSOM Strep A Test (50 Tests)

No. 141E-20 - OSOM Strep A Test (20 Tests)

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on *In vitro* Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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1. Youmans, G.P., Paterson, P.Y. and Sommers, H.M. Upper Respiratory Tract Infection: General Considerations, in The Biologic and Clinical Basis of Infectious Diseases, W. B. Saunders Co., Philadelphia, 177-183, 1980.
2. Faklam, R.R. and Washington, J.A., Streptococcus and Related Catalase-Negative Gram-Positive Coccii, in Manual of Clinical Microbiology, 5th Edition, Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D., and Shadomy, H.J., Eds., Am. Society of Microbiology, Washington, D.C., 238-257, 1991.
3. CDC, Biosafety in Microbiological and Biomedical Laboratories, 2nd Ed., HHS Publication No. 8808395, 4-6, 1988.
4. American Academy of Pediatrics. Summaries of Infectious Disease. In: Pickering LK, ed. 2000 Red Book: Report of the Committee on Infectious Diseases. 25th ed. Elk Grove Village IL: American Academy of Pediatrics; 2000 528.
5. Lauer, B.A., Reller, L.B., and Mirrett, S., Effect of Atmosphere and Duration of Incubation on Primary Isolation of Group A Streptococci from Throat Cultures, *J. Clin. Microb.*, 17: 338-340, 1983.
6. Wannamaker, L.W., Differences Between Streptococcal Infections of the Throat and of the Skin, *N. Eng. J. Med.*, 282: 23-31, 78-85, 1970.

LOT

REF

Manufacturer	Batch code	Catalogue number
Производител	Партиден номер	Каталожен номер
Producent	Batchkode	Katalognummer
Fabrikant	Batchcode	Catalogusnummer
Valmistaja	Eränumero	Luettelon numero
Fabricant	Code de lot	Référence de catalogue
Hersteller	Chargen-Code	Katalognummer
Fabbricante	Codice del lotto	Numero di catalogo
Gamintojas	Partijos kodas	Katalogo numeris
Producent	Kod partii	Numer katalogowy
Fabricante	Código de lote	Número de referência
Producător	Nr. lot	Nr. catalog
Производитель	Номер партии	Номер по каталогу
Fabricante	Código de lote	Número de catálogo
Tillverkare	Batchkod	Katalognummer
Üretici	Seri kodu	Katalog numarası

Rx ONLY

Caution: Federal law restricts this device to sale by or on the order of a physician

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OSOM and the QC Inside logo are registered U.S. trademarks of SEKISUI Diagnostics, LLC.
Licensed under U.S. Patent Nos. 5,714,389; 5,989,921; 6,485,982 and 6,979,576 and related non-U.S. patents and patent applications.
OSOM och QC Inside-logotypen är registrerade varumärken i USA som tillhör SEKISUI Diagnostics, LLC.
OSOM ve QC Inside logosu SEKISUI Diagnostics, LLC'nin tescilli ABD ticari markalarıdır.
SEKISUI Diagnostics, LLC
MDSS GmbH