

The FastPack® IP αGST Immunoassay (International Only)



FAST / SIMPLE / ACCURATE

The FastPack® IP αGST Immunoassay

The FastPack® IP αGST Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 12 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



SEKISUI
DIAGNOSTICS

Because every result matters™

The FastPack® IP αGST Immunoassay (International Only)

Simple Process



1
Dispense sample into the FastPack® Injection Port



2
Place FastPack® into the FastPack® IP Instrument



3
Press the START button

The FastPack® IP Free αGST Assay Specifications:

- Quantitative measurement of Alpha Glutathione S-Transferase
- Lower Limit of 0.5 ng/mL
- Upper Limit of 200 ng/mL; up to 500,000 ng/mL with dilution
- Assay time: 12 minutes

ORDERING INFORMATION	CONFIGURATION	CATALOG NUMBER
FastPack® IP αGST (International Only)	30 Tests	25000047
FastPack® IP αGST (Research Use Only – US)	30 Tests	25000051
FastPack® IP αGST Diluent	30/Kit	25000050
FastPack® IP αGST Diluent (Research Use Only – US)	30/Kit	25000054
FastPack® αGST Calibrator	Kit	25000048
FastPack® αGST Calibrator (Research Use Only – US)	Kit	25000053
FastPack® αGST Control	Kit	25000049
FastPack® αGST Control (Research Use Only – US)	Kit	25000052

THE AMERICAS

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent ME16 0LQ, UK
Phone: +44 1622 607800
Fax: +44 1622 607801
Email: info@sekisui-dx.com

www.sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS

Because every result matters™