


Automated Coagulation System CP3000

Offering the benefits of a **small footprint** with **high throughput**, combined with reagents demonstrating **excellent performance** and **stability**, the CP3000 system is an exciting new choice for laboratories of all sizes.






CP3000

Easy-to-Use

- Intuitive touchscreen software, ideal for all levels of operator
- Colored icons show the status of instrument, reagents and samples at-a-glance
- Easy to load barcoded reagents with internal barcode reader for positive reagent identification
- Step-by-step software guided maintenance

Flexible

- Closed tube sampling option
- Able to handle multiple vials and lots of each reagent
- Multiple calibration curves are stored and managed automatically by the software
- User definable settings for repeats, dilutions and reflex testing
- Novel mixing test functionality




Automated Coagulation System
CP3000

High Productivity

- High throughput of up to 400 PT tests per hour
- Small footprint requiring no bench space makes the CP3000 ideal for laboratories of all sizes
- Continuous rack sample loading with easy access to samples when required
- Dedicated STAT lanes, with time to result from just 2.5 minutes
- Capability to link to laboratory automation systems using point-in-space sampling



Efficient

- Enhanced on board stability due to excellent refrigerated reagent carousel
- Tapered reagent vials minimise reagent waste
- Minimal wash solution requirements: 1x 250mL wash solution is all that is required for up to 1 year
- Cuvette loading can occur whilst instrument is in use



With the Quality and Reliability You Expect

- Optical detection system with clotting, chromogenic and immunoturbidimetric measuring channels
 - Light scattering photometric clotting test measurement at 660nm which minimises interference from hemolysis, icterus and lipemia
- Sample assessment checks; sample volume and turbidity checks give more reliable test results
- Multi-level password protection
- 2.5 year data storage





Automated Coagulation System
CP3000

SEKISUI
DIAGNOSTICS

Because every result matters™

THE AMERICAS

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com
ISO 13485 Certified

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent, ME16 0LQ, UK
Phone: +44 (0)1622 607800
Fax: +44 (0)1622 607801
Email: info@sekisui-dx.com
ISO 9001 Certified, ISO 13485 Certified

