

# The FastPack® IP System Specifications



## Analyzer Specs

- Technology: Paramagnetic Particle Chemiluminescence Immunoassay
- Size: 13 x 9 x 12 inches (33 x 23 x 30 cm)
- Weight: 28.64 lbs (12.99 kg)
- Electrical Requirements: 100-240  $\pm$  10% VAC, 1.9A
- Power Consumption: 100 Watts (maximum)
- Ambient Operating Temperature: 15°C (59°F) to 32°C (90°F)
- Operating Humidity: 10% to 80% relative humidity
- Operating Altitude:  $\leq$  2000 m (6562 ft.)
- Display Output: 2-line alphanumeric display

## Printer Specs

- Size: 7.1 x 5.1 x 7.3 inches (18 x 13 x 18.5 cm)
- Weight: 6 lbs (2.7 kg)
- Electrical Requirements: 90-264 V AC, 50/60 Hz (input); 24V DC  $\pm$  5%, 2.0A (output)

### THE AMERICAS

Sekisui Diagnostics, LLC  
4 Hartwell Place  
Lexington, MA 02421  
Phone: 800 332 1042  
Fax: 800 762 6311  
Email: [questions@sekisui-dx.com](mailto:questions@sekisui-dx.com)

### INTERNATIONAL

Sekisui Diagnostics (UK) Limited  
Liphook Way, Allington  
Maidstone, Kent ME16 0LQ, UK  
Phone: +44 1622 607800  
Fax: +44 1622 607801  
Email: [info@sekisui-dx.com](mailto:info@sekisui-dx.com)

[www.sekisuidiagnostics.com](http://www.sekisuidiagnostics.com)

**SEKISUI**  
DIAGNOSTICS

*Because every result matters™*