

Asahi Kasei Pharma Corporation  
Fine Chemicals & Diagnostics Division  
9-1 Kanda Mitoshiro-cho  
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department  
Phone +81-(0)3-3259-5775  
Fax +81-(0)3-3259-5741  
diagnostics@om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41H04-85

April 1, 2005

[ 3<sup>rd</sup> edition ]

Product Name: 3 $\alpha$ -HYDROXYSTEROID DEHYDROGENASE [3 $\alpha$ -HSD II]

(3 $\alpha$ -Hydroxysteroid : NAD oxidoreductase, EC 1. 1. 1 . 50)

Product Number:	T-58
Origin:	<i>Pseudomonas</i> sp.
Description:	White lyophilized powder
Activity:	$\geq 30$ U/mg
Appearance of solution: (5U/ml in buffer solution)	Clear A400nm $\leq 0.010$ A500nm $\leq 0.010$ A650nm $\leq 0.010$
Contaminating enzymes: NADHoxidase	$\geq 90$ % (NADH Residual)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*

*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*