

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41A05-85
April 1, 2005
[4th edition]

Product Name: ADP-HEXOKINASE [ADP-HKP II]
(D-Hexose-6-phosphotransferase, EC 2. 7. 1)

Product Number: T-92

Origin: *Pyrococcus furiosus*

Description: White lyophilized powder

Activity: ≥ 30 U/mg

Appearance of solution: Clear
(1mg/ml in buffer solution)
A400nm ≤ 0.015
A500nm ≤ 0.015
A650nm ≤ 0.015

Contaminating enzymes:
ATPase $\leq 0.01\%$ (U/U)
NADPH oxidase $\leq 0.01\%$ (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T.OTANI
Plant Director
Ohito diagnostics plant



H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant