

Asahi Kasei Pharma Corporation  
Fine Chemicals & Diagnostics Division  
9-1 Kanda Mitoshiro-cho  
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department  
Phone +81-(0)3-3259-5775  
Fax +81-(0)3-3259-5741  
diagnostics@om.asahi-kasei.co.jp

## SPECIFICATION

ISO: 75F41A02-85  
April 7, 2005  
[ 3<sup>rd</sup> edition ]

Product Name: ACYL-COA OXIDASE [ACOD]  
(Acyl-CoA : oxygen oxidoreductase, EC 1. 3. 3. 6)

Product Number: T-17

Origin: *Arthrobacter* sp.

Description: Yellowish lyophilized powder

Activity: 20~40 U/mg

Appearance of solution: Clear  
(1mg/ml in buffer solution) A400nm  $\leq$ 0.080  
A500nm  $\leq$ 0.040  
A650nm  $\leq$ 0.020


Contaminating enzymes:  
Catalase  $\leq$ 1.00 (U/U)  
Glucose oxidase  $\leq$ 0.01% (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.*

*Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



-----  
T. OTANI  
Plant Director  
Ohito diagnostics plant



-----  
H. NAGAI  
Quality control Assistant Manager  
Ohito diagnostics plant