

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41A08-85

April 1, 2005

[3rd edition]

Product Name: ALCOHOL OXIDASE [ALOD]

(Alcohol:oxygen oxidoreductase , EC 1. 1. 3. 13)

Product Number: T-38

Origin: *Candida* sp.

Description: Light yellowish lyophilized powder

Activity: 7~20 U/mg


Appearance of solution: Clear
(12U/ml in buffer solution)
A400nm \leq 0.200
A500nm \leq 0.100
A650nm \leq 0.030

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T.OTANI
Plant Director
Ohito diagnostics plant



H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant