

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41B01-85
November 27, 2006
[1st edition]

Product Name: **BILIRUBIN OXIDASE [BOD II]**

 (bilirubin:oxygen oxidoreductase, EC 1. 3. 3. 5)

Product Number: **T-188**

Origin: *Bacillus subtilis*

Description: Light blue to blue solution

Activity: ≥ 100 U/ml

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Tobe

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Manager
Ohito diagnostics plant