



SPECIFICATION

ISO: 75F41C04-85
April 16, 2013
【 6th edition 】

Product Name:	CHOLESTEROL ESTERASE [CEN] (Sterol-ester acylhydrolase, EC 3. 1. 1. 13)
Product Number:	T-18
Origin:	<i>Pseudomonas</i> sp.
Description:	White to pale brownish lyophilized powder
Activity:	100~250 U/mg
Activity/Protein:	≥ 100 U/mgp
Appearance of solution: (150U/mL in buffer solution)	Clear A400nm ≤0.500 A500nm ≤0.250 A650nm ≤0.125
Contaminating enzymes: Catalase	≤1.0 % (U/U)

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

T. Matsuoka

T.MATSUOKA
Plant Director
Ohito diagnostics plant

Y. Arai

Y.ARAI
Quality control Manager
Ohito diagnostics plant