

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41C03-85

April 1, 2005

[3rd edition]

Product Name: CHOLESTEROL ESTERASE [CEBP-M]
(Sterol-ester acylhydrolase, EC3. 1. 1. 13)

Product Number: T-98

Origin: *Microorganism*

Description: White to off white lyophilized powder

Activity: ≥ 10.0 U/mg


Appearance of solution: Clear
(1mg /ml in buffer solution)
A400nm ≤ 0.030
A500nm ≤ 0.015
A650nm ≤ 0.010

Contaminating enzymes:
Catalase $\leq 1.0\%$ (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



.....
T. OTANI
Plant Director
Ohito diagnostics plant



.....
H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant