

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41C06-85

April 1, 2005

[3rd edition]

Product Name: CHOLINE OXIDASE [COD]
(Choline: oxygen 1-oxidoreductase, EC 1. 1. 3. 17)

Product Number: T-05

Origin: *Arthrobacter globiformis*

Description: Yellowish lyophilized powder

Activity: 8~20 U/mg


Appearance of solution: Clear
(3U /ml in buffer solution)
A400nm \leq 0.050
A500nm \leq 0.030
A650nm \leq 0.015

Contaminating enzymes:
Catalase \leq 10.0 (U/U)
Glucose oxidase \leq 0.01% (U/U)

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*



T.OTANI
Plant Director
Ohito diagnostics plant



H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant