

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41G04-85
July 14, 2006
[1st edition]

Product Name: GLYCEROPHOSPHORYLCHOLINE PHOSPHODIESTERASE [GPCP]

(sn-Glycero-3-phosphocholine glycerophosphohydrolase, EC 3. 1. 4. 2)

Product Number:	T-33
Origin:	<i>Gliocladium roseum</i>
Description:	White to brownish lyophilized powder
Activity:	≥ 10.0 U/mg
Appearance of solution: (50U/ml in buffer solution)	Clear A400nm ≤ 0.100 A500nm ≤ 0.050 A650nm ≤ 0.050

*Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.*

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Manager
Ohito diagnostics plant