

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41G06-85

April 1, 2005

[3rd edition]

Product Name: L- α -GLYCEROPHOSPHATE OXIDASE [GPOM]

(L-Glycero-3-phosphate:oxygen oxidoreductase, EC 1. 1. 3. 21)

Product Number:	T-107
Origin:	
Description:	Yellowish lyophilized powder
Activity:	≥ 15.0 U/mg
Appearance of solution: (50U/ml in buffer solution)	Clear A400nm ≤ 0.150 A500nm ≤ 0.060 A650nm ≤ 0.015

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant