

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41N01-85
April 1, 2005
[3rd edition]

Product Name: NAD synthetase [NADS II]
(Deamido-NAD:Ammonia ligase(AMP-forming), EC 6. 3. 1 . 5)

Product Number: T-67

Origin: *Bacillus stearothermophilus*


Description: White lyophilized powder

Activity: ≥ 1.0 U/mg

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



T. OTANI
Plant Director
Ohito diagnostics plant



H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant