



SPECIFICATION

ISO: 75F41P14-85

June 3, 2010

【 3rd edition 】

Product Name: **PHOSPHOLIPASE A₂ [PLA₂ II L]**
(Phosphatidylcholine 2-acylhydrolase, EC 3. 1. 1. 4)

Product Number: **T-194**

Origin: *Streptomyces violaceoruber*

Description: Colorless to brownish solution

Activity: ≥ 200 U/mL

Appearance of solution: Clear
(10U/mL in buffer solution)
A400nm ≤ 0.040
A500nm ≤ 0.020
A650nm ≤ 0.020

Storage: Store in tightly closed containers, desiccated, protected from light, at 2~8°C.

Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.

M. Toke

M.TOBE
Plant Director
Ohito diagnostics plant

H. Nagai

H.NAGAI
Quality control Assistant Manager
Ohito diagnostics plant