

Asahi Kasei Pharma Corporation
Fine Chemicals & Diagnostics Division
9-1 Kanda Mitoshiro-cho
Chiyoda-ku, Tokyo 101-8481 Japan

Diagnostics Marketing Department
Phone +81-(0)3-3259-5775
Fax +81-(0)3-3259-5741
diagnostics@om.asahi-kasei.co.jp

SPECIFICATION

ISO: 75F41P06-85
April 1, 2005
[3rd edition]

Product Name: PHOSPHOLIPASE-D [PLDP]

(Phosphatidylcholine phosphatidohydrolase, EC 3. 1. 4. 4)

Product Number: T-39

Origin: *Streptomyces* sp.

Description: Off white lyophilized powder

Activity: ≥ 100 U/mg

Appearance of solution:
(100U/ml in buffer solution) Clear


A400nm	≤ 0.100
A500nm	≤ 0.050
A650nm	≤ 0.025

Moisture $\leq 8.5\%$

Storage: Store in tightly closed containers, desiccated, protected from light, at -20°C.
Hazard Statement: Non-Hazardous. And for In Vitro Diagnostic Use.



.....
T. OTANI
Plant Director
Ohito diagnostics plant



.....
H. NAGAI
Quality control Assistant Manager
Ohito diagnostics plant